15



Stirred Reactors and Pressure Vessels



Designing and Building Quality Pressure Apparatus for 120 Years



Founded in 1899 by University of Illinois Professor S.W. Parr, Parr Instrument Company has consistently strived to provide for its customers the very best in product, service and support.



Welcome to the Fifteenth Edition of the Parr Stirred Reactor and Pressure Vessel Catalog. We proudly present here our latest catalog describing the continued expansion of this product line. New products, new designs, and expanded options are offered; all backed by Parr's continuing dedication to product quality, safety, and customer satisfaction. We hope you will find this new catalog to be a useful tool.

211 Fifty-third Street

Moline, Illinois 61265 USA

Phone: 1-309-762-7716

1-800-872-7720

Fax: 1-309-762-9453

Email: parr@parrinst.com

web: www.parrinst.com

1

Table of Contents

Introduction	 	 	2
The User's Responsibility.	 	 	6

1 : Design Features

Design Codes and Certification8
Materials of Construction
Pressure and Temperature Limits 13
High Torque Magnetic Drives14
Split-Ring Closures
Gaskets and Seals 18
Mounting Styles
Warranty

2 : Stirred Reactors

Reactor Selection Guide
How to use Ordering Guide
Series 4520 1-2 L Bench Top Reactors 28
Series 4530 1-2 L Floor Stand Reactors32
Series 4540 600-1200 mL HP Reactors 36
Series 4550 1-2 Gallon Reactors 40
Series 4555 2.6-5 Gallon Floor Stand Reactors
Series 4560 Mini Reactors
Series 4570-4580 HT/HP Reactors52
Series 4590 Micro Reactors56
Series 5100 Glass Reactors60
Series 5500 HP Compact Reactors64

3 : Multiple Reactors

System (MRS)68
Series 2500 Micro Batch Reactor System72
Parallel Reactor Systems74
4 : Tubular Reactors 77
Series 5400 Continuous Flow Tubular Reactors
Fluidized Bed Reactors

67

5 : Specialty and Custom Reactor Systems

Introduction to Specialty and Custom Reactors
Biofuels and Alternative Fuels Research Systems
Horizontal Reactors
Rotating Cage Apparatus92
A1000G Tri-Axial Core Holder93
Supercritical Fluids94
GTL (Gas-To-Liquid) System95
Other Specialty and Custom Reactor Systems 96

99

115

6 : Controllers

Controller Overview100
Series 4848 Reactor Controllers 102
Series 4848T Touchscreen Reactor Controller 105
Series 4838 Temperature Controller 107
SpecView Software for 4848 & 4838 Controllers
Series 4871 Process Controllers 109
Series 4875 Power Controller 113
Series 4876 Power Controllers113
4877 Valve Controller
A2200E Flow Controller114
A2110E Motor Controller114

7 : Optional Accessories

Heaters
Stirring Motors and Drives118
Stirrer Options 120
Gas Entrainment 121
Catalyst Baskets122
Condensers 123
Safety Rupture Discs124
Pressure Relief Valves and Gages 125
Gas Measurement Systems 126
Liquid Charging Systems
Solids/Slurry Addition Devices 128
Cooling Coils
Cylinder Liners
Liquid Sampling130

Sample Collection Vessels	130
4878 Automated Liquid Sampler	130
Filtered Dip Tubes	130
Catalyst Testing System	131
Bottom Drain Valves	131
Valves and Fittings	132
Pressure Hoses	133
Thermocouples	133
Equipment for Use in Potentially Ignitable Atmospheres	
Insulated Electrical Glands	
Spare Parts Kits	137
Temperature Limits	138
External Valves and Fittings	
8: Non-Stirred Vessels 1	39

8: Non-Stirred Vessels

Introduction 140
Series 4700 Pressure Vessels, 22 & 45 mL
Series 4740 HP/HT Pressure Vessels, 25 & 75 mL
Series 4750 General Purpose Pressure Vessels, 125 & 200 mL
Series 4790 Pressure Vessels, 25-100 mL
<u>Series 4760 Pressure Vessels,</u> 100-600 mL152
Series 4600 Pressure Vessel Systems, 1 & 2 Liter 156
Series 4625 High Pressure Vessel Systems, 600 & 1200 mL 160
Series 4650 High Pressure / High Temperature Vessel Systems
Series 4660 Pressure Vessel Systems, 1 & 2 Gallon
Series 4680 High Pressure / High Temperature Vessel Systems, 1 & 1.8 Liters
Series 4670 HP-HT Vessel Systems, <u>1 and 1.5 Gallon168</u>
Series 4675 General Purpose Vessel Systems, 2.6 & 5 Gallon
Gage Block Assemblies
Coned Pressure Fittings

Some of the things that

make Parr Instrument Company the world leader in laboratory equipment cannot be easily shown in a catalog. We refer here to the long-established practices and policies of Parr that have helped to build our reputation as a reliable supplier in our specialized field, such as:





A technical sales and support staff with an average of over 20 years of experience in this specialized field. We have "been there and done that" and we are looking forward to helping you solve your unique requirements.

Communication

A commitment to listen to our customers and a willingness to make the additions or changes in a reactor or pressure vessel that the customer may want or need. Approximately 40 percent of all the reactors and vessels we ship are modified in some way over-and-above the options listed in this catalog. Specialty valves, custom head layouts, electrical leads, special stands, additional openings, unique motors, and non-standard materials of construction are just a few of the modifications we work with every day.





www.parrinst.com is continually updated to provide you with 24 hours access to a more extensive discussion of our products and capabilities. Look for our online resources to expand as we strive to continue to provide you with the world class experience you have come to expect from Parr Instrument Company.

2



nvestment

An investment in modern machine tools and superior production management methods enables us to make available not only all of the options listed on the following pages, but also to handle custom orders for one-of-a-kind designs along with regular production orders on a routine basis.

Service

A commitment to maintain commonly used replacement and service parts in stock for next day shipment.





Research

A commitment to new materials, manufacturing methods, research disciplines, and technological advances to continue our leadership position.

Delivery

A commitment to prompt and dependable delivery

schedules for not only catalog items, but also for custom equipment as well. Typical delivery time for catalog reactors constructed of stainless steel is about six weeks. For reactors made of special alloys or including custom features, lead time will be extended based on the complexity of the customization. Of equal importance over the life of your reactor is our commitment to maintain commonly used replacement and service parts in stock for next day shipment.





At Parr Instrument Company,

we never rest on our past successes. We are constantly working to further improve our products and to streamline our processes. We continue to lead the industry in innovation, safety, and reliability. Researchers from around the world rely on Parr's scientific and engineering design expertise. We are continually expanding our expertise for designing and developing custom systems. From single units to complex multi-reactor systems, we can help find a solution for your research needs.

Our innovative horizontal stirred reactors can be designed to load vertically and operate in either the vertical or horizontal position. The tilt mechanism can be manual or automated. The options are only limited by your design requirements.



Six reactor 4590-based multi-reactor with automated sampling. See page 76.



Temperature, pressure, stirred or non-stirred, options such as windows for measurment or observation, we can build one for you.

Hinged Split-Rings. Available on many of our larger Reactor Systems throughout this catalog.

122

The User's Responsibility

All Parr reactors and pressure vessels are designed and manufactured with great care to ensure safe operation when used within their prescribed temperature and pressure limits. But, the basic responsibility for safety when using this equipment rests entirely with the user; who must:

1. Select a reactor or pressure vessel which has the capability, pressure rating, corrosion resistance and design features that are suitable for its intended use. Parr engineers will be glad to discuss available equipment and material options with prospective users, but the final responsibility for selecting a reactor or pressure vessel that will perform to the user's satisfaction in any particular reaction or test must rest with the user — not with Parr.

In exercising the responsibility for the selection of pressure equipment, the prospective user is often faced with a choice between over or under-designed equipment. The hazards introduced by under-designed pressure vessels are readily apparent, but the penalties that must be paid for over-designed apparatus are often overlooked.

Recognizing these criteria, Parr reactors and pressure vessels are offered in several different styles, each designed for convenient use in daily operation within certain temperature and pressure limits, using gaskets, closures and other elements carefully selected for safe operation within the limits specified for that design. But in order to preserve the validity of these designs, all temperature and pressure limits must be observed, and no attempt should be made to increase these limits by making alterations or by substituting components which are not recommended by Parr Instrument Company.

2. Install and operate the equipment within a suitable barricade, if required, with appropriate safety accessories and in full compliance with local safety codes and rules.

All standard Parr reactors and pressure vessels are provided with either a suitable relief device or a means to attach one (typically in the form of a plugged opening). When a pressure vessel is delivered without a pressure relief device, it is the customer's responsibility to provide pressure relief in order to protect the operator and the equipment from destructive high pressures. If you need more information or need help in selecting a proper relief device, please contact Parr Instrument Company.

- **3. Establish training procedures** to ensure that any person handling the equipment knows how to use it properly.
- **4. Maintain the equipment** in good condition and establish procedures for periodic testing to be sure the vessel remains structurally sound.

6

Design Features Chapter 1



Design Codes Certification Quality Assurance: - ISO 9001: 2015 - ASME - "R" Symbol - PED

- PARR
- <u>CSA</u>
- <u>TPS ACT</u>
- <u>CE</u>

Materials of Construction

- <u>Alloy Designations</u> <u>Pressure and Temperature Limits</u> <u>Magnetic Drives</u> <u>Split-Ring Closures</u> <u>Gaskets & Seals</u> <u>Mounting Styles</u>

Design Codes and Certification



Commitment to a Worldwide Market

In 1973, Parr Instrument Company made a commitment to serve customers on a worldwide basis in a reliable and continuing manner. Parr products were exhibited for the first time at Achema in Frankfurt, Germany that year. An extensive effort was initiated to select and train a network of local agents to provide technical sales and service support as well as import services in each of the countries where modern chemical research and development is conducted.

In 1982, Parr Instrument Company was awarded the President's E Award for Excellence in Exporting in recognition of our successful efforts in this area. In 1988, Parr was awarded the U.S. Department of Commerce's highest level award the President's E Star Award for Excellence in Exporting in recognition of our continuing and accelerating success in servicing the world market.

Today, Parr pressure reaction equipment is in use in over 80 countries and active exclusive agents are operating in 28 of these. Today, over half of all Parr pressure reaction equipment is exported from the USA. As a part of this expansion in the world market, international considerations are designed into, not added onto, all Parr equipment. Design criteria specified in the ASME Code for Unfired Pressure Vessels are closely observed in the manufacture and testing of all Parr pressure equipment. Cylinders for each reactor and pressure vessel of standard design are machined from solid, hot-rolled or forged bars of selected corrosion resistant alloys with no seams or welds in the vessel as potential sources of weakness or localized corrosion. Each individual cylinder is tested hydrostatically. For vessels bearing ASME Certification Mark with "U" Certification designator, the minimum hydrostatic test pressure is 1.3 times the rated working pressure corrected for temperature. For CE marked pressure vessels the minimum hydrostatic test pressure is the higher of 1.43 times the rated working pressure or 1.25 times the rated working pressure corrected for temperature. Each complete reactor and vessel assembly is further tested with nitrogen to be sure that it is leak-free and operating properly.

Protection against equipment damage and possible personal injury in case of an accidental over-pressurization is typically provided by a safety rupture disc installed in the head of a reactor or in a gage block assembly mounted on the head. A description of these safety discs and rules to be observed in their selection and use are given on page 124.

Quality Assurance

Parr Instrument Company has designed, installed and operates under a Quality Assurance Program which ensures that all aspects of the design, materials selection and procurement, manufacture, testing and certification of its pressure vessels are performed in accordance with accepted codes and practices. Currently this Quality Assurance Program has been certified to be in compliance with ISO 9001: 2015 for Total Quality Procedures, ASME, CSA (Canadian Standards Association), and all applicable European Directives (CE) which include but are not limited to PED (Pressure Equipment Directive), LVD (Low Voltage Directive), EMC (Electromagnetic compatibility), and Machinery Safety.

ISO 9001: 2015 Certification

Parr Instrument Company's overall Quality Assurance System has been certified to be in compliance with ISO 9001: 2015 by TÜV SÜD. ISO 9001: 2015 covers the overall quality assurance and management compliance aspects of Parr's activities as opposed to the certification of an individual product.





8

ASME Certification

Parr Instrument Company holds a Certificate of Authorization issued by the Boiler and Pressure Vessel Committee of the American Society of Mechanical Engineers (ASME) and the National Board of Boiler & Pressure Vessel Inspectors (NBBI) as an approved facility for manufacturing unfired pressure vessels. If required, any Parr reactor or pressure vessel can be certified to the ASME BPVC, Section VIII, Division 1, which involves:

- Vessel inspection and tests by an Authorized Inspector of the National Board of Boiler and Pressure Vessel Inspectors.
- Issuance of Form U-1 or U-1A, "Manufacturer's Data Report for Pressure Vessels", which provides all relevant documentation of the finished vessel including certification of the chemical analysis and physical properties of materials used in the vessel.
- Application of the ASME Certification Mark with "U" Certification Designator to the vessel by the Inspector.
- 4. Registration of the vessel with the National Board.

It should be noted that unless specifically requested, ASME certification provided by Parr will cover the vessel only and not the relief device. Parr can furnish certain relief devices with ASME certification upon request.

Parr is also certified under Section IX of the ASME Code for welding. Normally welding is done only to attach jackets or fittings to the primary vessel.

Parr Instrument Company is authorized to apply the "R" Symbol in accordance with the provisions of the National Board of Boiler and Pressure Vessel Inspectors for metallic repairs and/or alterations at its Moline, Illinois headquarters.

PED Certification

Parr Instrument Company has implemented, operates and maintains a quality assurance system as described in the Pressure Equipment Directive Annex III, Module H/H1 for the scope of Design and Manufacture of Pressure Reactors and Assemblies for Laboratory Applications.

Parr Certification

If requested, Parr will furnish a signed certificate listing the materials of construction used in the manufacture of an individual reactor or pressure vessel, the pressure tests applied to that reactor, material certificates and reference layout drawings.

CSA Certification

Where appropriate, Parr reactors are manufactured and certified to the electrical code established by the Canadian Standards Association. Identification of those units for which CSA Certification has been received will be provided upon request. The CSA logo is shown on the nameplate of each CSA certified unit.

TPS ACT Certification (CUE Mark)

Where appropriate, Parr reactors will carry the CUE mark of TPS (TÜV SÜD Product Service division) ACT (Acceptance of Client's Testing) certifying the manufacture and testing in compliance with mutually agreed product categories and/or standards, including the requirements of the current revision of ISO/IEC 17025, appropriate to applicable Canadian, US, and European market regulations.

CE Certification

Where appropriate, Parr reactors will carry the CE Mark certifying compliance with all applicable European Community Directives.

Other National or Local Codes

Parr regularly works with other national, state or international authorities to obtain individual approval for specific vessels. Parr has obtained CRN Approval for pressure vessels in Canadian provinces. The internationally recognized Quality Assurance Program in place at Parr and the experience of the Parr engineering department in working with these authorities makes it possible to obtain these approvals with little difficulty. It is the user's responsibility to identify any such applicable code so that these requirements can be met before the vessel is fabricated and delivered. Parr's network of international distributors are familiar with the applicable codes for pressure vessels within their countries of responsibility.





Materials of Construction



Parr reactors are normally made of Type 316 Stainless

Steel, but they can also be made of other alloys as well. Available construction materials include:

- Type 316/316L Stainless Steel
- Alloy A-286
- Alloy 20
- Alloy 230
- Alloy 400
- Alloy 600
- Alloy 625
- Alloy B-2/B-3
- Allov C-276
- Nickel 200
- Titanium Grades 2, 3, 4, & 7
- Zirconium 702 & 705
- Other material may be available upon request

Alloy Designation

Parr uses alloy designation numbers to identify the various corrosion resistant alloys available for use in Parr reactors and pressure vessels. These allovs can also be identified by trade names and by ASTM, ASME, DIN and other specification numbers. Many of the high nickel alloys were originally patented and sold under trade names, such as Monel¹, Incoloy¹, Carpenter Alloy 20², Hastelloy³, etc. Most of the original patents have expired and these alloys are now available from other suppliers, as well as from the owners of the original trade names.

Among the many corrosion resistant alloys now available, there may be two or three with very similar compositions and intended for use in the same corrosive environment. In these cases, Parr will select and offer the most widely used alloy in each of the basic corrosion resistance categories, rather than catalog and stock all three.

Each of these alloys has its own physical strength and temperature characteristics as well as its own unique resistance to certain corrosive materials. All of these factors must be considered when making a selection, with cost and availability also becoming factors in the final choice.

The basic composition of these alloys is listed in Table I. Corrosion resistance information can be obtained from various corrosion handbooks and metallurgical publications. Helpful information can also be obtained from the individual alloy manufacturers.

Any abridged listing of corrosion resistance of various metals and alloys can be potentially misleading since it cannot possibly deal with all of the effects of concentration, temperature, pressure and the presence of additional ions, all of which have a significant effect upon the ability of a reactor to withstand corrosion. In addition, the vulnerability

MONEL, INCONEL and INCOLOY are Registered Trademarks of Special Metals Corp. ² CARPENTER 20 is a Registered Trademark of Carpenter Technology Corporation. ³ HASTELLOY is a Registered Trademark of Haynes International, Inc.

Parr Instrument Company will attempt to answer questions regarding corrosion resistance and mav suggest materials that might be suitable for specific applications, but no guarantee can be made that any particular alloy will be fully resistant to a prescribed set of corrosive conditions.

of any material to stress corrosion cracking, intergranular corrosion and pitting must also be considered when judging the suitability of a material for a particular application.

The principal characteristics of the construction materials offered by Parr are summarized on the following pages. These listings are intended to serve only as a starting point for any study of comparative corrosion resistance and physical properties. Material manufacturers booklets on each alloy are available on our website at www.parrinst.com. Additional details may also be obtained from other sources.

Type 316/316L Stainless Steel

Type 316 Stainless Steel is an excellent material for use with most organic systems. A few organic acids and organic halides can, under certain conditions, hydrolyze to form inorganic halogen acids which will attack T316SS. Acetic, formic and other organic acids are routinely handled in T316SS. T316SS is not normally the material of choice for inorganic acid systems. At ambient temperatures it does offer useful resistance to dilute sulfuric, sulfurous, phosphoric and nitric acids, but sulfuric, phosphoric and nitric acids readily attack T316SS at elevated temperatures and pressures. Halogen acids attack all forms of stainless steel rapidly, even at low temperatures and in dilute solutions.

Halogen salts can cause severe pitting in all stainless steels. Chlorides can cause stress corrosion cracking, but many other salt solutions can be handled in stainless vessels, particularly neutral or alkaline salts.

At moderate temperatures and pressures, T316SS can be used with most commercial gases. In scrupulously anhydrous systems even hydrogen chloride, hydrogen fluoride and chlorine can be used in stainless steel. Essentially all of the T316SS produced today also meets the specifications for 316L, low carbon stainless steel.

Alloy 20

Alloy 20 is an enriched grade of stainless steel, designed specifically for use with dilute (up to 30 percent by weight) sulfuric acid at elevated temperatures. It can also be used for nitric and phosphoric acid systems as well as for all systems for which T316SS is suitable.

Alloy 400

Alloy 400 is an alloy comprised essentially of twothirds nickel and one-third copper. For many applications it offers about the same corrosion resistance as nickel, but with higher maximum working pressures and temperatures and at a lower cost because of its greatly improved machinability.

Alloy 400 is widely used for caustic solutions because it is not subject to stress corrosion cracking in most applications. Chloride salts do not cause stress corrosion cracking in Alloy 400. It is also an excellent material for fluorine, hydrogen fluoride and hydrofluoric acid systems. Alloy 400 offers some resistance to hydrochloric and sulfuric acids at modest temperatures and concentrations, but it is seldom the material of choice for these acids. As would be expected from its high copper content, Alloy 400 is rapidly attacked by nitric acid and ammonia systems.

Continued, next page >

Table I Nominal Chemical Composition of Pressure Vessel Materials

Major Elements (Percent)							
Material	Typical Trade Name	Fe	Ni	Cr	Мо	Mn	Other
T316	Stainless Steel	65	12	17	2.5	2.0	Si 1.0
Alloy 20	Carpenter 20	35	34	20	2.5	2.0	Cu 3.5, Cb 1.0 max
Alloy 230	Haynes 230	3	52	22	2	0.7	Co-5, W-14, Si-0.5
Alloy 400	Monel 400	2.5 max	63 min			2 max	Cu 31.5
Alloy 600	Inconel 600	8	76	15.5			
Alloy 625 Gr 1	Inconel 625 Gr 1	5	58	22	9	0.5	Nb+Ta 3.7
Alloy A-286	Alloy 286	53	25	15	1.2	2	Si-1, Ti-2.1, Al-0.35
Alloy B-2/B-3	Hastelloy B-2/B-3	2	66	1	28	1	Co 1.0
Alloy C-276	Hastelloy C-276	6.5	53	15.5	16	1	W4.0, Co 2.5
Nickel 200			99				
Titanium Grade 2, 3, 4	Commercially pure titaniumn						Ti 99 mi
Titanium Grade 7							Ti 98.8 min, 0.15 Pd
Zirconium Grade 702							Zr + Hf 99.2 min, Hf 4.5 max
Zirconium Grade 705							Zr + Hf 95.5 min, Hf 4.5 max, Nb 2.5

All Parr raw

Materials of Construction

< Continued, from previous page

Alloy 600

Alloy 600 is a high nickel alloy offering excellent resistance to caustics and chlorides at high temperatures and high pressures when sulfur compounds are present. In caustic environments, Alloy 600 is unexcelled. It also is often chosen for its high strength at elevated temperatures. Although it can be recommended for a broad range of corrosive conditions, its cost often limits its use to only those applications where its exceptional characteristics are required.

Alloy B-2/B-3

Alloy B-2/B-3 is an alloy, rich in nickel and molybdenum, which has been developed primarily for resistance to reducing acid environments, particularly hydrochloric, sulfuric and phosphoric. Its resistance to these acids in pure forms is unsurpassed, but the presence of ferric and other oxidizing ions in quantities as low as 50 ppm can dramatically degrade the resistance of this alloy.

Alloy C-276

Alloy C-276 is a nickel chromium-molybdenum alloy having perhaps the broadest general corrosion resistance of all commonly used alloys. It was developed initially for use with wet chlorine, but it also offers excellent resistance to strong oxidizers such as cupric and ferric chlorides, and to a variety of chlorine compounds and chlorine contaminated materials. Because of its broad chemical resistance, Alloy C-276 is the second most popular alloy, following T316SS, for vessels used in research and development work.

Nickel 200

Nickel 200 is one of the designations of commercially pure nickel. It offers the ultimate in corrosion resistance to hot caustic environments, but its applications are severely restricted because of its poor machinability and resultant high fabrication costs.

Titanium

Titanium is an excellent material for use with oxidizing agents, such as nitric acid, aqua regia and other mixed acids. It also offers very good resistance to chloride ions. Reducing acids, such as sulfuric and hydrochloric, which have unacceptably high corrosion rates in their pure form can have their corrosion rates in titanium reduced to acceptable levels if relatively small quantities of oxidizing ions, such as cupric, ferric, nickel or even nitric acid are present to act as corrosion inhibitors.

This phenomenon leads to many successful applications for titanium in the hydrometallurgy field where acids, particularly sulfuric acid, are used to leach ores. In these operations, the extracted ions act as corrosion inhibitors.

Prospective users must remember that titanium will burn vigorously in the presence of oxygen at elevated temperatures and pressures. While there have been many successful applications in hydrometallurgy where oxygen and sulfuric acid are handled in titanium equipment, the danger of ignition is always present and must be protected against whenever titanium and oxygen are used together.

Commercially pure titanium is available in several grades. Grade 2 is the material most commonly used for industrial equipment since it can be fabricated by welding and can be used to make vessels compliant to the PED and ASME codes. Grade 4, which has slightly higher trace levels of iron and oxygen, has higher strength than Grade 2 but it is not suitable for welding and it is not covered by the PED or ASME Codes.

Since most Parr vessels are not welded, Grade 4 can be used to obtain higher working pressures than can be obtained with Grade 2. Grade 7, containing small amounts of palladium, and Grade 12 containing small amounts of nickel and molybdenum, offer enhanced resistance to certain environments and can be used for Parr reactors and pressure vessels if suitable billets can be obtained.

Zirconium

Zirconium offers excellent resistance to hydrochloric and sulfuric acids, however, as with Alloy B-2/B-3, oxidizing ions such as ferric, cupric and fluorides must be avoided. Zirconium also offers good resistance to phosphoric and nitric acids, and to alkaline solutions as well. Two different grades are available: Grade 702 which contains hafnium is the standard commercial grade offering the best resistance to most corrosive agents. Grade 705 contains small amounts of both hafnium and niobium which increases the strength characteristics and allows for higher maximum working pressures for a vessel. Grade 702 typically offers better corrosion resistance than Grade 705. Grade 702 is also more widely available from commercial stocks of raw materials.

High Temperature / High Strength Alloys

In addition to the metals chosen for their corrosion resistance Parr also offers some alloys that are selected for their outstanding strength values, their high temperature strengths, or both.

Alloy 230

Alloy 230 is an alloy approved for ASME pressure vessel design for temperatures up to 980 °C. It is an alloy high in nickel, chromium, tungsten, and cobalt. While it has resistance similar to Alloy 600, it is normally selected for its high strengths at very high temperatures. It is sometime selected as a bolting material.

Alloy 625

Alloy 625 is an alloy with chemical resistance similar to Alloy C-276, but with much greater strength. We use this alloy to obtain additional pressure ratings for high temperature applications.

Alloy A-286

Alloy A-286 is an alloy of the Stainless Steel family with very high strengths up to 371 °C. It is commonly used as a bolting material.

Pressure and Temperature Limits

The maximum pressure and temperature at which any reactor or pressure vessel can be

used will depend upon the design of the vessel, its material of construction, and other components integral to its design. Since all materials lose strength at elevated temperatures, any pressure rating must be stated in terms of the temperature at which it applies. The listings shown in this catalog show the maximum allowable working pressure (MAWP) for each vessel in pounds per square inch (psi) and in bar at the maximum rated temperature for that particular design when that vessel is constructed of Type 316 Stainless Steel. Maximum pressure and temperature limits for vessels constructed of other alloys are computed and assigned by the Parr Engineering Department in accordance with all applicable regulations.

Lower operating temperatures sometimes permit higher working pressures. For example, the 4560HT High Temperature reactors are rated at 2000 psi (138 bar) maximum pressure and 500 °C maximum temperature. Standard 4560 reactors are rated at 3000 psi (200 bar) maximum pressure at 350 °C maximum temperature.

One should not assume that any vessel being operated at a lower temperature can be used at pressures exceeding the rated MAWP. Factors other than the material strength of the vessel wall may well be the constraint controlling the rating. Other factors that can limit the pressure and temperature ratings are the closures design, the magnetic drive, the type of seal, the choice of other components used, as well as the material of construction. The maximum operational temperature of some materials is much lower than what is permissible with stainless steel as shown in Table II. Users are encouraged to contact the Parr Customer Service Department with any questions. Multiple factors are involved in safely calculating the maximum working pressures and temperatures of Parr Reactors and Pressure Vessels. Please contact Parr Customer Service Department for more information and to assist you in making the correct purchase decision.

Table II

Maximum Allowable Temperatures


Materials of Construction	Maximum Temperature
T316/316L Stainless Steel	800 °C
Alloy 20	426 °C
Alloy 230	980 °C
Alloy 400	482 °C
Alloy 600	625 °C
Alloy 625 Gr 1	648 °C
Alloy A-286	371 °C
Alloy B-2/B-3	426 °C
Alloy C-276	625 °C
Nickel 200	315 °C
Titanium Grade 2, 3, 4, 7	315 °C
Zirconium Grade 702, 705	371 °C



Maximum pressure and temperature ratings for vessels constructed of alloys other than T316SS will be computed and assigned by the Parr Engineering Department.



High Torque Magnetic Drives



All Parr stirred reactors are equipped with a magnetic drive to provide a trouble-free linkage to an internal stirrer, thereby avoiding the leakage problems which can arise with a packed gland stirrer drive. With a Parr magnetic drive there are no rotating seals. The drive turns freely and the system remains gas-tight, permitting long, continuous runs at pressures up to 5000 psi (345 bar) with little or no attention to the seal and drive.

Parr drives are assembled with specially designed permanent magnets which have excellent temperature stability and can be depended upon to operate for long periods with little or no flux degradation. Magnets for the inner rotor to which the stirrer shaft is attached are enclosed in a stainless steel (or other alloy) housing, permanently sealed by welding and supported by graphite-filled, PTFE bushings to provide a long life, chemically inert stirring system. Magnets for the outer drive are also fully enclosed and supported by twin, high quality sealed ball bearings for smooth operation and long life. A water cooling sleeve attached to each drive protects the components from excessive heat arising from the reactor.

Four Sizes

Parr magnetic drives are made in four sizes (as shown in the chart on the facing page), designed to match the full range of Parr reactor sizes and to provide alternate drives for high viscosity loads, higher stirring speeds and other special requirements. Each drive is assembled in a sealed housing which threads directly into the reactor head.

Parr uses neodymium-iron-boron magnets with 25% more coupling force than samarium-cobalt magnets. With very few exceptions involving gear reduction drives, the magnetic stirrers fitted to reactors have higher coupling torques than the stall conditions of standard motors. Today magnetic drives are used with confidence for high viscosity polymerization reactions.

Parr magnetic drives are supported with three graphite-filled PTFE bushings and quality internal ball bearings. They routinely deliver 2000 hours of operation without service.



The General Purpose and Heavy Duty models are the standard units normally furnished with the reactor sizes listed in the adjoining table. The Extra Heavy Duty model is a special high torque drive intended primarily for heavy loads and high viscosity applications. When it is used to replace a standard drive, the standard motor and drive system may have to be modified to provide the higher torque which the Extra Heavy Duty drive is capable of transmitting.

The 5500 Series Compact Reactors have a smaller magnetic drive that is intended for low viscosity applications.

Two Styles Available

Parr offers a choice of two styles of magnetic drives. The first operates with small diameter stirrer shafts which require a lower support bracket or "foot" bearing to stabilize the stirrer shaft. These drives are intended for high speed stirring for applications involving liquid-liquid or gas-liquid mixing.

The second are footless magnetic drives and employ a larger diameter stirrer shaft

Parr Magnetic Drive Series

designed to operate without this lower support bracket or "foot" bearing. These drives were originally designed for digesting ores where the abrasive solids would get caught in the PTFE foot bearing and wear away the stirrer shaft. They are also recommended for slower-speed mixers such as the anchor, paddle, or spiral stirrers. It is also important to select the appropriate motor which can handle the increased drag associated with the larger diameter shafts.

Alternate Packed Gland Drive

For rare circumstances where a direct mechanical drive is preferable to a magnetically coupled system, Parr can furnish a self-sealing packed gland which will maintain a reliable seal on the stirrer shaft at working pressures up to 2000 psig (138 bar). These glands are made to a Parr design which uses a combination of cones and O-rings in conjunction with pressure from within the vessel to maintain a positive seal on the rotating shaft.

Magnetic Drive	Coupling Torque (in-lb)*	Shaft Diameter (inches)	Foot Bushing Required? (see text above)	Ordering Guide Abbreviation	Supplied as Standard Mag Drive on Reactors with these Volumes:	
General Purpose	16	3/16"	Yes	М	25 mL to 2 L	
	16	3/8"	No (Footless)	FMD1		
Heavy Duty	60	3/8"	Yes	HD	1 to 5 Gallon	
	60	5/8"	No (Footless)	FMD2		
Extra Heavy Duty	120	3/8"	Yes	XHD		
	120	5/8"	No (Footless)	FMD3	1 to 5 Gallon	
Compact	2.5	3/16"	Yes	N/A	5500 Compact Reactor	

* in-lb = 0.11 Newton Meter

Please refer to the Parr Magnetic Drives Operating Manual (234M) at https://www.parrinst.com/download/37494



Split-Ring Closures



Split-Ring with Compression Bolts for Moveable Vessels



Split-Ring with Captive Compression Ring and Bolts for Fixed Head Vessels



Split-Ring for Self-Sealing O-ring Closures

Split-Ring U.S. Patent No. 2625296

We have all heard the story of the man who resigned from the Patent Office around 1900 because he believed that everything that could be invented already had been. That is not the philosophy that has driven Parr Instrument Company over its more than a century of service to the chemical industry. The following prediction first appeared in February 1930 in the "Forward" of our bulletin describing our catalytic hydrogenation apparatus.

"The greatest chemical discoveries are yet to be made, untold numbers of which will result from the new studies in catalysis and catalytic reactors. These studies have already unearthed a mine of information, but with a more exact knowledge of catalytic processes, much more will be learned in the future."

We are delighted with the foresight of our predecessors here at Parr and are dedicated to continuing this outlook for the future as we conduct our second century of service to our friends and customers in this critical industry.

Easy Access to Pressure Vessels

Parr reactors and pressure vessels are equipped with a unique split-ring cover clamp which adds greatly to the


convenience of the equipment and the ease with which it can be handled. This is an exclusive Parr design (see sidebar) which allows easy access to a pressure vessel without using a heavy screw cap, cumbersome cover clamps, or a wide flange for moveable bolts. Instead, the head is clamped to the cylinder by either a Type 4140 steel or a Type 316 Stainless Steel ring which has been split into two sections. These sections slide into place from the sides without interfering with any fittings attached to the head.

Self-Sealing and Compression Bolt Designs

For reactors and vessels in which a self-sealing O-ring is used as the main head seal, there are no compression bolts in the split-ring sections. The vessel is closed by simply sliding the two ring sections into place and locking them with either an encircling drop band or with attached, quick-opening latches. This convenient closure can be used on most small and mid-sized stirred reactors and pressure vessels, provided the intended operating temperature does not exceed the allowable working temperature limit for the O-ring seal. Self-sealing O-rings are not recommended for low pressure applications. See <u>page</u> 18 for additional information.

Split-ring closures for reactors with PTFE, flexible graphite, metal or other contained, flat gaskets have a set of compression bolts in the rings which must be tightened to develop the compressive force required to seal the gasket. These split-rings are locked together with either a drop band or quick opening latches.

With either style, the closure parts come completely away from the vessel so that the head with all of its fittings can be lifted from the cylinder or the cylinder can be dropped away without disrupting any attached fittings. When compression bolts are used, they are simply tightened or loosened, they are never completely removed from the split-ring or drop band. This saves time both in opening and closing the vessel and in looking for lost parts.




Split-Ring with Compression Bolts and Drop Band

Split-ring closures add many attractive features to Parr reactors and pressure vessels.

- The reactor or vessel can be opened and closed without disturbing any connections or fittings attached to the head.
- The full inside diameter of the vessel is exposed when the head is removed.
- A maximum area is exposed on the head for attaching valves and fittings.
- There are no cumbersome bolt flanges or threaded studs to interfere with operations.
- No delicate threads on the cylinder to gall or to be damaged in handling.

Screw Cap Closures

Parr uses screw cap closures on small vessels where enough sealing force can be developed by simply tightening the main screw cap. This design can be made more compact than the split-ring closure and is used primarily on general purpose vessels with volumes of less than 100 mL.



Split-Ring with Drop Band and No Compression Bolts




Screw Cap Closure


Split-Ring with Compression Bolts, Latches, and No Drop Band 0 0 0

Hinged Split-Rings


17

Gasket and Seals






O-ring FKM Seal



Flat Flexible Graphite Gasket



High Pressure Metal Gasket

There are four different types of gasket materials for the main head seal in Parr reactors

and pressure vessels, each with its own advantages and limitations. Some of these are recent additions which have significantly expanded the choices a user can consider when selecting a closure and gasket material for the intended operating conditions.

Confined and Contained Flat PTFE Gaskets for Temperatures to 350 °C

The traditional and most popular main head gasket for Parr vessels is a flat gasket made of a PTFE fluoropolymer. In Parr flat gasket closures, the gasket is held in a recess in the vessel cover. The mating lip on the cylinder closes the recess, leaving the gasket completely confined with only a small inside edge exposed to the reactants within the vessel. This combination of complete gasket containment and the exceptional properties of PTFE materials produces a reliable closure for working temperatures up to 350 °C.

Flat contained gaskets require an initial loading pressure in order to develop and to maintain a tight seal. In Parr designs this is produced by tightening a ring of compression bolts in a split-ring cover clamp. Fortunately PTFE is slightly "plastic" and will flow under pressure, producing a seal that improves with each use as the gasket is forced into the faces on the head and cylinder. It also is a very forgiving seal which does not require the special care needed to achieve a uniform loading, which is essential when working with a metal or other non-plastic gasket material.

An equally important advantage of the PTFE gaskets is their essentially universal chemical resistance.

Self-Sealing O-rings

Parr has greatly expanded its offerings of reactors and vessels which feature self-sealing O-ring closures. In these designs the sealing force on the gasket is developed from pressure within the vessel itself, eliminating the need for compression bolts in the split-ring to pre-load the seal. In these self-sealing closures the split-ring sections simply lock the head and cylinder together.

Operating reactors and vessels that incorporate a self-sealing O-ring at pressures less than ~20 percent of their maximum working pressure will need to be modified since diminished force is available to fully energize the seal.

Users who select the self sealing O-ring design must consider two important characteristics of elastometric materials. First, they will not withstand operating temperatures as high as the PTFE gaskets. Secondly, none of these materials offers the universal chemical resistance of PTFE polymers. The chemical resistance is especially important since the O-ring is directly exposed to the contents of the vessel.

Although there are a number of available O-ring materials, the real choice comes down to two. Fluoroelastomer (FKM) O-rings, such as Viton, are a first choice for Parr self-sealing closures. They have good chemical resistance and a working temperature up to 225 °C. Perfluoroelastomer (FFKM) O-rings, such as Kalrez[®], have extremely broad chemical resistance and can be used at working temperatures up to 300 °C. Unfortunately, this material should probably be considered an "exotic" because it costs approximately 80 times



O-ring Main Head Seal with Retaining Lip

as much as an FKM O-ring. While it will raise the allowable working temperature to 300 °C, as a practical matter, most users intending to work at this temperature level would be well advised to choose a closure with a flat PTFE gasket and a 350 °C temperature limit. Other exotic O-ring materials are available, and there are economically priced materials such as ethylene-propylene that will resist some materials that cause FKM to fail, with only slight sacrifices in operating temperatures.

Contained Flat Flexible Graphite Gaskets for Temperatures to 600 °C

For operating temperatures above 350 °C, Parr uses a flexible form of graphite, called Grafoil[®], which has proven to be an excellent high temperature sealing material. It consists of flexible layers of graphite bonded together to produce a gasket that is almost as easy to seal as a flat, PTFE gasket, but with an almost unlimited temperature range and excellent chemical resistance.

Parr has converted all of its standard designs to accept a flat, Grafoil gasket whenever operating temperatures above 350 °C are required, replacing the metal gaskets formerly used for high temperatures. These flexible graphite gaskets are held in grooves identical to the ones used for PTFE gaskets and sealed with the same split-ring closures. This makes it possible to substitute a PTFE gasket whenever the vessel is to be used at temperatures below 350 °C. Grafoil gaskets can be reusable, but their service life is shorter than can be obtained with a PTFE gasket.

Metal Gaskets

Metal gaskets have traditionally been the only gaskets available for use at temperatures above 350 °C. Parr has designs for diamond cross-section metal gaskets which can be furnished for special applications, but we would recommend the flexible graphite gaskets described above for most applications.

Trademarks of Sealing Materials

A number of gasket materials have so dominated their product categories that their Trade Names have become more common than the actual material designation itself. In an attempt to respect the value of these Trade Names and their proper usage and to minimize the disruptions in our descriptions, we have adopted the following generic material descriptions and designations for use in this catalog. Where available we have selected the ASTM material designation.

Common or Trade Name	Material	Designation
Viton®	fluoroelastomer	FKM
Kalrez®	perfluoroelastomer	FFKM
Teflon®	tetrafluoroethylene polymer	PTFE
Grafoil ®	flexible graphite	FG

Viton®, Kalrez®, Teflon® are Registered Trademarks of DuPont. Grafoil® is a Registered Trademark of UCAR Carbon Co. Inc.



Two Mounting Styles



Model 4523 Bench Top Reactor, 1 liter, Fixed Head

Fixed Head Reactors

Parr offers most of its laboratory reactors in a fixed head design. This includes all reactors with volumes from 25 mL to 20 liters. In these reactors the head of the vessel may remain fixed in the reactor support stand. All attachments to the head: gas and liquid feed and discharge lines, cooling water, vapor take-off and condenser, thermocouple and any electrical leads can remain permanently in place. The reactor is opened by simply removing the split-ring and lowering the cylinder away from the head, leaving all of the attachments undisturbed. If desired the complete vessel assembly can be removed from the head support plate.

The support frames for fixed head reactors vary based on the sizes and weights of the vessels:

- For vessels with volumes of 25 mL to 600 mL, the cylinder may be lowered and removed manually by the operator
- Vessels with volumes of 1 to 2 liters include a manual lift mechanism to raise or lower the cylinder
- Vessels with volumes greater than two liters or where weight dictates include a pneumatic cylinder lift mechanism.




Close up view of a fixed head reactor with the vessel and heater lowered

83

Moveable Vessel Reactors

As an alternate to the fixed head designs described on the previous page, all Parr reactors can be furnished in designs which allow the entire vessel to be removed as a complete assembly from the support stand allowing for charging, product recovery, and cleaning. In the smaller and mid-sized models the entire vessel is simply lifted out of the stand by hand. In the larger models, the vessel can be opened and closed with the cylinder remaining in the heater, but the head must be lifted off manually. In the largest (10 Liter and 20 Liter) models, a chain hoist is provided for lifting the heavier head and cylinder.

These moveable reactors will be attractive to users that intend to operate fairly simple batch systems rather than continuous flow arrangements, and that want to be able to remove the vessel for product recovery, charging or cleaning. They will also appeal to those that need to prepare the vessel in a special atmosphere, or want to clean the cylinder and stirrer at a site away from the reactor stand and heater. There is also the added advantage of being less expensive than the fixed head models since they do not require the more elaborate head or the cylinder and heater lift mechanism.



High Pressure Moveable Vessel, 600 mL, for use to 5000 psi



Parr Instrument Company Warranty

Parr Instrument Company (Parr) combustion bombs, calorimeters, reactors, pressure vessels and associated products are designed and manufactured only for use by or under the direct supervision of trained professionals in accordance with specifications and instructions for use supplied with the products. For that reason, Parr sells only to professional users or distributors to such users. Parr produces precision equipment and associated products which are not intended for general commercial use.

Exclusive Warranty. To the extent allowed by law, the express and limited warranties herein are the sole warranties. Any implied warranties are expressly excluded, including but not limited to implied warranties of merchantability or fitness for a particular purpose.

Warranty Conditions:

- **1. Non-assignable.** The warranties herein extend only to the original purchaser/user and to the distributors to such users. These warranties or any action or claims based thereon are not assignable or transferable.
- 2. Use of product. The warranties herein are applicable and enforceable only when the Parr product:
 - (a) Is installed and operated in strict accordance with the written instructions for its use provided by Parr.
 - (b) Is being used in a lawful manner.
 - (c) Has not been modified by any entity other than Parr Instrument Company.
 - (d) Has been stored or maintained in accordance with written instructions provided by Parr, or if none were provided, has been stored and maintained in a professionally reasonable manner.
- **3. The user's responsibility.** Parr engineers and sales personnel will gladly discuss available equipment and material options with prospective users, but the final responsibility for selecting a reactor, pressure vessel or combustion bomb which has the capacity, pressure rating, chemical compatibility, corrosion resistance and design features required to perform safely and to the user's satisfaction in any particular application or test must rest entirely with the user not with Parr. It is also the user's responsibility to install the equipment in a safe operating environment and to train all operating personnel in appropriate safety, operational and maintenance procedures.
- **4. Warranty period.** Unless otherwise provided in writing by Parr, the warranties herein are applicable for a period of one year from date of delivery of the product to the original purchaser/user. Note, however, that there is no guarantee of a service life of one year after delivery.
- **5. Notification.** To enforce any express warranty created herein, the purchaser/user must notify Parr in writing within thirty (30) days of the date any defect is detected. Upon request of Parr, the part or product involved must be returned to Parr in the manner specified by Parr for analysis and non-destructive testing.

Express Warranties.

Subject to the above Conditions, Parr expressly warrants that its products:

- 1. Are as described in the applicable Parr sales literature, or as specified in Parr shipping documents.
- 2. Will function as described in corresponding Parr sales bulletins or, for specially engineered assemblies, as stated in the sales proposal and purchase agreement.
- 3. Will remain free from defects in materials and workmanship for the Warranty Period.

Limitations on the Parr Warranty.

As to the original purchaser/user and to the distributors to such users, Parr limits its liability for claims other than personal injury as follows:

- 1. Replacement or repair. With respect to express warranties herein, Parr's only obligation is to replace or repair any parts, assemblies or products not conforming to the warranties provided herein.
- 2. Disclaimer of consequential damages. In no event shall Parr be liable for consequential commercial damages, including but not limited to: damages for loss of use, damages for lost profits, and damages for resulting harm to property other than the Parr product and its component parts.

Indemnity and Hold Harmless.

Original purchaser/user agrees to indemnify and hold Parr harmless for any personal injuries to original purchaser/user, its employees and all third parties where said injuries arise from misuse of Parr products or use not in accordance with specifications and instructions for use supplied with the Parr products.

Stirred Reactor Systems Chapter 2



Inside this chapter you will find:

4520 Bench Top, 1000 & 2000 mL

4520 HP Bench Top, 970 & 1900 mL

4530 Moveable Cart or Floor Stand, 1000 & 2000 mL

4530 HP Moveable Cart or Floor Stand, 970 & 1900 mL

4540 HP Bench Top, Floor Stand, or Moveable Cart, 600 & 1200 mL

4550 Moveable Cart or Floor Stand, <u>1 & 2 Gallon (3.75 & 7.99 L)</u>

4555 Floor Stand, 2.6 & 5 Gallon (10 & 18.75 L)

4560 Mini, Bench Top, 100-600 mL

4560 HT Mini, Bench Top, <u>300-600 mL</u>

4570 HP/HT, Bench Top, Moveable Cart, or Floor Stand, 250-1800 mL

4580 HP/HT, Moveable Cart or Floor Stand, <u>1 & 1.5 Gallon (3.75 & 5.5 L)</u>

4590 Micro, Bench Top, 25-100 mL

4590 HP Micro, Bench Top, 25-100 mL

<u>4590 HP/HT Micro, Bench Top,</u> <u>25-100 mL</u>

5100 Low Pressure Glass or Metal, Bench Top, 160-1500 mL

5500 Compact Bench Top, <u>25-600 mL</u>

Reactor Selection Guide

It is possible to convert most of these reactors from one size to another within the same series. This is done by substituting a longer or shorter cylinder with corresponding internal fittings including the stirrer shaft, thermowell or thermocouple, dip tube and cooling coil (if installed). In some cases, the heater will also need to be changed. The Parr **Customer Service** department will be happy to provide a list of the appropriate conversion parts for any contemplated conversion.

The selection process starts with establishing the Four Basic Specifications discussed below. Having set these requirements, the user can then identify a suitable series group from the Guide to

Establish 4 Basic Specifications

1. Maximum Operating Pressure

Parr offers a number of operating pressures. Vessels 600 mL and under in volume are typically rated for 3000 psi, and vessels 1 L or larger are typically rated for 1900 psi. We also offer 5000 psi versions of models up to approximately 2 L, and 3000 psi versions of models up to 5.5 L. Some pressure limits can be increased with custom designs. However, higher pressure vessels generally require thicker walls, which can make temperature control more difficult, and larger volume, high pressure vessels are heavier and more difficult to handle. We also have a line of 5100 Series Low Pressure Reactors which can handle pressures up to 150 psi with a glass cylinder or up to 1000 psi with a metal cylinder.

2. Maximum Operating Temperature

Traditionally, the choices here have been up to 350 °C for vessels with PTFE gaskets and up to 500 °C for flexible graphite (previously metal) gaskets. Parr added the option of a self-sealing O-ring closure for general purpose vessels over a full range of volumes. These quick closing designs are limited to 225 °C, although this can be raised to 300 °C with special O-ring materials.

Moveable and Fixed Head Designs



The moveable head is best when the user prefers to remove the entire reactor in one piece after running the operation.



Parr Stirred Reactors. The user should then review the standard fittings. Finally, a list of Secondary Specifications should determine some of the finer details of the system.

3. Vessel Size

Standard Parr stirred reactors are offered in many sizes ranging from 25 mL to 18.75 liters (5-gallon). It should be noted that these volumes refer to the free space in the vessel, and for safe operation the maximum liquid charge held in the vessel typically should not exceed two-thirds of the available free space in sealed batch operations. Generally, several vessel volumes are offered within most series, and reactor sizes can be reconfigured with conversion parts.

4. Material of Construction

Parr reactors are normally made of Type 316 Stainless Steel, but they can be made of other alloys as well.

The list of available materials of construction includes:

- Type 316/316L Stainless Steel
- Alloy A-286
- Alloy 20
- Alloy 230
- Alloy 400
- Alloy 600
- Alloy 625
- Alloy B-2/B-3
- Alloy C-276
- Nickel 200
- Titanium Grades 2, 3, 4, & 7
- Zirconium 702 & 705

Other materials may be available upon request.

The majority of organic reactions can be handled in a standard T316 Stainless Steel vessel, but other corrosion resistant alloys are available to provide vessels suitable for use with a wide range of corrosive acids, bases, salts and gases. Special alloy construction can be provided for both the internal parts of the vessel and the external valves and fittings. However, there are considerable cost savings if the user can accept standard external parts made of stainless steel instead of a special alloy.

There is more detailed information on special materials in the "Materials of Construction" section of Chapter 1, starting on <u>page 10</u>, of this catalog.

2 Sele	ect Appropriate Series			
Series No.	Reactor Type	Nominal Size	Maximum Pressure psi (bar)	Maximum Temperature °C
<u>4520</u>	Bench Top	1000 and 2000 mL	1900 (131)	225-350
<u>4520 HP</u>	High Pressure Bench Top	970 and 1900 mL	2900 (200)	350
<u>4530</u>	Moveable Cart or Floor Stand	1000 and 2000 mL	1900 (131)	225-350
<u>4530 HP</u>	High Pressure Moveable Cart or Floor Stand	970 and 1900 mL	2900 (200)	350
<u>4540</u>	High Pressure Bench Top, Floor Stand or Cart	600 and 1200 mL	5000 (345)	350
<u>4550</u>	Moveable Cart or Floor Stand	1 and 2 gallon (3.75 and 7.99 L)	1900 (131)	225-350
<u>4555</u>	Floor Stand	5 and 2.6 gallon (18.75 and 10 L)	1900 (131)	225-350
<u>4560</u>	Mini, Bench Top	100-600 mL	3000 (207)*	225-350
<u>4560 HT</u>	High Temperature Mini, Bench Top	300-600 mL	2000 (138)	500
<u>4571-4572</u> 4577-4578	High Pressure/High Temperature, Moveable Cart or Floor Stand	1000 and 1800 mL	5000 (345)	500
<u>4575A-4576A</u> <u>4575B-4576B</u>	High Pressure/High Temperature, Bench Top	250 and 500 mL	5000 (345)	500
<u>4581-4584</u>	High Pressure/High Temperature, Moveable Cart or Floor Stand	1 and 1.5 gallon (3.75 L and 5.5 L)	3000 (207)*	500
<u>4590</u>	Micro, Bench Top	25-100 mL	3000 (207)*	225-350
<u>4590 HP</u>	High Pressure Micro, Bench Top	25-100 mL	5000 (345)	350
<u>4590 HP/HT</u>	High Pressure Micro, Bench Top, High Temperature (Fixed Head Only)	25-100 mL	5000 (345)	500
<u>5100</u>	Low Pressure Glass or Metal Reactors	160 mL-1.5 L	150 (10.3) Glass 1000 (69) Metal	225 225-300
<u>5500</u>	Mini or Micro, Bench Top, Compact Stand	25-600 mL	3000 (207)*	225-350
				0 bar for CE ord

Custom options available. Contact Parr's Customer Service Department.

Size	Pressures (psi)	Temperatures
25 mL 18.75 L	1900 3000 5000	225 °C 500 °C
Parr Instrument Company offers	Maximum Allowable Working Pressures (MAWP) are	Operating Temperatures have

Parr Instrument Company offers laboratory reactors and pressure vessels in sizes from 25 mL to 18.75 L. Generally it is best to select a size that will allow for 1/3 free space. This allows for some liquid expansion during the heating phase of a reaction. Maximum Allowable Working Pressures (MAWP) are determined by adherence to the ASME pressure vessel design criteria. There are three ranges; up to 1900 psi, up to 5000 psi, and in most cases reactors in the 25 mL to 600 mL range can be rated to 3000 psi. Parr Instrument Company recommends research be restricted to 90% of any vessel's MAWP rating because of the standard rupture disc installed. A vessel's pressure capacity is directly tied to the temperature the vessel will be operating.

Operating Temperatures have traditionally been limited to 350 °C using a flat PTFE gasket. Flexible graphite gaskets have now extended the range to 500 °C. New quick opening designs have been added to the Parr

designs have been added to the Parr product line that use O-rings for the seal and normally limit the temperature to 225 °C (Viton®) or 300 °C (Kalrez®).



Reactor Selection Guide

Confirm Standard Fittings

Pressure Gage, analog type, which shows the pressure within the vessel at all times.

Liquid Sampling Valve for withdrawing liquid samples through the dip tube shared with the gas inlet valve. Incoming gas can be used to clear the dip tube between liquid samples.

Gas Release Valve to release gas from the reactor during or at the completion of a run.

Thermocouple or Thermowell for measuring the temperature within the vessel. In small reactors, a thermocouple encased in a metal sheath extends directly into the vessel. In larger reactors, and in vessels constructed of special alloys, the control thermocouple is inserted into a thermowell which extends to a point near the bottom of the vessel.

Internal Stirring System consists of a motor drive magnetically coupled to an internal stirrer shaft with attached turbine-type impeller(s).

Dip Tube

Cooling Coil

Safety Rupture Disc to protect the vessel and the operator from dangerous pressures beyond the rated limit for the vessel. **Gas Inlet Valve** for charging gas into the reactor. This valve and the liquid sampling valve are connected to a dip tube which extends to the bottom of the vessel.

Guide or Foot Bearing with a PTFE bushing to support and stabilize lower stirrer shaft. (*Not required with footless magnetic drive*)

Specify Detailed Options

There are a number of options to be considered and selections to be made in order to complete the specifications for a reactor. You will need to choose from the following:

1. Sealing Style

Parr reactors may be provided with a flat gasket seal or an O-ring seal. Most commonly, a flat PTFE gasket is used. These are generally good for temperatures up to 350 °C. A graphite seal for temperatures up to 500 °C is used for high temperature systems.

O-ring seals allow for a convenient, boltless closure, but the material must be checked carefully against the intended process as there are many solvents and gases which attack O-rings.

Our typical O-ring seal is designed to be self-sealing at high pressures. However, this design needs to be modified when the reactor is going to be operated at or below 20% of its high pressure rating.

2. Mounting Style

Parr reactors may be provided with a moveable vessel mounting or fixed head mounting.

In the moveable vessel style the complete vessel assembly (only the head in the one gallon or larger) is removed from the heater for charging, product recovery, and vessel cleaning. In the fixed head style the head remains in the mounting fixture and the reactor heater, then the cylinder drops away to cool and open the vessel. This is useful for users who wish to leave inlet and feed lines, discharge and vent lines, condensers, and similar head connections undisturbed between runs. The head is easily removable by sliding the head out of its slot on the stand.

3. Support Stand

Parr stirred reactors are most commonly mounted on a bench top stand for sizes of 2 liter or smaller. Vessels 1 liter or larger may be mounted on either moveable carts or fixed floor stands. Care should be taken to check stand dimensions with the available space at the user's facility.

4. Stirrer Type

Parr reactors have a magnetic drive equipped for most low viscosity fluids and light slurries. Drives with higher coupling torque and more powerful motors are available for higher viscosity applications.

Also, many reactors require a lower support bracket to stabilize the stirring shaft. In applications with abrasive particles or heavy mixing, a "footless" magnetic drive with a larger diameter shaft may be provided.

5. Cooling Coil

An internal cooling coil can be installed in all reactors (except the Micro sizes) to remove the heat during a reaction and/ or to cool the vessel at the end of a run. In some reactors a cooling coil is furnished as a standard fitting. In others, a coil can be added as an option in either serpentine style or alternately in a spiral style for selected models. For the micro sizes, with their limited dimensions, a cold finger may be added as an option.

6. Gage and Rupture Disc Ranges

The pressure gage range must be selected to provide the resolution desired, while having a range high enough to handle the maximum pressure the reactor will experience. One would be tempted to simply get a gage with the largest range possible, but this can reduce the resolution to unacceptable levels. One good rule of thumb for good resolution is to operate with the gage at half of its maximum range.

The rupture disc is most commonly matched to the Maximum Allowable Working Pressure of the vessel, ensuring the disc will burst before reaching a pressure beyond the MAWP range. Care must be taken not to take the vessel above 90% of its burst pressure under normal operation. This will help protect the disc from the effects of repeated pressure cycles.

7. Controller Options

Parr controllers are typically set up to control temperature, but may also be set up with digital displays of pressure, stirring speed, secondary temperature, and other displays and controls.

These may be set up with data logging via PC as well. For more information on these options, see the "<u>Reactor</u> <u>Controllers Section</u>" of this catalog.

8. Certification

ASME, CE, and PARR Certifications are available for users who require these recognized quality assurance certifications.

9. Custom Options

A wide range of custom options such as special openings in the head or cylinder, high viscosity stirrers, windows, special valves, heaters or jackets, explosion proof wiring, and volume modifications are available.

Various accessories, such as glass or PTFE liners, condensers, catalyst holders, and alternate stirrers are available to further adapt these reactors to the individual user's applications.

Detailed information for these and other options are found in Optional Accessories, <u>Chapter 7, starting on page 115</u> of this catalog.

How to use the Ordering Guide:

The last page of each Reactor Series in this chapter is an Ordering Guide. A composite identification/order number can be generated by combining the individual symbols from the separate sections. Omit any symbols not desired for the system.

Example: Using the Ordering Guide from the Series 4520 Stirred Reactors we have chosen to order a 1000 mL fixed head reactor, PTFE Flat Gasket, the vessel manufactured out of T316 Stainless Steel, a general purpose magnetic drive made of T316 Stainless Steel, 115V power supply, 1/4 hp explosion-proof motor, 3000 psi pressure gage, no bottom drain valve, and a 4848 controller.

The order number for that particular configuration would appear as follows:



<u>No. 4523-T-SS-M-SS-115-XP.25-3000-4848</u>

Α.	С.	D.	Е.	F.	G.	н.	Ι.	L.
Model	Gasket	Material	Stirrer	Drive MOC	Voltage	Motor	Gage	Controller
No. 4523	-T	-SS	-M	-SS	-115	-XP.25	-3000	-4848



Type: General Purpose

Stand: Bench Top

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: 1000-2000 970-1900 HP

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temp., °C:

225 w/ FKM O-ring 300 w/ FFKM O-ring 350 w/ PTFE Flat Gasket

High Pressure (HP) MAWP Rating, psi (bar): **2900 (200)**

Maximum Operating Temperature, °C, at High Pressure (HP): **350** @ **2900 psi**

Series 4520 Bench Top Reactors, 1 & 2 L

These are the largest of the Parr Reactors that can be handled on a bench top.

These reactors are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

An optional HP (High Pressure) flat gasket version has been added for maximum allowable working pressure of 2900 psi (200 bar) at 350 °C. Both fixed head and moveable vessel designs are available. An optional pneumatic lift is available for the heavier 2 liter cylinder and heater. It should be noted, however, that the 2 liter. fixed head model is tall and may not be convenient to operate on a standard height bench top. We recommend using a floor stand. (See Series 4530)

With their 4-inch internal diameter, these mid-size reactors have sufficient space for special modifications, such as an internal cooling coil, bottom drain valve (not on 2 liter bench top), ball valve for a solids charging port, catalyst addition devices, condensers, electrical feed-throughs and more. Details are provided in Chapter 7 of this catalog, starting on <u>page 115</u>.

The standard magnetic stirrer drive on these 4520 Models works well for reaction mixtures with viscosities up to 25,000 centipoise. For heavier stirring loads, these reactors can be equipped with larger magnetic drives, more powerful motors, and drive trains capable of delivering additional stirring torque.



Model 4523 Bench Top Reactor, 1000 mL, Fixed Head, open to show Internal Fittings, and a 4848 Controller shown with optional Expansion Modules.

4520

	in series.						
Model Number	4523	4524	4525	4526			
Approximate Volume, mL	1000	2000	1000	2000			
Approximate Volume for HP Models, mL	970	1900	970	1900			
Maximum Pressure (MAWP)	-		i (131 bar)				
HP Maximum Pressure (MAWP)		2900 ps	i (200 bar)				
Maximum Temperature							
with FKM 0-ring		225 °C					
with FFKM 0-ring			0° 00				
with PTFE Flat Gasket		35	50 °C				
Reactor Details							
Mounting Style	Fixed	Head	Move	able			
Stand Type			ch Top				
Closure		•	Gasket, no Compression Bo	0.			
High Pressure Closure	Split-Ring (12 0	•	t Gasket, no Compression B	olts for O-ring)			
Valve Connections			lale NPT				
Magnetic Stirrer, Model No.		A1120HC6 General Pur	pose or A2140HC Footless				
Maximum Torque			n-Pounds				
Impeller(s), 6 blades		2 (2.28" dia.	/ HP: 2.00" dia.)				
Stirrer Motor		1/8 hp vai	riable speed				
Pressure Gage, Size		4.5 i	nches				
Range		0-2000 psi (138 bar) /	HP: 0-3000 psi (207 bar)				
Temperature Measurement		Ther	mowell				
Cooling Coil (optional)		Serp	pentine				
Bottom Drain Valve (BDV) (optional)							
Standard Pressure (1900 psi)	3/8" NPT	N/A	3/8" NPT	N/A			
High Pressure (2900 psi)	1/4" NPT	N/A	1/4" NPT	N/A			
Lift Mechanism (optional)	Pneu	matic	N/	A			
Heater Style		Ca	alrod				
Heater Power, Watts	1000 / 1500	1500 / 2000	1000 / 1500	1500 / 2000			
Electrical Supply							
Volts, AC		115	i / 230				
Maximum Load, amps, 115 / 230		1	2/9				
Cylinder Dimensions	-						
I.D. x Depth, inches	4.0 x 5.4	4.0 x 10.5	4.0 x 5.4	4.0 x 10.5			
I.D. x Depth, inches for HP	3.75 x 5.4	3.75 x 10.5	3.75 x 5.4	3.75 x 10.5			
Vessel Assembly Weight, pounds	33	38	30	36			
Vessel Assembly Weight for HP, pounds	36	44	35	43			
Cylinder Weight, pounds	8.9	13.7	8.9	13.7			
Cylinder Weight for HP, pounds	11.3	18.4	11.3	18.4			
Reactor Dimensions							
Width x Depth w/o Controller, inches		17	x 24				
Height, inches	38 (48 with BDV)	48	38 (48 with BDV)	48			
Weight, pounds	105	120	100	115			
Spare Parts Kit		45	i29M				

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.



Series Number:

Index to Options

Option	Page
Certification	8
Materials of Construction	n 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fitting	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4520 Bench Top Reactors, 1 & 2 L



Model 4526 Bench Top Reactor, 2000 mL, Moveable Head, and a 4848 Controller shown with optional Expansion Modules.

Series 4520 Ordering Guide

The Order No. for the Base System is: 452_-T-SS-M-115-VS.12-2000-4848

A composite identification number to be used when ordering a 4520 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, <u>please see page 27</u>.

A Base Mode		
Model No.	Size	Vessel Style
4523	1000 mL	Fixed
4524	2000 mL	Fixed
4525	1000 mL	Moveable
4526	2000 mL	Moveable
	s are 970 mL & 1900 mL	movedbio
B High Press		000 nai (121 har)
-No Symbol -HP	Standard Configuration (1 2900 psi (200 bar)	900 psi / 131 bar)
	• • •	
C Gasket / Ma	aximum Temperatur	9
-0V	FKM 0-ring, 225 °C	
-0K	FFKM 0-ring, 300 °C	
-T	PTFE Flat Gasket, 350 °C	
D Vessel Mat	erial of Construction	1
-SS	T316 Stainless Steel	
-M0	Alloy 400	
-IN	Alloy 600	
-HB	Alloy B-2 / B-3	
-HC	Alloy C-276	
-CS	Alloy 20	
-TI2	Titanium Grade 2	
-TI4	Titanium Grade 4	
-ZR702	Zirconium Grade 702	
-ZR705	Zirconium Grade 705	
	or complete list of availabl	e alloys.
E Magnetic S	tirrer Drive	
-M	General Purpose, 16 in-lb	
-HD	Heavy Duty, 60 in-lb	
-FMD1	Footless, General Purpos	a 16 in-lb
-FMD2	Footless, Heavy Duty, 60 i	
\sim		
	Material of Constru	
-MOC Symbol	Indicate Material of Cons	truction
G Electrical S	upply	
-115	115 VAC	
-230	230 VAC	
H Motor Opti	on	
-VS .12	Variable Speed, 1/8 hp	
-VS .25	Variable Speed, 1/4 hp	
-VS .50*	Variable Speed, 1/2 hp (23	30V only)
-XP .25	Explosion Proof, Variable	· ·
-XP .50*	Explosion Proof, Variable	
-AM .25	Air Motor, 1/4 hp	
-AM .50	Air Motor, 1/2 hp	
* For use with Hea	vy Duty Magnetic Stirrer	
I Pressure G		
· Fressure G		

Pressure G	age
-3000	3000 psi / 207 bar
-2000	2000 psi / 138 bar
-1000	1000 psi / 69 bar
-600	600 psi / 40 bar
-200	200 psi / 14 bar
-100	100 psi / 7 bar

-BDV 1/4" NPT for HP sys -4848 PID Control, Ramp Control, and Data la cable not included) display modules. -4848B Same as above but additional display not controller	
-BDV Bottom Drain Valve 1/4" NPT for HP sys Controller -4848 (included in base system) PID Control, Ramp Control, and Data la cable not included) display modules. -4848B Same as above but additional display modules	tems (4523 & 4525 only) & Soak Programming, Motor Speed ogging with Software. (RS-485 to USB For use with up to three additional for use with up to six nodules.
-BDV Bottom Drain Valve 1/4" NPT for HP sys L Controller -4848 (included in base system) PID Control, Ramp Control, and Data la cable not included) display modules. -4848B Same as above but additional display not -A2110E Motor Controller	tems (4523 & 4525 only) & Soak Programming, Motor Speed ogging with Software. (RS-485 to USB For use with up to three additional for use with up to six nodules.
-4848 (included in base system)PID Control, Ramp Control, and Data la cable not included) display modules4848BSame as above but additional display m -A2110E	ogging with Software. (RS-485 to USB For use with up to three additional for use with up to six nodules.
-4040 (included in base system)Control, and Data la cable not included) display modules4848BSame as above but additional display in -A2110E	ogging with Software. (RS-485 to USB For use with up to three additional for use with up to six nodules.
-4046D additional display n -A2110E Motor Controller	nodules.
-4871 Process Controller	(for enhanced control options)
See <u>Chapter 6</u> for a complete list o	f controllers and options.
M 4848 Controller Option	IS
-TDM Tachometer Displa	
-MCM Motor Control Mod	•
-PDM Pressure Display N	1odule
-HTM High Temperature	Cut Off Module
-ETLM External Temperatu	ıre Limit Module
-MTM Motor Torque Mod	ule (<i>must be installed with the MCM</i>)
-SVM Solenoid Valve Mo	dule (for cooling control)
-A1925E4 RS-485 to USB Cab (required for data la	le for 4848 Controller ogging)
-A1925E6 RS-485 to USB Con	verter, isolated, 30-ft
-A2208E RS-485 Daisy Chair (Must be used with	n for Multiple Controllers 1 1925E6)
-A3504HC SpecView Software	e Package for 4838/4848 Controller
N Custom Options (List A	All Desired)
-AS Anchor Stirrer	
-PS Paddle Stirrer	
-SA Spiral Stirrer	
-GE Gas Entrainment St	tirrer
-BF Removable Baffle S	Set
-SB Static Catalyst Bas	ket
-DB Dynamic Catalyst E	Basket
-CAD Internal Catalyst Ac	
-XCAD External Catalyst A	
-SCP Solids Charging Po	rt (Ball Valve)
-RC Reflux Condenser	
-RTC Reflux/Take-Off Co	ndenser
-WJ Welded Jacket	
	eater with Cooling Channels
-PL Pneumatic Lift	6
See <u>Chapter 7</u> for a complete list o	i upuunai accessories.
0 Certifications	
-ASME ASME Documentat	ion
-CE CE Documentation	
-P Parr Certification	

Please note that all options and combinations are not compatible with all models.

Spare Parts Kit for 4520 Series

P Spare Parts Kit

-4529M





Type: General Purpose

Stand: Floor Stand or Cart

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: 1000-2000 Standard 970-1900 HP

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temp., °C:

225 w/ FKM O-ring 300 w/ FFKM O-ring 350 w/ PTFE Flat Gasket

High Pressure (HP) MAWP Rating, psi (bar): **2900 (200)**

Maximum Operating Temperature, °C, at High Pressure (HP): **350** @ **2900** psi

Series 4530 Floor Stand Reactors, 1 & 2 L

Series 4530 1 & 2 L Floor Stand Reactors can be moved, modified to fit your needs, and configured for polymer studies. They are also available in high pressure versions.

Reactors configured for high torque applications. The reactors in this series are mounted on a sturdy floor stand which will accommodate larger stirrer drive motors and stronger power trains that are generally not available for the bench top models. With these options it is possible to provide the high torque and low stirring speeds required for work with polymers, slurries, and other mixtures with viscosity of 1 million centipoise and more. These mountings are also better suited for the installation of a drain valve in the cylinder bottom for convenient removal of the reaction mixture while it is still hot.

Reactors requiring extensive modifications. Floor stand mountings provide a good base for reactor modifications and for the addition of accessories, such as: condensers, packed columns, special motors, custom heaters, jacketed vessels, automatic valves or regulators, and many other fittings. Both fixed head and moveable vessel designs are available. An optional pneumatic lift is available for fixed head models

Reactors that need to be moved. These reactors are designed for use in an area where a bench top is not available. Users who wish to move the reactor to storage when not in use, will appreciate the moveable cart design of the Series 4531 and 4532.

These reactors are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

- **HP versions.** These reactors are also available in a HP (High Pressure) version up to 2900 psi (200 bar) at 350 °C.

4534 Floor Stand Reactor, 2000 mL, Fixed Head, with Heavy-Duty Footless Mag Drive, Serpentine Cooling Coil, optional Solids Charging Device, Cylinder equipped with Welded Jacket, and Bottom Drain Valve.



4530

Shaded bar indicates specifications that change w	vithin series.						
Model Number	4531	4532	4533	4534	4535	4536	
Approximate Volume, mL	1000	2000	1000	2000	1000	2000	
Approximate Volume for HP Models, mL	970	1900	970	1900	970	1900	
Maximum Pressure (MAWP)	1900 psi (131 bar)						
Maximum Pressure for HP (MAWP)	2900 psi (200 bar)						
Maximum Temperature							
with FKM O-ring	225 °C						
with FFKM 0-ring	300 °C						
with PTFE Flat Gasket			35	0 °C			
Reactor Details							
Mounting Style	Move	eable	Fixed	Head	Mov	eable	
Stand Type	Movea	ble Cart	Floor	Stand	Floor	Stand	
Closure	Split	t-Ring (6 Compres	sion Bolts for Flat	Gasket, No Compi	ession Bolts for O	-ring)	
High Pressure Closure	Split	-Ring (12 Compres	sion Bolts for Flat	Gasket, No Comp	ression Bolts for O	-ring)	
Valve Connections				ale NPT			
Magnetic Stirrer, Model No.		A112)HC General Purp	ose or A2140HC Fo	ootless		
Maximum Torque				-Pounds			
Impeller(s), 6-Blade				HP: 2.00" dia.)			
Stirrer Motor				iable Speed			
Pressure Gage, Size			•	nches			
Range		Standard		ar) / HP: 0-3000 ps	si (207 har)		
Temperature Measurement		Otandure	-	nowell	(207 Bul)		
Cooling Coil (optional)				entine			
Bottom Drain Valve (BDV) (optional)				enune			
Standard Pressure (1900 psi)			3/8	NPT			
High Pressure (2900 psi)				NPT			
Lift Mechanism (optional)	N	/A		atic Lift	N	/A	
Heater Style	N	/A			I	/A	
-	10	100	1		1000 / 1500	1500 / 2000	
Heater Power, Watts	13	00	1000 / 1500	1500 / 2000	1000 / 1500	1500 / 200	
Electrical Supply				/ 000			
Volts, AC			-	/ 230			
Maximum Load, amps, 115 / 230			12	2/9			
Cylinder Dimensions		4.0.405	40.54	4.0.405	10.54	4.0 40.5	
I.D. x Depth, inches	4.0 x 5.4	4.0 x 10.5	4.0 x 5.4	4.0 x 10.5	4.0 x 5.4	4.0 x 10.5	
I.D. x Depth, inches for HP	3.75 x 5.4	3.75 x 10.5	3.75 x 5.4	3.75 x 10.5	3.75 x 5.4	3.75 x 10.5	
Vessel Assembly Weight, pounds	30	35	33	38	301	36	
Vessel Assembly Weight for HP, pounds	35	43	36	44	35	43	
Cylinder Weight, pounds	8.9	13.7	8.9	13.7	8.9	13.7	
Cylinder Weight for HP, pounds	11.3	18.4	11.3	18.4	11.3	18.4	
Reactor/Stand Dimensions							
Width x Depth w/o Controller, inches	35 x	x 18	21	x 28	21	x 28	
Height, inches	5	0	(3	E	3	
Weight, pounds	200	215	225	240	225	240	

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.



Series Number:

Index to Options

Option	Page
Certification	8
Materials of Construction	1 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4530 Floor Stand Reactors, 1 & 2 L



4532 Moveable Head on Moveable Cart, 2000 mL, and a 4848 Temperature Controller with optional Expansion Modules.



2000 mL Moveable Vessel



Series 4530 Ordering Guide

The Order No. for the Base System is: 453_-T-SS-M-115-VS.25-2000-4848

A composite identification number to be used when ordering a 4530 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Mo	odel				K Intern	al Cooling Coil		
Model No.	o. Size Vessel Style		-SC Serpentine Coil					
4531	1000 mL	Moveable	Head, Mo	veable Cart	Bottom Drain Valve			
4532	2000 mL	Moveable	Head, Mo	veable Cart		Bottom Drain Valve, 3/8" NPT for standard systems,		
4533	1000 mL	Fixed Hea	d, Stationa	ary Floor Stand	-BDV	1/4" NPT for HP systems (1000 / 970 mL only)		
4534	2000 mL	Fixed Hea	d, Stationa	ary Floor Stand				
4535	1000 mL	Moveable	Head, Sta	tionary Floor Stand	M Contro	bller		
4536	2000 mL	Moveable	Head, Sta	tionary Floor Stand	-4848	PID Control, Ramp & Soak Programming, Motor Speed		
High Pressure	sizes are 970 m	nL and 1900 ml	_ respectiv	vely.	(included in	Control, and Data logging with Software. (RS-485 to USI cable not included) For use with up to three additional		
B High Pr	essure Opti	on			base system	display modules.		
-No Symbol			00 nsi / 13	1 har)	•	1 7		
-No Symbol Standard Configuration (1900 psi / 131 bar) -HP 2900 psi (200 bar)		-4848B	Same as above but for use with up to six additional display modules.					
\sim					Δ.2110Γ	1, ,		
C Gasket / Maximum Temperature					-A2110E	Motor Controller Process Controller (for enhanced control options)		
-0V	FKM 0-ring,				-4871			
-0K	FFKM 0-ring	-			See <u>Chapter 6</u> for a complete list of controllers and options.			
-T	PTFE Flat Ga	sket, 350 °C			_ <u>N</u> 4848 (Controller Options		
D Materia	l of Constru	uction			-TDM	Tachometer Display Module		
-SS	T316 Stainle	ss Steel			-MCM	Motor Control Module w/Tachometer		
-M0	Alloy 400				-PDM	Pressure Display Module		
-IN	Alloy 600				-HTM	High Temperature Cut Off Module		
-HB	Alloy B-2 / B	-3			-ETLM	External Temperature Limit Module		
-HC	Alloy C-276				-MTM*	Motor Torque Module		
-CS	Alloy 20				-SVM	Solenoid Valve Module (for cooling control)		
-TI2	Titanium Gra	de 2			-A1925E4	RS-485 to USB Cable for 4848 Controller		
-TI4	Titanium Gra	de 4			-ATJZJL4	(required for data logging)		
-ZR702	Zirconium G	rade 702			-A1925E6	RS-485 to USB Converter, isolated, 30-ft		
-ZR705 Zirconium Grade 705		-A2208E	RS-485 Daisy Chain for Multiple Controller					
See <u>page 10</u> oi	r <u>24</u> for comple	te list of availa	ble alloys		AZZOOL	(must be used with A1925E6)		
E Magnet	ic Stirrer D	rivo			-A3504HC	SpecView Software Package for 4838/4848		
-M		pose, 16 in-lb			* The MTM	must be installed in conjunction with the MCM.		
-HD	Heavy Duty, 60 in-lb		O Custo	m Options (List All Desired)				
-XHD		Duty, 120 in-lb	1		-AS	Anchor Stirrer		
-FMD1	Footless, General Purpose, 16 in-lb		-PS	Paddle Stirrer				
-FMD2	Footless, Heavy Duty, 60 in-lb		-SA	Spiral Stirrer				
-FMD3		tra Heavy Duty			-GE	Gas Entrainment Stirrer		
\sim					-BF	Removable Baffle Set		
	rive Materia				-SB	Static Catalyst Basket		
-MOC Symbol	Indicate Ma	terial of Constr	ruction		-DB	Dynamic Catalyst Basket		
G Electric	al Supply				-CAD	Internal Catalyst Addition Device		
-115	115 VAC				-XCAD	External Catalyst Addition Device		
-230	230 VAC				-SCP	Solids Charging Port (Ball Valve)		
H Motor C	Intion				-RC	Reflux Condenser		
	Variable Spe	od 1/9 hr			-RTC	Reflux/Take-Off Condenser		
-VS .12					-WJ	Welded Jacket		
-VS .25	Variable Spe				-ABH	Aluminum Block Heater with Cooling Channels		
-VS .50		ed, 1/2 hp (230		hn	-PL	Pneumatic Lift		
-XP .25		oof, Variable S	-			7 for a complete list of optional accessories.		
-XP .50		oof, Var. Spee	u, 1/2 np (
-AM .25	Air Motor, 1/	•				cations		
-AM .50	Air Motor, 1/	•		4	-ASME	ASME Documentation		
-GDD	Geared Dire	ct Drive (Fixed		,	-CE	CE Documentation		
I Geared	Direct Driv	e Ratio*	J Pre	ssure Gage	-P	Parr Certification		
-3D 3:1, (0-600 rpm		-3000	3000 psi / 207 bar	0 Spare	Parts Kit		
-5D 5:1, 0	0-360 rpm		-2000	2000 psi / 138 bar	-4539M	Spare Parts Kit for 4530 Series		
-10D 10:1,	, 0-180 rpm		-1000	1000 psi / 69 bar				
	otor Option Only	/.	-600	600 psi / 40 bar	Please note that all options and combinations are not compatible with all models.			
	. ,		-200	200 psi / 14 bar				



Type: High Pressure

Stand: Bench Top, Floor Stand, or Cart

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: 600 and 1200

High Pressure MAWP Rating, psi (bar): **5000 (345)**

Standard Maximum Operating Temperature, °C: **350 w/ PTFE Flat Gasket** Series 4540 High Pressure Reactors, 600 & 1200 mL

This series of stirred reactors has been designed for users who need higher operating pressures than the 1900 to 2900 psi offered by the General Purpose Reactors, but do not require the high operating temperatures provided by the Series 4570 High Pressure, High Temperature Reactors. The 4540 series reactors offer working pressures to 5000 psi (345 bar) at temperatures to 350 °C.

These vessels have been designed with outside dimensions comparable to the Series 4520 and 4530 Reactors so that they can use the same support system, stirrer drive and heater as these popular general purpose models. This provides not only an attractively priced high pressure/ moderate temperature system, but also reactors that can be interchanged with the 1 and 2 liter sizes. The thicker walls required for higher operating pressures reduce the volumes of these reactors to 600 and 1200 mL.

These reactors can be used in either the bench top or floor stand mountings. While the 1200 mL reactor is offered as a fixed head bench top model; it is too tall and too heavy to be handled comfortably on a standard height bench top. It is recommended that the floor stand support option should be selected unless the user has an adjustable bench top which will accommodate the overall height of the 1200 mL systems.

Model 4544 High Pressure Bench Top Reactor, 600 mL, Moveable Head, with heater lowered, and a 4848 Controller shown with optional Expansion Modules.



4540

2

Shaded bar indicates specifications that cha	ange within series.						
Model No Moveable Vessel	4544	—	4544A	_	4544C	4546C	
Model No Fixed Head	4545	4547	4545A	4547A			
Approximate Volume, mL	600	1200	600	1200	600	1200	
Maximum Pressure (MAWP)			5000 psi	(345 bar)			
Maximum Temperature							
with PTFE Flat Gasket	350 °C						
Reactor Details							
Mounting Style	Moveable or Fixed Head						
Stand Type	Bend	ch Top	Floor	Stand	Moveable Cart		
Closure	Sp	lit-Ring (8 Compres	sion Bolts for Flat (Gasket, no Compre	ssion Bolts for O-r	ing)	
Valve Connections			1/4" Ma	le NPT			
Magnetic Stirrer, Model No.		A1120	HC6 General Purpo	ose or A2140HC Fo	otless		
Maximum Torque	16 Inch-Pounds						
Impeller(s), 6-Blade	2 (1.75" dia.)						
Stirrer Motor	1/8 hp 1/8 hp 1/4 hp				hp		
Pressure Gage, Size		4.5 inches					
Range		0-5000 psi (345 bar)					
Temperature Measurement			Therm	owell			
Cooling Coil (optional)			Serpentin	e or Loop			
Bottom Drain Valve (BDV) (optional)	Ν	I/A	1" NPS		N/A		
Lift Mechanism	N/A	Pneumatic	N/A	Pneumatic	N/A	N/A	
Heater Style			Calrod				
Heater Power, Watts	1000 / 1500	1500 / 2000	1000 / 1500	1500 / 2000	1300	1300	
Electrical Supply							
Volts, AC		115 / 230					
Maximum Load, amps, 115 / 230			12	/ 9			
Cylinder Dimensions							
I.D. x Depth, inches	3.25 x 4.7	3.25 x 9.8	3.25 x 4.7	3.25 x 9.8	3.25 x 4.7	3.25 x 9.8	
Vessel Assembly Weight, pounds	4544: 50 4545: 53	65	4544A: 50 4545A: 53	65	50	61	
Cylinder Weight, pounds	18	29	18	29	18	29	
Reactor Dimensions							
Width x Depth w/o Controller, inches	17	x 24	21 >	: 28	35 :	k 18	
Height, inches	38	38 48 63 50				0	
Cylinder Weight, pounds	125	135	175	185	175	185	
Spare Parts Kit			454	9M			

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.



Index to Options

•	
Option	Page
Certification	8
Materials of Construction	n 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4540 High Pressure Reactors, 600 & 1200 mL



Model 4546C Moveable Head Reactor on a Moveable Cart, 1200 mL, and a 4848 Reactor Controller shown with optional Expansion Modules.

Series 4540 Ordering Guide

The Order No. for the Base System is: 454_-T-SS-M-115-VS.12-2000-4848

A composite identification number to be used when ordering a 4540 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, <u>please see page 27</u>.

A Base Model				
Model No.	Size	Vessel Style		
4544	600 mL	Moveable Head, Bench Top		
4544A	600 mL	Moveable Head, Floor Stand		
4544C	600 mL	Moveable Head, Moveable Cart		
4545	600 mL	Fixed Head, Bench Top		
4545A	600 mL	Fixed Head, Stationary Floor Stand		
4546C	1200 mL	Moveable Head, Moveable Cart		
4547	1200 mL	Fixed Head, Bench Top		
4547A	1200 mL	Fixed Head, Stationary Floor Stand		

B Gasket / Maximum Temperatur

	B Gaske	et / Maximum Temperature
	-T	PTFE Flat Gasket, 350 °C
ļ		
	C Mater	ials of Construction
	-SS	T316 Stainless Steel
	-M0	Alloy 400
	-IN	Alloy 600
	-HB	Alloy B-2 / B-3
	-HC	Alloy C-276
	-CS	Alloy 20
	-TI2	Titanium Grade 2*
	-TI4	Titanium Grade 4*
	-NI	Nickel 200*
	-ZR702	Zirconium Grade 702*
	-ZR705	Zirconium Grade 705*

* Maximum Pressure and/or Temperature are limited.

See page 10 or 24 for complete list of available alloys.

000 <u>page 10</u>	
D Magne	etic Stirrer Drive
-M	General Purpose, 16 in-lb
-HD	Heavy Duty, 60 in-lb
-XHD	Extra Heavy Duty, 120 in-Ib
-FMD1	Footless, General Purpose, 16 in-lb
-FMD2	Footless, Heavy Duty, 60 in-lb
-FMD3	Footless, Extra Heavy Duty, 120 in-lb
E Mag.	Drive Meterial of Construction
	Drive Material of Construction
-MOC Symbol	Indicates Material of Construction
e y iniser	
F Electr	ical Supply
-115	115 VAC
-230	230 VAC
G Motor	
	Option
-VS .12	Variable Speed, 1/8 hp
-VS .25	Variable Speed, 1/4 hp
-XP .25	Explosion Proof, Variable Speed, 1/4 hp
-AM .25	Air Motor, 1/4 hp
H Press	ure Gage
-10000	10000 psi / 690 bar
-7500	7500 psi / 517 bar
-5000	5000 psi / 345 bar
-3000	3000 psi / 207 bar

-SC	Serpentine Coil
-CL	Single Cooling Loop
-BDV	Bottom Drain Valve, 1" NPS
K Control	ler
-4848 (included in base system)	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software. (RS-485 to USB cable not included) For use with up to three additional display modules.
-4848B	Same as above but for use with up to six additional display modules.
-A2110E	Motor Controller
-4871	Process Controller (for enhanced control options)
See <u>Chapter 6</u>	for a complete list of controllers and options.
L 4848 C	ontroller Options
-TDM	Tachometer Display Module
-MCM	Motor Control Module w/Tachometer
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM*	Motor Torque Module
	Motor Torque Module Solenoid Valve Module (for cooling control)
-MTM*	•
-MTM* -SVM	Solenoid Valve Module (for cooling control) RS-485 to USB Cable for 4848 Controller

* The MTM must be installed in conjunction with the MCM.

M Custom	o Options (List All Desired)		
-AS	Anchor Stirrer		
-PS	Paddle Stirrer		
-GE	Gas Entrainment Stirrer		
-SP	Spiral Stirrer		
-BF	Removeable Baffle Set		
-CAD	Internal Catalyst Addition Device		
-XCAD	External Catalyst Addition Device		
-SCP	Solids Charging Port (Ball Valve)		
-RC	Reflux Condenser		
-RTC	Reflux/Take-Off Condenser		
-WJ	Welded Jacket		
-ABH	Aluminum Block Heater w/ Cooling Channels		
See <u>Chapter 7</u> for a complete list of optional accessories.			

N Certifi	estions
\sim	
-ASME	ASME Documentation
-CE	CE Documentation
-P	Parr Certification
O Spare	Parts Kit
-4549M	Spare Parts Kit for 4540 Series

Please note that all options and combinations are not compatible with all models.



Type: General Purpose

Stand: Floor Stand or Cart

Mounting Style: Moveable or Fixed Head

Vessel Sizes, Gallons (Liters): 1 (3.75) and 2 (7.99)

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

Series 4550 Floor Stand Reactors, 1 & 2 Gallon

These 4550 Reactors extend the size range of the 4530 Series to 1 and 2 gallon (3.75 and 7.99 liter) sizes, providing excellent options for pilot plant studies.

These reactors are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

Both fixed head and moveable head designs are available. The moveable head. Model 4551 and 4552, cart reactors are designed so that they can be opened or closed conveniently without removing the cylinder from the heater and without auxiliary handling equipment. The split-ring sections can be moved into place from the sides, and the compression bolts can be tightened with the vessel in place in its heater. The fixed head versions. Model 4553 and 4554, allow the head (and the connections made to it) to remain in place while the cylinder and heater are lowered with the aid of the included pneumatic lift system.

The 1 gallon size is usually recommended for high viscosity polymer studies. An optional bottom drain valve may be added for convenient product recovery. As with the smaller floor stand models, these larger, self- contained systems can be equipped with a variety of attachments, such as condensers, solids charging ports, bottom drain valves, special motors, custom heaters, jacketed vessels and automatic valves and regulators. Because of the higher wattage heaters for these reactors, all models in the 4550 Series require a 230 volt power supply.



Model 4554 Floor Stand Reactor, Two Gallon, Fixed Head, Pneumatic Lift, Hinged Split-Rings, opened to show Internal Fittings and Serpentine Cooling Coil, with 4848 Reactor Controller shown with optional Expansion Modules.

4550

2

Shaded bar indicates specifications that change with	iin series.					
Model Number	4551	4552	4553	4554		
Approximate Volume, Gallon (Liter)	1 (3.75)	2 (7.99)	1 (3.75)	2 (7.99)		
Maximum Pressure (MAWP)		1900 psi	(131 bar)			
Maximum Temperature						
with FKM 0-ring		225 °C				
with FFKM 0-ring		300 °C				
with PTFE Flat Gasket		350 °C				
Reactor Details						
Mounting Style	Mo	veable	Fixed	Head		
Stand type	Move	able Cart	Floor	Stand		
Closure	Split-Ring (10	Compression Bolts for Flat	Gasket, No Compression I	Bolts for O-ring)		
Valve Connections		1/4" NPT Male				
Magnetic Stirrer, Model No.	A1180HC General Purpose or A2160HC Footless					
Maximum Torque		60 Inch-Pounds				
Impeller(s), 6-Blades		2 (3.85" dia.)				
Stirrer Motor	1/2 hp variable speed					
Pressure Gage, Size	4.5 inches					
Range	0-2000 psi (138 bar)					
Temperature Measurement	Thermowell					
Cooling Coil (optional)		Serpe	entine			
Bottom Drain Valve (BDV) (optional)		3/8"	NPT			
Lift Mechanism (optional)		N/A	Pneu	matic		
Heater Style		Cal	rod			
Heater Power, Watts	2250	2700	2250	2700		
Electrical Supply						
Volts, AC		23	30			
Maximum Load, amps, 230	13	15	13	15		
Cylinder Dimensions						
I.D. x Depth, inches	6.0 x 8.6	6.0 x 17.2	6.0 x 8.6	6.0 x 17.2		
Vessel Assembly Weight, pounds	96	116	111	131		
Cylinder Weight, pounds	33	52	33	52		
Reactor/Stand Dimensions						
Width x Depth w/o Controller, inches	3!	5 x 18	20 ×	: 31		
Height, inches		50	7	5		
Weight, pounds	265	295	315	345		
Spare Parts Kit	4559M					

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.

Index to Options

Page
8
n 10
14
18
118
124
125
129
131
102
109
114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4550 Floor Stand Reactors, 1 & 2 Gallon



The innovative Parr Hinged Split-Rings on the 4553 and 4554 add to a safe vessel removal routine. Simply loosen the compression bolts, unlatch the split-ring closures, and pivot the split-rings out of the way.



4553 Floor Stand Reactor, 1 Gallon, Fixed Head, with Gear Drive, Hinged Split-Rings, and heater engaged.

Series 4550 Ordering Guide

The Order No. for the Base System is: 455_-T-SS-HD-230-VS.50-2000-4848

A composite identification number to be used when ordering a 4550 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, <u>please see page 27</u>.

A Base I		Vessel Style		
Nodel No.	Size	Vessel Style		
1551	1 Gallon	Moveable Head, Moveable Cart		
1552	2 Gallon	Moveable Head, Moveable Cart		
1553	1 Gallon	Fixed Head, Stationary Floor Stand		
1554	2 Gallon	Fixed Head, Stationary Floor Stand		
		um Temperature		
-0V	FKM 0-ring			
-0K	FFKM 0-rir			
-T	PTFE Flat G	asket, 350 °C		
C Mater	ials of Co	nstruction		
-SS	T316 Stainl	ess Steel		
-M0	Alloy 400			
-IN	Alloy 600			
HB	Alloy B-2 /	B-3		
HC	Alloy C-276			
CS	Alloy 20			
TI2	Titanium G			
TI4	Titanium G	rade 4		
ZR702		Zirconium Grade 702		
ZR705	Zirconium	Grade 705		
See <u>page 10</u>	or <u>24</u> for cor	nplete list of available alloys.		
D Magne	etic Stirre	r Drive		
HD	Heavy Duty			
-XHD		y Duty, 120 in-Ib		
-FMD2		eavy Duty, 60 in-lb		
FMD3	Footless, E	xtra Heavy Duty, 120 in-Ib		
E Mag. I	Drive Mat	erial of Construction		
MOC Symbol	Indicate M	aterial of Construction		
F Electr	ical Suppl	У		
-230	230 VAC			
G Motor	Option			
-VS .50		beed, 1/2 hp		
-XP .50	Explosion F	Proof Variable Speed, 1/2 hp		
AM .50	Air Motor,	1/2 hp		
-VS .75		Variable Speed, 3/4 hp		
-GDD	Geared Dir	ect Drive (Fixed Head Units Only)		
H Geare	d Direct D	Prive Ratio*		
-3D	3:1, 0-600 r			
-5D	5:1, 0-360 r	om		
-10D	10:1, 0-180			
		Motor Option Only.		
D Pressi	ure Gage			
-2000	2000 psi / 1	38 bar		
1000	1000	0.1		

-			
-2000	2000 psi / 138 bar		
-1000	1000 psi / 69 bar		
-600	600 psi / 40 bar		
-200	200 psi / 14 bar		
-100	100 psi / 7 bar		

J Interna	I Cooling Coil
-SC	Serpentine Coil
V Rottom	Drain Valve
-BDV	Bottom Drain Valve 3/8" NPT
-00 v	
L Contro	ller
-4848 (included in base system)	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software. (RS-485 to USB cable not included) For use with up to three additional display modules.
-4848B	Same as above but for use with up to six additional display modules.
-A2110E	Motor Controller
-4871	Process Controller (for enhanced control options)
See <u>Chapter (</u>	for a complete list of controllers and options.
M 4848 C	ontroller Options
-TDM	Tachometer Display Module
-MCM	Motor Control Module w/Tachometer
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM*	Motor Torque Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A3504HC	SpecView Software Package for 4838/4848
* The MTM m	ust be installed in conjunction with the MCM.
	Options (List All Desired)
-AS	Anchor Stirrer
-PS	Paddle Stirrer
-SA	Spiral Stirrer
-GE	Gas Entrainment Stirrer
-BF	Removeable Baffle Set
-CAD	Internal Catalyst Addition Device
-XCAD	External Catalyst Addition Device
CUD	Solids Charging Port (Ball Valve)
-SCP	
-RC	Reflux Condenser
-RC -RTC	Reflux/Take-Off Condenser
-RC	

O Certific	ations

-ASIVIE	ASIVIE Documentation
-CE	CE Documentation
-P	Parr Certification
_	

P Spare Parts Kit

-4559M Spare Parts Kit for 4550 Series

Please note that all options and combinations are not compatible with all models.



Type: General Purpose

Stand: Floor Stand

Mounting Style: Moveable or Fixed Head

Vessel Sizes, Gallons (Liters): **2.6 (10) and 5 (18.75)**

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

Series 4555 Floor Stand Reactors, 2.6 & 5 Gallon

The Model 4555 and 4557 Reactors with their 5 gallon (18.75 L) capacity are the largest stirred reactors

offered by Parr. The 4556 and 4558 Models are similar units with a 10 liter (2.6 gal) volume which falls between the larger 5 gallon design and the smaller 1 and 2 gallon models in the 4550 Series. Vessel styles are offered in a moveable head or fixed head design.

These reactors are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

In the moveable head design (see page 21), the vessel is held in a support system which minimizes the physical effort required to handle these heavy components. The hoist is attached to a support column which provides a convenient means for lifting the head and cylinder out of the stand. These components may be transferred to the holding position on the right side of the stand. Vessels equipped with a bottom drain valve will probably remain in the heater most of the time, but can be lifted out when necessary.

The fixed head support stand features hinged split-rings that swing to either side allowing the head to remain fixed to the stand while a pneumatic lift allows the cylinder to be raised and lowered. When lowered, the cylinder can be slid forward for cleaning and servicing.

These reactors are generally used for pilot plant or for custom chemical production purposes, usually with a variety of attachments added to the basic units. Various heaters, larger motors, heavier stirrer drives and remote controls appropriate to the size of these reactors have been designed and are available.

Modified versions of these units are available with higher working pressures and temperatures.



Hinged Split-Rings open to reveal Serpentine Cooling Coil, with Heater and Vessel lowered via Pneumatic Lift.



4557 Floor Stand Reactor, 5 Gallon, Fixed Head, 3-zone Band Heater, with Split-Rings and Pneumatic Lift.

4555

Shaded bar indicates specifications that change	within series.			
Model Number	4555	4556	4557	4558
Approximate Volume, Gallon (Liter)	5 (18.75)	2.6 (10)	5 (18.75)	2.6 (10)
Maximum Pressure (MAWP)		1900 psi	(131 bar)	
Maximum Temperature				
with FKM 0-ring		225 °C		
with FFKM 0-ring		300	٥° C	
with PTFE Flat Gasket		350	0°C	
Reactor Details				
Mounting Style	Mov	eable	Fixed I	Head
Stand Type	Table Fl	oor Stand	Floor S	Stand
Closure	Split-Ring (12 (Split-Ring (12 Compression Bolts for Flat Gasket, no Compression Bolts for O-ring)		
Valve Connections		3/8" NPT Male		
Magnetic Stirrer, Model No.		A1750HC Heavy Duty or A2160HC Footless		
Maximum Torque		60 Inch-Pounds		
Impeller(s), 6-Blades	2 (5.25" dia.)	2 (3.85" dia.)	2 (5.25" dia.)	2 (3.85" dia.)
Stirrer Motor		3/4 hp variable speed		
Pressure Gage, Size		4.5 inches		
Range		0-2000 psi (138 bar)		
Temperature Measurement		Thermowell		
Cooling Coil		Serpentine		
Bottom Drain Valve (BDV) (optional)		1.0"	NPT	
Lift Mechanism	Manu	al Hoist	Pneumatic	
Heater Style	Ceramic, 3-zone	Band Heater, 3-zone	Band Heat	er, 3-zone
Heater Power, Watts	4500	4250	5500	3300
Electrical Supply				
Volts, AC		230-1P or 400-415V 3-Phase "Y"		
Maximum Load, amps		1P-30 amps / 3P-15 amps/leg		
Cylinder Dimensions				
I.D. x Depth, inches	9.5 x 16.3	7.75 x 12.2	9.5 x 16.3	7.75 x 12.2
Vessel Assembly Weight, pounds	354	206	355	209
Cylinder Weight, pounds	157	97	157	97
Reactor Dimensions				
Width x Depth w/o Controller, inches	63	63 x 25 31 x 43		
Height, inches		91	95	87
Weight, pounds	1000	900	1000	900
Spare Parts Kit		4559	PCM	

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.

Power Requirements

Typical power requirements for Parr's large, electrically heated reactors are 40 amp single phase or 3-Phase power sources. Users are advised to have a qualified electrician determine and install an appropriate mains power supply for the large reactor system.

Large reactor systems with lower electrical power requirements, such as low temperature applications are available for use with typical 20 amp, 230 volt sources.

Contact Parr Technical Sales staff for assistance with electrical specifications.



Series 4555 Floor Stand Reactors, 2.6 & 5 Gallon

Index to Options

Option	Page
Certification	8
Materials of Construction	1 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fitting	gs 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	e 128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126



4555 Table Floor Stand Reactor, 5 Gallon Moveable Vessel with Head removed from Vessel, Manual Hoist, and a 4848M Controller.

Series 4555 Ordering Guide

The Order No. for the Base System is: 455_-T-SS-HD-230-VS.75-2000-SC-C3-4848-ASME

A composite identification number to be used when ordering a 4555 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Mo	del							
Model No.	Size	Vessel Style						
4555	5 Gallon	Moveable He	ad					
4556	10 Liter	Moveable He	ad					
4557	5 Gallon	Fixed Head						
4558	10 Liter	Fixed Head						
\sim	Maximum Tempe	rature						
-0V		FKM O-ring, 225 °C						
-0K		FFKM O-ring, 300 °C						
-T PTFE Flat Gasket, 350 °C								
C Material	s of Construction							
-SS	T316 Stainless Steel							
-M0	Alloy 400							
-IN	Alloy 600							
-HB	Alloy B-2 / B-3							
-HC	Alloy C-276							
-CS	Alloy 20							
-TI2	Titanium Grade 2							
-TI4	Titanium Grade 4							
-ZR702	Zirconium Grade 702							
-ZR705	Zirconium Grade 705							
See <u>page 10</u> or	24 for complete list of a	vailable alloys.						
	Ctimer Drive							
	C Stirrer Drive							
-HD	Heavy Duty, 60 in-lb	1						
-XHD	Extra Heavy Duty, 120							
-FMD2	Footless Magnetic Dri							
-FMD3	Footless Magnetic Dri	ve, 120 m-10						
E Mag. Dri	ve Material of Co	nstruction						
-MOC Symbol	Indicate Material of C	onstruction						
	I Cumulu							
F Electrica								
-230	230 VAC, Single Phase	9						
-400	400-415V 3-Phase "Y"							
G Motor O	ption							
-VS .75	Variable Speed, 3/4 hp)						
-XP .75	Explosion Proof Varial		пр					
-AM .100	Air Motor, 1 hp							
-GDD	Geared Direct Drive (F	ixed Head Only	()					
	Direct Drive Ratio	-	ssure Gage					
	-600 rpm	-2000	2000 psi / 138 bar					
	-360 rpm	-1000	1000 psi / 69 bar					
	D-180 rpm	-600	600 psi / 40 bar					
* for -GDD Moto	or Option Only	-200	200 psi / 14 bar					
		-100	100 psi / 7 bar					
J Internal	Cooling Coil							
-SC	Cooling Coil Serpentine Coil							
	Drain Valve							
-BDV	Bottom Drain Valve, 1	" NPT						

L Heater	Options
-C3	Ceramic, 3-zone (4555 only)
-B3	Band Heaters, 3-zone
-WJ	Welded Jacket
-FM	Flexible Mantle Heater, 1-zone, 2500 W, 200 °C Max.
	ift Mechanism (Table Floor Stand only)
-EH	Electric Hoist (115 VAC)
N Contro	ller
-4848 (included in base system)	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software. (RS-485 to USB cable not included) For use with up to three additional display modules.
-4848B	Same as above but for use with up to six additional display modules.
-A2110E	Motor Controller
-4871	Process Controller (for enhanced control options)
See <u>Chapter (</u>	for a complete list of controllers and options.
0 4040 0	
	ontroller Options
	Techemoter Dianlay Medule

4040 0	ontroller Options
-TDM	Tachometer Display Module
-MCM	Motor Control Module w/Tachometer
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM*	Motor Torque Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E	RS-485 Daisy Chain for Multiple Controller (must be used with A1925E6)
-A3504HC	SpecView Software Package for 4838/4848
* The MTM n	nust be installed in conjunction with the MCM.

P Custo	m Options (List All Desired)
-AS	Anchor Stirrer
-PS	Paddle Stirrer
-SA	Spiral Stirrer
-GE	Gas Entrainment Stirrer
-BF	Removeable Baffle Set
-XCAD	External Catalyst Addition Device
-SCP	Solids Charging Port (Ball Valve)
-RC	Reflux Condenser
-RTC	Reflux/Take-Off Condenser
-WJ	Welded Jacket
-ABH	Aluminum Block Heater
See <u>Chapter</u>	<u>7</u> for a complete list of optional accessories.

Q Certific	cations
-ASME	ASME Documentation
-CE	CE Documentation
-P	Parr Certification
	~

R Spare Parts Kit

-4559PCM Spare Parts Kit for 4555 Series

Please note that all options and combinations are not compatible with all models.



Type: Mini

Stand: Bench Top

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: 100 - 600

Standard Temperature MAWP Rating, psi (bar): 3000 (207)

High Temperature (HT) MAWP Rating, psi (bar): 2000 (138)

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

High Temperature (HT) Maximum Operating Temperature, °C: 500 w/ FG Flat Gasket

Series 4560 Mini Reactors, 100-600 mL

These are the most popular of all Parr Stirred

Reactors. Although they are called "Mini" reactors, they offer a range of sizes large enough to work with significant sample sizes, yet small enough to be handled with ease.

They are offered in both fixed head and moveable vessel styles and are available with an FKM O-ring seal for operating temperatures to 225 °C, an FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures to 350 °C.

Choosing the high temperature option (HT) boosts the maximum temperature to 500 °C, but reduces the MAWP to 2000 psi (138 bar). Standard Mini reactors can be converted to high temperature reactors (500 °C max temperature and 2000 psi MAWP) by changing the head assembly (contains cone connections, high temperature valves, and graphoil gasket) replacing the heater with a ceramic fiber heater and replacing the split-ring. Contact Parr Instrument Company for details.



Although internal and external space is limited in these small vessels, gas entrainment impellers, catalyst baskets, condensers and other options are available. All reactors in this series can also be easily converted from one size to another by simply substituting a longer or shorter cylinder with the corresponding internal fittings and heaters.

The support system for these Mini Reactors is designed specifically to provide stability at stirring speeds up to 1700 rpm, in a compact mounting small enough to fit into a laboratory hood. The support system can also be adapted to accept any of the smaller vessels from the 4590 Micro Reactor Series.

> 4560 Reactor Vessels from left to right, 300 mL, 100 mL, 450 mL, 160 mL, and 600 mL.

4566 Bench Top Mini Reactor, 300 mL, Fixed Head, Drop-Band Split-Ring, and a 4848 Controller shown with optional Expansion Modules.

4560

Shaded bar indicates specifications th				AEGA	AFGE	AEGG	4567	4560	AFCCD	AFGGO
Model Number	4561	4562	4563	4564	4565	4566	4567	4568	4566B	4566C
Approximate Volume, mL	300	450	600	160	100	300	450	600	160	100
Maximum Pressure (MAWP)	3000 psi (207 bar, 200 bar for CE orders)									
HT Max. Pressure (MAWP)	2000 psi (138 bar)									
Maximum Temperature										
with FKM 0-ring		225 °C								
with FFKM O-ring					300	0°C				
with PTFE Flat Gasket					350	0°C			1	-
with Flexible Graphite Flat Gasket (HT)		500 °C		N	/A		500 °C		N	/A
Reactor Details										
Mounting Style			Moveable					Fixed Head		
Stand Type					Benc	h Top				
Closure		Spli	t-Ring (6 Cor	npression B	olts for Flat (Gasket, no C	ompression	Bolts for O-	ring)	
Valve Connections					1/8" NP	T Male			_	
Magnetic Stirrer, Model No.					A112	0HC6				
Maximum Torque					16 Inch	-Pounds				
Impeller(s), 4-Blade, 1.38" dia.	1	1	2	1	1	1	1	2	1	1
Stirrer Motor					1/8 hp vari	able speed				
Pressure Gage, Size						ches				
Range, Standard Temperature	0-3000 psi (207 bar)									
Range, High Temperature		0-2000 psi (227 bar)								
Temperature Measurement			Fi	xed Thermo	couple (The		special allov	(s)		
Cooling Coil		Included			cluded		Included		Not In	cluded
Style					Loop (Optior	ı nal Sernentir				
Bottom Drain Valve, optional				-	T (Not avail	-				
Heater Style		Mantle			ip-on		Mantle		Clam	ip-on
Heater Power, Watts	510	590	780		25	510	590	780		25
Heater Style w/ HT Option		Ceramic Fibe			/A		Ceramic Fibe			/A
Heater Power, Watts	800	800	1100		/A	800	800	1100		/A
Electrical Supply	000	000	1100	14	/^	000	000	1100		
Volts, AC					115	/ 230				
Maximum Load, amps, 115 / 230				10 /	5 (14 / 7 Hig		uro)			
Cylinder Dimensions				107	5 (14/7 Hig	gir reinperat				
I.D x Depth, inches	2.5 x 4.0	2.5 x 6.0	2.5 x 8.0	2.5 x 2.0	2.1 x 2.0	2.5 x 4.0	2.5 x 6.0	2.5 x 8.0	2.5 x 2.0	2.1 x 2.0
• *										
Vessel Assembly Weight, pounds	17	21	21	17	18	18	21	22	18	19
Weight of Cylinder, pounds	3.7	4.9	6.2	2.4	3.3	3.7	4.9	6.2	2.4	3.3
Reactor Dimensions					10	. 10				
Width x Depth w/o Controller, inches							00			
Height, inches	28	32	36	28	28	28	32	36	28	28
Weight, pounds	60	63	66	60	60	62	65	68	60	60
Spare Parts Kit				4569M	(4569HT for	High Tempe	rature)			



Index to Options

Page
8
ı 10
14
18
118
124
125
129
131
102
109
114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4560 Mini Reactors, 100-600 mL



Model 4563 Bench Top Mini Reactor, 600 mL, Moveable Head, PTFE Flat Gasket Seal, with vessel disassembled, and a 4848 Controller shown with optional Expansion Modules.

Series 4560 Ordering Guide

The Order No. for the Base System is: 456_-T-SS-M-115-VS.12-2000-4848

A composite identification number to be used when ordering a 4560 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Mo		14 16 i		I Cooling Coil
Vodel No.	Size	Vessel Style	-CL	Internal Cooling Loop
561	300 mL	Moveable Head	K Bottom	Drain Valve
562	450 mL	Moveable Head	-BDV	Bottom Drain Valve (Not Available on High Temperature Vessels)
1563	600 mL	Moveable Head		Bottom Brain Valve (Not Available on right temperature vessels)
1564	160 mL*	Moveable Head	L Control	ller
4565	100 mL *	Moveable Head	-4848	PID Control, Ramp & Soak Programming, Motor Speed Control,
4566	300 mL	Fixed Head	(included in	and Data logging with Software. (RS-485 to USB cable not include
4567	450 mL	Fixed Head	base system)	For use with up to three additional display modules.
4568	600 mL	Fixed Head	-4848B	Same as above but for use with up to six additional display modul
4566B	160 mL*	Fixed Head	-A2110E	Motor Controller
4566C	100 mL*	Fixed Head	-4871	Process Controller (for enhanced control options)
* High tempera	ature option not ava	ilable.	See <u>Chapter 6</u>	for a complete list of controllers and options.
B High Te	mperature Opt	tion	M 4848 C	ontroller Options
-No Symbol		uration: 3000 psi / 207 ba		Tachometer Display Module
-HT	2000 psi / 138 ba	•	-MCM	Motor Control Module w/Tachometer
			-PDM	Pressure Display Module
\sim	/ Maximum Te		-HTM	High Temperature Cut Off Module
-0V	FKM 0-ring, 225		-ETLM	External Temperature Limit Module
-0K	FFKM 0-ring, 300		-MTM*	Motor Torque Module
-T	PTFE Flat Gasket	,	-SVM	Solenoid Valve Module (for cooling control)
-FG	Flexible Graphite	e Flat Gasket, 500 °C	-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging
D Materia	l of Constructi			
\smile			-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-SS	T316 Stainless S	teel	-A2208E	RS-485 Daisy Chain for Multiple Controller (must be used with A1925E6)
-MO	Alloy 400*		40504110	
-IN	Alloy 600		-A3504HC	SpecView Software Package for 4838/4848
-HB	Alloy B-2 / B-3*		* The MITM m	ust be installed in conjunction with the MCM.
-HC	Alloy C-276		N Custom	n Options (List All Desired)
-CS	Alloy 20*	×	-AS	Anchor Stirrer
-TI2	Titanium Grade 2*		-PS	Paddle Stirrer
-TI4	Titanium Grade		-SA	Spiral Stirrer
-ZR702	Zirconium Grade		-GE	Gas Entrainment Stirrer
-ZR705	Zirconium Grade	/05*	-BF	Removeable Baffle Set
* Not available			-SB	Static Catalyst Basket
See <u>page 10</u> of	r <u>24</u> for complete lis	st of available alloys.	-CAD	Internal Catalyst Addition Device
E Magnet	ic Stirrer Drive	•	-XCAD	External Catalyst Addition Device
-M	General Purpose	e, 16 in-Ib	-SCP	Solids Charging Port (Ball Valve)
-FMD1	Footless, 16 in-lb	l	-RC	Reflux Condenser
	Duites Madanial	- 6 0	-RTC	Reflux/Take-Off Condenser
\smile		of Construction	-WJ	Welded Jacket
-MOC Symbol	Indicate Materia	l of Construction	-ABH	Aluminum Block Heater w/ Cooling Channels
G Electric	al Supply			7 for a complete list of optional accessories.
-115	115 VAC			
-230	230 VAC		O Certific	ations
~			-ASME	ASME Documentation
H Motor C			-CE	CE Documentation
-VS .12	Variable Speed,		-P	Parr Certification
-VS .25	Variable Speed,	1/4 hp		
-XP .25	Explosion Proof Variable Speed, 1/4 hp			Parts Kit
-AM .25	Air Motor, 1/4 hp	1	-4569M	Spare Parts Kit for 4560 Series
			-4569HT	Spare Parts Kit for 4560HT Series
I Pressur	e Gage		Diana and di	hat all antions and combinations are not competible with "
			Please note th	hat all options and combinations are not compatible with all models.
-3000 3000	psi / 207 bar	-600 600 psi / 40 ba		·····
	psi / 207 bar psi / 138 bar	-600 600 psi / 40 ba -200 200 psi / 14 ba	T	



Series Number: **4570/80**

Series 4570/80 HP/HT Reactors, 250-5500 mL

Type: High Temperature / High Pressure

Stand: Bench Top, Floor Stand or Cart

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: **250-5500**

Standard Pressure MAWP Rating, psi (bar): **4570: 5000 (345) 4580: 3000 (207)**

Maximum Operating Temperature, °C: **500 w/ FG Flat Gasket** Parr offers a number of different reactors in the

4570 and 4580 Series in a variety of sizes and with different mountings for operating pressures to 5000 psi (345 bar) and temperatures to 500 °C.

4571 - 4572 Reactors

These are 1000 and 1800 mL moveable head reactors mounted in moveable carts.

4577 - 4578 Reactors

These are 1000 and 1800 mL fixed head reactors mounted on a floor stand with pneumatic lift.

4575 - 4576A Reactors

These are 500 and 250 mL bench top reactors with 5000 psi (345 bar) and 500 °C ratings. They are available as either moveable or fixed head vessels.

4581 - 4582 Reactors

These 1 gallon (3.75 L) and 1.5 gallon (5.5 L) reactors are designed for 3000 psi (207 bar) maximum pressure with a 500 °C maximum working temperature and mounted on moveable carts.

4583 - 4584 Reactors

These are 1 gallon (3.75 L) and 1.5 gallon (5.5 L) fixed head reactors mounted on floor stands with pneumatic lift.

All of these reactors have been designed with suitable valves, fittings and modifications to the vessel itself to operate

at the higher temperature and pressure ratings.

Flat, flexible graphite (FG) gaskets are used for temperatures in the higher range to 500 °C. These can be replaced with PTFE gaskets for procedures in which the temperature will not exceed 350 °C.



The Series 4570/80 High Temperature/High Pressure Reactors come in three configurations, from left to right: Bench Top (250 & 500 mL only), Floor Stand, and Cart.

Model 4578 High Temperature/High Pressure Reactor, Floor Stand, 1800 mL, Fixed Head, with Heater and Pneumatic Lift.

4570/80

Model No. – Moveable Vessel	4571	4572	4581	4582	4575B	4576B
Model No Fixed Head	4577	4578	4583	4584	4575A	4576A
Approximate Volume, mL	1000	1800	3750	5500	500	250
Maximum Pressure, MAWP	5000 psi	(345 bar)		0 psi ar for CE orders)		(345 bar)
Maximum Temperature			(20) 541, 200 54			
with Flexible Graphite Flat Gasket			50	0 °C		
Reactor Details						
Stand Type	Cart / Flo	oor Stand	Cart / Fl	oor Stand	Benc	h Top
Closure (Compression Bolts)	Split-Ring (12 Co for Flat	ompression Bolts Gasket)	Split-Ring (16 Co for Flat	ompression Bolts : Gasket)	Split-Ring (8 Compression Bol for Flat Gasket)	
Valve Connections			1/4" NP	T Female		
Magnetic Stirrer, Model No.	A11	20HC	A11	80HC	A11	20HC
Maximum Torque	16 Inch	-Pounds	60 Inch	-Pounds	16 Inch	-Pounds
Impeller(s), 6-Blade (4-blade)	2 (2'	ˈdia.)	2 (3.5	5" dia.)	2 (.81" / 1	I.38" dia.)
Stirrer Motor, Variable Speed	1/4	l hp	1/2	2 hp	1/8 hp	
Pressure Gage, Size			4.5 inches			
Range	0-5000 psi (345 bar) 0-3000 psi (207 bar)				0-5000 psi (345 bar)	
Temperature Measurement			Thermowell			
Cooling Coil			Incl	uded		
Style	Serp	Serpentine Serpentine Single Loc				e Loop
Bottom Drain Valve			1" NPS (Floor stand supports only)			
Heater Style		Ca	lrod: Fixed Head,	Ceramic: Removea	ble	
Heater Power, Watts	28	300	28	800	1800	1400
Electrical Supply						
Volts, AC	2	30	2	30	115 / 230	
Maximum Load, amps, 115 / 230	1	15		15	9 /	15
Cylinder Dimensions						
Inside Diameter, inches	3.75 x 6.2	3.75 x 10.5	5.5 x 9.7	5.5 x 15.1	2.5 x 6.6	2.5 x 3.2
Moveable Vessel Assembly Weight, pounds	83	100	124	143	44	38
Fixed Vessel Assembly Weight, pounds	92	109	136	155	44	38
Cylinder Weight, pounds	36	52	53	72	15	9
Reactor Dimensions						
Moveable Width x Depth x Height, in.	Cart: 35 x 18 x 50 Cart: 35 x 18 x 50			17 v 21 v 12		
Fixed Width x Depth x Height, in.	Floor Stand	: 20 x 31 x 75	Floor Stand	: 20 x 31 x 75	Bench Top: 17 x 24 x 4	
Moveable Weight, pounds	245	265	325 355		125	120
Fixed Weight, pounds	265	285	345	375	127	122
Spare Parts Kit	457	79M	45	89M	457	/9M

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.

Series Number: **4570/80**

Series 4570/80 HP/HT Reactors, 250-5500 mL

Model 4581 High Temperature/High Pressure Reactor, with 1 Gallon Moveable Vessel, and a 4848 Reactor Controller with optional Expansion Modules.

Index to Options

Option	Page
Certification	8
Materials of Construction	n 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
 Stirrers	120
 Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

54



Model 4576A HP/HT Bench Top Reactor, 250 mL Vessel with Fixed Head.

Series 4570/80 Ordering Guide

The Order No. for the Base System is: 457_-FG-SS-M-115-VS._-5000-4848 or 458_-FG-SS-HD-230-VS.50-3000-4848

A composite identification number to be used when ordering a 4570 or 4580 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Mo			H Pressur	
Model	Size	Vessel Style	-5000	5000 psi / 345 bar (4570 Series Only)
4571	1000 mL	Moveable Head, Moveable Cart	-3000	3000 psi / 207 bar
4572	1800 mL	Moveable Head, Moveable Cart	-2000	2000 psi / 138 bar
4577	1000 mL	Fixed Head, Stationary Floor Stand	-1000	1000 psi / 69 bar
4578	1800 mL	Fixed Head, Stationary Floor Stand	Bottom	Drain Valve
4575B	500 mL	Moveable Head, Bench Top	-BDV	Bottom Drain Valve (Not available on 250 mL or Bench
4576B	250 mL	Moveable Head, Bench Top	-001	Top Systems)
4575A	500 mL	Fixed Head, Bench Top		
4576A	250 mL	Fixed Head, Bench Top	J Control	
4504	100 1		-4848 (included in	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software. (RS-485 to USB cable not includ
4581	1.0 Gallon	Moveable Head, Moveable Cart	base system)	For use with up to three additional display modules.
4582	1.5 Gallon	Moveable Head, Moveable Cart	-4848B	Same as above but for use with up to six additional display modu
4583	1.0 Gallon	Fixed Head, Stationary Floor Stand	-A2110E	Motor Controller
1584	1.5 Gallon	Fixed Head, Stationary Floor Stand	-4871	Process Controller (for enhanced control options)
B Gasket	/ Maximum	Temperature		for a complete list of controllers and options.
FG		ohite Gasket, 500 °C		
·Т	PTFE Compr	ession Gasket, 350 °C		ontroller Options
			-TDM	Tachometer Display Module
	Is of Constr		-MCM	Motor Control Module w/Tachometer
SS	T316 Stainles	SS STEEL	-PDM	Pressure Display Module
-M0	Alloy 400*		-HTM	High Temperature Cut Off Module
-IN	Alloy 600		-ETLM	External Temperature Limit Module
-HB	Alloy B-2 / B	-3*	-MTM*	Motor Torque Module
-HC	Alloy C-276		-SVM	Solenoid Valve Module (for cooling control)
-CS	Alloy 20*		-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data loggin
-TI2	Titanium Grade 2*		-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-TI4	Titanium Grade 4*		-A2208E	RS-485 Daisy Chain for Multiple Controller
-ZR702	Zirconium Grade 702*		, LLLOOL	(must be used with A1925E6)
-ZR705	Zirconium G		-A3504HC	SpecView Software Package for 4838/4848
*Maximum ter	perature and/o	r pressure limited. <i>e list of available alloys.</i>	* The MTM m	ust be installed in conjunction with the MCM.
~ ~ ~	_ ·	,		Options (List All Desired)
	ic Stirrer Dr		-AS	Anchor Stirrer
-M		oose, 16 in-lb	-PS	Paddle Stirrer
-HD	Heavy Duty,		-SA	Spiral Stirrer
-XHD		Duty, 120 in-lb *	-GE	Gas Entrainment Stirrer
-FMD1		neral Purpose,16 in-lb.	-BF	Removeable Baffle Set
-FMD2		avy Duty, 60 in-lb **	-SB	Static Catalyst Basket
-FMD3	Footless, Ext	ra Heavy Duty, 120 in-lb *	-DB	Dynamic Catalyst Basket
 Not available 	e on 4575-76 (25	0 mL & 500 mL)	-CAD	Internal Catalyst Addition Device
E Mag. Dr	vivo Matoria	of Construction	-XCAD	External Catalyst Addition Device
-MOC Symbol		terial of Construction	-SCP	Solids Charging Port (Ball Valve)
			-RC	Reflux Condenser
F) Electrica	al Supply		-RTC	Reflux/Take-Off Condenser
-115	115 VAC (457	'5-4576 Only)	-WJ	Welded Jacket
-230	230 VAC		See <u>Chapter 7</u>	for a complete list of optional accessories.
	Option		M Certific	ations
G Motor C	perom	ed, 1/8 hp (4575-4576 Only)	-ASME	ASME Documentation
G Motor C	Variable Sno			CE Documentation
-VS .12			-CE	
-VS .12 -VS .25	Variable Spe	ed, 1/4 hp		Parr Certification
-VS .12 -VS .25 -VS .50*	Variable Spe Variable Spe	ned, 1/4 hp ned, 1/2 hp	-P	Parr Certification
VS .12 VS .25 VS .50* XP .25	Variable Spe Variable Spe Explosion Pr	ned, 1/4 hp ned, 1/2 hp oof Variable Speed, 1/4 hp	-P N Spare F	Parr Certification Parts Kit
VS .12 VS .25 VS .50* XP .25 XP .50*	Variable Spe Variable Spe Explosion Pr Explosion Pr	red, 1/4 hp red, 1/2 hp oof Variable Speed, 1/4 hp oof Variable Speed, 1/2 hp	-P	Parr Certification Parts Kit Spare Parts Kit for 4570 Series
VS .12 VS .25 VS .50* XP .25 XP .50* AM .25	Variable Spe Variable Spe Explosion Pr Explosion Pr Air Motor, 1/	red, 1/4 hp ned, 1/2 hp oof Variable Speed, 1/4 hp oof Variable Speed, 1/2 hp 4 hp (4575-4576 Only)	-P N Spare F	Parr Certification Parts Kit
-VS .12 -VS .25	Variable Spe Variable Spe Explosion Pr Explosion Pr Air Motor, 1/ Air Motor, 1/	red, 1/4 hp ned, 1/2 hp oof Variable Speed, 1/4 hp oof Variable Speed, 1/2 hp 4 hp (4575-4576 Only)	-P • Spare F -4579M -4589M	Parr Certification Parts Kit Spare Parts Kit for 4570 Series



Type: Micro

Stand: Bench Top

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: **25, 50, 100**

Standard Pressure MAWP Rating, psi (bar): **3000 (207)**

High Pressure (HP) MAWP Rating, psi (bar): 5000 (345)

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

High Temperature (HT) / High Pressure (HP) Maximum Operating Temperature, °C: 500 w/ FG Flat Gasket (Fixed Head Only)



Smaller vessel sizes are available in our **Series 2500 Micro Batch System** which includes three reactors in 5 mL or 10 mL. <u>See</u> page 72.

Series 4590 Micro Stirred Reactors, 25-100 mL

These are Parr's smallest individual

Stirred Reactors. They are a good choice for chemists working with very expensive materials, materials only available in small amounts, or those chemists working with hazardous materials. Use of the Series 4590 Reactors helps minimize the quantities of waste products created which may require special disposal procedures. The 4590 systems are now available in high pressure (5000 psi) or high pressure/high temperature (500 °C) versions.

These micro reactors have been designed to provide as many of the features of the larger vessels as possible in the limited space available. All of the standard head fittings are provided, with optional cooling now available with aluminum block heaters with cooling capability, welded jackets or internal cold fingers.

These systems are offered in both fixed head and moveable head were a solf

vessel styles with choices of a self sealing FKM O-ring for temperatures up to 225 °C or with FFKM O-ring for temperatures to 300 °C or with a flat PTFE gasket for temperatures to 350 °C, or with an Flexible Graphite gasket for temperatures to 500 °C with the HP/HT option (for fixed head versions only). A split-ring closure is standard.

All three volumes use the same ceramic fiber heater.

These micro reactors can be easily converted from one size

to another by simply substituting a larger or smaller cylinder and the corresponding internal fittings. The support system can also be readily adapted to accept any of the vessels from the 4560 Mini Reactor Series. The opportunity to modify these small reactors is restricted because of the limited head space available.

4598 HP/HT Micro Reactor, 100 mL with Fixed Head.

4590

Model Number	4591	4592	4593	4596	4597	4598
Approximate Volume, mL	25	50	100	25	50	100
Maximum Pressure (MAWP)		3	, 000 psi (207 bar, 20	0 bar for CE order	s)	
HP, HP/HT Max. Pressure (MAWP)			5000 psi	(345 bar)		
Maximum Temperature						
with FKM 0-ring			225	o °C		
with FFKM O-ring			300	0°C		
with PTFE Flat Gasket			350	0°C		
with FG Flexible Graphite Flat Gasket (HP/HT)			500) °C		
Reactor Details				·		
Mounting Style		Moveable			Fixed Head	
Stand type			Benc	h Top		
Closure	Spli	it-Ring (6 Compres	sion Bolts for Flat (Gasket, No Compre	ession Bolts for O-	ring)
Valve Connections		1/8" NPT Male	e, 1/4" NPT Male fo	r HP, 1/4" NPT Fen	nale for HP/HT	
Magnetic Stirrer, Model No.		A1120HC6				
Maximum Torque	16 Inch-Pounds					
Impeller, 4-Blade			1 (.81	" dia.)		
Stirrer Motor			1/8 hp Vari	able Speed		
Pressure Gage, Size			3.5 in	iches		
Range		0-3000	psi (207 bar), 0-500	00 psi (345 bar) HP	, HP/HT	
Temperature Measurement		Fixed T	hermocouple (The	rmowell for specia	al alloys)	
Cooling (optional)			Cold I	Finger		
Bottom Drain Valve		1/4" N	IPT (Not available	on HP or HP/HT m	odels)	
Heater Style			Ceramic Fi	ber Heater		
Heater Power, Watts			70	00		
Electrical Supply						
Volts, AC			115,	/ 230		
Maximum Load, amps, 115 / 230		7 / 4			7 / 4	
Maximum Load, amps, HP/HT, 115 / 230		N/A			7 / 4	
Cylinder Dimensions						
I.D. x Depth, inches	1.0 x 2.0	1.3 x 2.3	1.3 x 4.6	1.0 x 2.0	1.3 x 2.3	1.3 x 4.6
Vessel Assembly Weight, pounds	18	18	19	19	19	20
Cylinder Weight, pounds	3.5	3.1	4.2	3.5	3.1	4.2
Reactor/Stands Dimensions						
Width x Depth w/o Controller, inches			12 :	k 18		
Height, inches*	2	27	35	2	27	35
Weight, pounds	Į	55	56	5	56	57
Spare Parts Kit		45	99M (4599HPHT fo	or High Temperatu	re)	

Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.



Index to Options

Option	Page
Certification	8
Materials of Construction	ו 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fitting	js 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	e 128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Page
129
137
133
125
127
126

Series 4590 Micro Stirred Reactors, 25-100 mL



Series 4590 Ordering Guide

The Order No. for the Base System is: 459_-T-SS-M-115-VS.12-3000-4848

A composite identification number to be used when ordering a 4590 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, <u>please see page 27</u>.

A Base M	lodel		Pressur	e Gage
Model	Size	Vessel Style	-5000	5000 psi / 345 bar
4591	25 mL	Moveable Head	-3000	3000 psi / 207 bar
4592	50 mL	Moveable Head	-2000	2000 psi / 138 bar
4593	100 mL	Moveable Head	-1000	1000 psi / 69 bar
4596	25 mL	Fixed Head	-600	600 psi / 40 bar
4597	50 mL	Fixed Head	-200	200 psi / 14 bar
4598	100 mL	Fixed Head	-100	100 psi / 7 bar
B Pressu	re / Tempe	erature Options	J Bottom	Drain Valve
-No Symbol		207 bar) / 350 °C	-BDV	Bottom Drain Valve (Not available on HP models)
-HP	5000 psi (3	345 bar) / 350 °C		
-HP/HT*	5000 psi (3	345 bar) / 500 °C	K Control	ler
* Fixed Head	Only		-4848	PID Control, Ramp & Soak Programming, Motor Speed Control.
C Gasket	/ Maximu	m Temperature	(included in	and Data logging with Software. (RS-485 to USB cable not included) For use with up to three additional display modules.
-0V	FKM 0-rir	-	base system)	
-0K		ing, 300 °C	-4848B	Same as above but for use with up to six additional display modules
-T		Gasket, 350 °C	-A2110E	Motor Controller
-FG		raphite, 500 °C (HP/HT option only)	-4871	Process Controller (for enhanced control options)
\sim		· · · · · · · · ·	See <u>Chapter 6</u>	for a complete list of controllers and options.
D Materia	al of Cons	truction nless Steel		ontroller Options
-33 -M0	Alloy 400*		-TDM	Tachometer Display Module
			-MCM	Motor Control Module w/Tachometer
-IN	Alloy 600		-PDM	Pressure Display Module
-HB	Alloy B-2 / B-3* Alloy C-276		-HTM	High Temperature Cut Off Module
-HC -CS	Alloy 20*		-ETLM	External Temperature Limit Module
-US -TI2		Crede 2*	-MTM*	Motor Torque Module
-112 -TI4	Titanium (Titanium (-SVM	Solenoid Valve Module (for cooling control)
-714 -ZR702		Grade 702*	-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging)
-ZR702		Grade 702	-A1925E6	RS-485 to USB Converter, isolated, 30-ft
		nd/or pressure limited.	-A2208E	RS-485 Daisy Chain for Multiple Controller (must be used with A1925E6)
See <u>page 10</u> c	or <u>24</u> for com	plete list of available alloys.	-A3504HC	SpecView Software Package for 4838/4848
E Stirrer	Drive			ust be installed in conjunction with the MCM.
-M	General P	urpose Magnetic		
F Mag. D	rivo Mato	rial of Construction	-GE	Options (List All Desired) Gas Entrainment Stirrer
		Vaterial of Construction	-BF	
			-BF -XCAD	Removeable Baffle Set
G Electric	cal Supply	/	-SCP	External Catalyst Addition Device
-115	115 VAC			Solids Charging Port (Ball Valve) Reflux Condenser
-230	230 VAC		-RC -RTC	Reflux/Take-Off Condenser
H Motor	Ontion		-WJ	Welded Jacket
-VS .12		peed, 1/8 hp	-ABH	Aluminum Block Heater with Cooling Channels
-VS .25		peed, 1/4 hp		for a complete list of optional accessories.
-V3 .25 -XP .25		Proof Variable Speed, 1/4 hp	~	
	ATEX Expl	osion Proof Variable Speed		
-XP .25X		D-1300 rpm)	-ASME	ASME Documentation
-AM .25	Air Motor,		-CE -P	CE Documentation Parr Certification
			O Spare P	
			-4599M	Spare Parts Kit for 4590 Series

-4599HPHT

Please note that all options and combinations are not compatible with all models.

Spare Parts Kit for 4590HPHT Series

Type: Glass or Metal

Stand: Bench Top

Mounting Style: Fixed Head

Vessel Sizes: **160 mL - 1.5 L**

Maximum Operating Pressure, MAWP Rating, psi (bar): Glass: 150 (10) Metal: 1000 (69)

Maximum Operating Temperature, °C: **225**





5100 REACTOR

The Parr Series 5100 Low Pressure Glass Reactors offer:

- 1. A system for running reactions similar to those that have been conducted for many years in the Parr shaker, but which offer stirring for better scalability, higher operating temperatures and pressures, and more extensive reactor controls and instrumentation.
- 2. Glass Reactors for use at elevated pressures that permit direct observation of mixing action, color changes, or changes of state.
- Reactors designed for convenient operation at moderate pressure and/or for corrosion resistance.

Convenient and Easy Sealing with O-rings and Split-Ring Closures

Parr has developed an O-ring and closure system to accommodate the requirements of this unique glass-to-metal seal and support, which is convenient to use. A face seal-type O-ring design is used with the proven and popular Parr split-ring closure. For this application a special gasket groove was designed to retain the O-ring on the head of the reactor when it is opened. A full range of O-ring materials is available for chemical compatibility with reactants, products and solvents.

The split-ring for the glass vessel is padded with high temperature cushions so the glass vessel does not come in direct contact with the metal split-ring. Six sealing screws are tightened finger tight to develop the seal on the O-ring. The split-ring snaps together with latches to provide a secure and positive closure.

The alternate metal cylinders use a different split-ring designed to handle the higher working pressure of the metal vessels.

Standard Reactor Fittings

The head of each glass reactor is equipped with:

- Pressure gage, 3-1/2 inch diameter, calibrated 0-200 psi and 0-14 bar
- Gas release valve





Split-Ring for Metal Reactors.



- Gas inlet valve
- Liquid sampling valve
- Internal thermocouple
- Internal cooling loop standard 300 mL to 1.5 L
- Internal stirrer with magnetic drive
- Heads intended for use with glass cylinders are equipped with spring-loaded relief valves adjustable between 50-150 psi.
- All heads are equipped with a rupture disc rated for 1000 psi.
- Internal fittings are T316 Stainless Steel with optional PTFE coating.

Materials of Construction

These reactors are a combination of a glass reaction vessel with a metal head, internal stirrer, dip tube, thermowell, cooling loop, and external valves and fittings, or alternately an all metal system.

The standard material of construction for the head is Type 316 Stainless Steel with PTFE coated T316SS internals. As an alternative the head and internal wetted parts can be provided in any of the standard Parr materials of construction. See the 5100 Ordering Guide, page 63.

Size

Series 5100 Reactors can be easily converted between the 160, 300, 450, and 600 mL sizes by simply changing the cylinders and internal parts. In a similar manner, 1 L and 1.5 L are interchangeable. While the 160-600 mL stand cannot be converted to hold 1 and 1.5 L vessels, the larger stand can be converted to accommodate the 160-600 mL vessels. If you plan to convert at a later time, be sure to order the stand for the largest size you plan to use so the shield and supports will not have to be replaced.

Heating and Temperature Control

These reactors are available with either jacketed or non-jacketed glass or metal cylinders allowing for heating by either a user supplied circulator or with a removable heating mantle respectively. While we would normally expect glass vessels to be equipped with circulating jackets to maintain their transparent feature, some users may not need to heat their reactions or may prefer to use removable heating mantles when they need to work at elevated temperatures. Although transparency is not an issue with metal vessels, users will generally want to select the same heating method for metal vessels as they use for glass vessels so they can utilize the same heating and control system for both.

Users who are using a circulating bath that has its own temperature control for use with these reactors will want to order the reactor without a heating mantle and may only need a motor controller for stirring speed to complete the system. Most laboratory circulators should be adequate for these small vessels. Moveable electric heating mantles are available for vessels that do not have attached circulating jackets. These mantles are for use with 115 or 230 VAC.

5110 Conversion Sets: Glass to Metal or Metal to Glass

Series 5100 Reactors can be easily converted between glass and metal cylinders. The conversion sets include the cylinder, closure, gage and safety relief devices for the "converted to" system. Sets do not include heater.

5110 & 5120 Conversion Sets

Glass to Metal

Catalan Number	C:	Commente From	Comunita To
Catalog Number	Size, mL	Converts From	Converts To
5110A	300	Glass Jacketed	Metal Jacketed
5110B	300	Glass	Metal
5110C	450	Glass Jacketed	Metal Jacketed
5110D	450	Glass	Metal
5110E	600	Glass Jacketed	Metal Jacketed
5110F	600	Glass	Metal
5110M	215	Glass Jacketed	Metal Jacketed
5110N	160	Glass	Metal
5120A	1000	Glass Jacketed	Metal Jacketed
5120B	1000	Glass	Metal
5120C	1500	Glass Jacketed	Metal Jacketed
5120D	1500	Glass	Metal

Metal to Glass

Catalog Number	Size, mL	Converts From	Converts To
5110G	300	Metal Jacketed	Glass Jacketed
5110H	300	Metal	Glass
51101	450	Metal Jacketed	Glass Jacketed
5110J	450	Metal	Glass
5110K	600	Metal Jacketed	Glass Jacketed
5110L	600	Metal	Glass
5110P	215	Metal Jacketed	Glass Jacketed
51100	160	Metal	Glass
5120E	1000	Metal Jacketed	Glass Jacketed
5120F	1000	Metal	Glass
5120G	1500	Metal Jacketed	Glass Jacketed
5120H	1500	Metal	Glass



Series 5100 Glass Reactors, 160-1500 mL

Shaded bar indicates specifications that chan	-			-		
Model Number	5101	5102	5103	5104	5111	5112
Approximate Volume, mL	300	450	600	160	1000	1500
Maximum Pressure, Glass				(10 bar)		
Maximum Pressure, Metal			1000 ps	i (69 bar)		
Maximum Temperature						
with FKM O-ring				5 °C		
with FFKM O-ring, Glass Vessel				5 °C		
with FFKM 0-ring, Metal Vessel			300) °C		
Reactor Details	1					
Mounting Style				Head		
Stand Type				h Top		
Closure	Gla	iss Vessels: 6 Thum			6 Compression B	olts)
Valve Connections			1/8" Ma	ale NPT		
Magnetic Stirrer, Model No.			A112	0HC9		
Maximum Torque		1	16 Inch	-Pounds		
Impeller(s)	1	2	2	1	2	2
Stirrer Motor			1/8 hp vari	able speed		
Pressure Gage, Size	3.5 inches					
Range, Glass Cylinder	0-200 psi (14 bar)					
Range, Metal Cylinder			0-1000 p:	si (69 bar)		
Temperature Measurement			Fixed, Type J,	Thermocouple		
Cooling Coil		Single Loop Include	d	None	Single Loo	p Included
Heater Style			Ma	ntle		
Heater Power Glass, Watts	510	590	780	400	400	550
Heater Power Metal, Watts	510	590	780	400	450	650
Electrical Supply						
Volts, AC			115	/ 230		
Maximum Load, amps, 115 / 230			9	/ 5		
Cylinder Dimensions						
I.D. x Depth, inches	2.5 x 4.0	2.5 x 6.0	2.5 x 8.0	2.5 x 2.0	4.0 x 6.0	4.0 x 8.0
Vessel Assembly Weight, Glass, pounds	15	15	16	14	29	28
Vessel Assembly Weight, Metal, pounds	18	19	21	16	33	36
Cylinder Weight, Glass, pounds	0.8	1.1	1.4	0.5	3.0	4.1
Cylinder Weight, Metal, pounds	3.4	4.6	5.9	2.2	8.3	10.8
Reactor/Stand Dimensions						
Width x Depth w/o Controller, inches		17 :	x 24		21 :	x 26
Height, inches		3	80		3	3
Weight, pounds	60	63	66	60	109	113
Spare Parts Kit			19M	1	511	9M

Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.

Series 5100 Ordering Guide

The Order No. for the Base System is: 51___-G-SS-4B-115-VS.12-200

A composite identification number to be used when ordering a 5100 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Mo	del	K C
Model	Size	
5101	300 mL	-4848
5102	450 mL	-4040
5103	600 mL	
5104	160 mL (215 mL Glass Jacketed)	-4848B
5111	1000 mL	-A2110
5112	1500 mL	-4871
B Cylinder	Type	See <u>Ch</u>
-GJ	Glass Jacketed	L 48
-G	Glass	-TDM
-MJ	Metal Jacketed	-MCM
-M	Metal	-PDM
C Gasket		-HTM
		-ETLM
-0V	FKM O-ring	-MTM*
-0K	FFKM 0-ring	-SVM
D Material	s of Construction	-A1925
-SS	T316 Stainless Steel	
-M0	Alloy 400	-A1925
-IN	Alloy 600	-A2208
-HB	Alloy B-2 / B-3	1004
-HC	Alloy C-276	-A3504
-CS	Alloy 20	* The I
-TI2	Titanium Grade 2	
-TI4	Titanium Grade 4	-PS
-ZR702	Zirconium Grade 702	-GE
-ZR705	Zirconium Grade 705	-XCAD
See <u>page 10</u> or	<u>24</u> for complete list of available alloys.	-SCP
E Magneti	c Stirrer Drive	-RC
-M	General Purpose, 16 in-lb	-RTC
-FMD1	Footless, General Purpose, 16 in-lb	See <u>Ch</u>
E Mag. Dri	ve Material of Construction	N C
-MOC Symbol	Indicate Material of Construction	-ASME
		-CE
G Electrica	I Supply	-P
-115	115 VAC	
-230	230 VAC	-5109M
H Motor O	ption	-5119N
-VS .12	Variable Speed, 1/8 hp	
-VS .25	Variable Speed, 1/4 hp	Please
-XP .25	Explosion Proof Variable Speed, 1/4 hp	with al
-AM .25	Air Motor, 1/4 hp	
Pressure	Corre	
-200 -100	200 psi / 14 bar 100 psi / 7 bar	
-100	100 psi / 69 bar (Metal Vessels Only)	
J Heater		
-MH	Mantle (Non-Jacketed Vessels Only)	
-NH	No Heater	

K Contro	ller
-4848	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software. (RS-485 to USB cable not included) For use with up to three additional display modules.
-4848B	Same as 4848 but with up to six additional modules
-A2110E	Motor Controller
-4871	Process Controller (for enhanced control options)
See <u>Chapter</u>	<u>6</u> for a complete list of controllers and options.
L 4848 C	Controller Options
-TDM	Tachometer Display Module
-MCM	Motor Control Module w/Tachometer
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM*	Motor Torque Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E	RS-485 Daisy Chain for Multiple Controller (must be used with A1925E6)
-A3504HC	SpecView Software Package for 4838/4848
* The MTM n	nust be installed in conjunction with the MCM.
M Custon	n Options (List All Desired)
-PS	Paddle Stirrer, 300 mL and Larger
-GE	Gas Entrainment Stirrer
-XCAD	External Catalyst Addition Device
-SCP	Solids Charging Port (Ball Valve)
-RC	Reflux Condenser
-RTC	Reflux/Take-Off Condenser
See <u>Chapter</u>	<u>7</u> for a complete list of optional accessories.
N Certific	cations
-ASME	ASME Documentation
-CE	CE Documentation
-P	Parr Certification
0 Spare	Parts Kit
-5109M	Spare Parts Kit for Models 5101, 5102, 5103, 5104
-5119M	Spare Parts Kit for Models 5111 and 5112
Please note t	hat all ontions and combinations are not compatible

e note that all options and combinations are not compatible all models.



Type: High Pressure, Compact

Stand: Bench Top

Mounting Style: Moveable

Vessel Sizes, mL: **25-600**

Standard Pressure MAWP, psi (bar): **3000 (207)**

Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

Series 5500 HP Compact Reactors, 25-600 mL

Familiar Parr Designs

All of the safety, convenience and reliability features which have been the hallmark of Parr pressure reaction equipment for more than 50 years have been incorporated into a new line of high pressure, compact laboratory reactors.

The Series 5500 Compact Reactors are based upon our popular micro and mini, Series 4590 & 4560 Reactors. There are several differences between these new reactors and their original counterparts. Primarily, these are:

- A smaller, more compact magnetic drive is installed.
- A smaller, more compact variable speed stirrer motor.
- The larger support stand, overarm and motor have been eliminated.
- An Aluminum Block Heater also serves as the support stand.

As a result of these changes, we are able to offer these reactors to users who do not require the wide variety of options and expandability provided by our Series 4560 and 4590 Reactors at a significantly lower cost. These designs will be welcomed by not only scientists with limited space or budgets, but also by those building multiple reactor installations for combinatorial chemistry or high throughput investigations.

The reaction vessels used in these reactors are identical to the ones furnished in the Series 4590 Micro and Series 4560 Mini Reactors and use the popular Parr split ring closure. These vessels are rated for a maximum working pressure of 3000 psi. The maximum operating temperature is dependent upon the seal selected, PTFE gasket for up to 350 °C; with FKM O-ring to 225 °C or FFKM O-ring to 300 °C.

The 25, 50, and 100 mL reactors are equipped with gas inlet and outlet valve, a liquid sampling valve, pressure gage, safety rupture disc, and internal thermocouple in addition to the internal stirrer. The 300, 450, and 600 mL reactors provide an

optional internal cooling loop in addition to these fittings. In addition to the standard Type 316 Stainless Steel, the vessels for these reactors can also be constructed from any of the standard Parr materials of construction.

These vessels are designed, built and can be certified to the ASME Pressure Vessel Code, European Community P.E.D. and other appropriate local codes. Electrical safety is certified by CSA and the EC-CE mark.

New Magnetic Drive

To take advantage of the new technology available in magnets today, Parr has designed a new compact, magnetically-coupled stirrer drive especially for these smaller vessels. Tests show that this new drive is sufficient to stir reaction mixtures with viscosities up to 10,000 centipoise in a 600 mL reaction vessel

Variable Speed Motor

A 1/17 hp variable speed motor provides stirring speeds adjustable from 0 to 1700 rpm. An optional tachometer pickup provides a signal to the optional tachometer display module which can be installed in the 4848 Controller.

New Heater / Reactor Support

A new heater that also serves as the vessel support has been designed for these reactors. This is an aluminum block style heater for excellent thermal uniformity. The cartridge heaters used in this heating block are easily replaced if required. A stainless steel heat shield is provided around the heating block. This style of heater/reactor support provides a very

small footprint, ideal for limited bench space.

Model 4848 Temperature Controller

The Series 4848 Controller used with the standard Parr line of medium and high pressure reactors is also furnished for use with these reactors. The 4848 offers the user options for redundant temperature sensor and alarm, digital

#11: 1

Parr Series 5500 HPCL Reactor and a 4848 Reactor Controller shown with optional Expansion Modules.

5500

Shaded bar indicates specifications that	it change with	in series.						
Model Number	5511	5512	5513	5521	5522	5523	5524	5525
Sizes, mL	25	50	100	300	450	600	160	100
Maximum Pressure			3000	psi (207 bar, 2	00 bar for CE or	ders)		
Maximum Temperature								
with FKM 0-ring				22	5 °C			
with FFKM O-ring				30	0 °C			
with PTFE Flat Gasket	350 °C							
Reactor Details								
Nounting Style				Mov	/eable			
Stand Type				Compact	Bench Top			
Closure		Split-Ring	(6 Compressio	n Bolts for Flat	Gasket, No Cor	npression Bolts	for O-ring)	
Valve Connections	1/8" Male NPT							
Maximum Torque	2.5 Inch-Pounds (0.28 Nm)							
Impeller(s), 4-blade	1 (0.81" dia.) 1 (0.81" dia.) 2 (1.38" dia.) 1 (0.81" dia.)							
Stirrer Motor, Variable Speed	1/17 hp							
Pressure Gage, Size	3.5 inches							
Range	0-3000 psi (207 bar)							
Temperature Measurement	Fixed Thermocouple							
Cooling Coil	Coldfinger (optional) Standard Single Loop Spiral (optional)							
Bottom Drain Valve	N/A							
Heater Style	Aluminum Block (External Cooling optional)							
Heater Power, Watts	7	00	1000	700	10	000	7(00
Electrical Supply								
Volts, AC				115	/ 230			
Maximum Load, amps, 115 / 230	8	/ 4	10 / 5	8 / 4	10	/ 5	8 /	/ 4
Cylinder Dimensions						`		
I.D. x Depth, inches	1.0 x 2.0	1.3 x 2.25	1.3 x 4.5	2.5 x 4.0	2.5 x 6.0	2.5 x 8.0	2.5 x 2.0	2.0 x 2.0
Vessel Assembly Weight, pounds	17	17	18	17	18	20	16	16
Cylinder Weight, pounds	3.5	3.1	4.2	3.7	4.9	6.2	2.4	3.3
Reactor/Stand Dimensions								
Width x Depth w/o Controller, inches				8.3	x 9.5			
Height, inches	23	23	25	23	25	27	23	23
Weight, pounds	25	25	25	26	28	30	25	25
Spare Parts Kit	5509M							

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information. Weights and dimensions are estimated from the base model. Final weights and dimensions will vary based on options selected.

pressure readout, stirring speed display or control, motor load, and bi-directional digital communication (RS-485).

Alternate Controllers Available

A single 4871 Process Controller can control up to eight high pressure, compact laboratory reactors.

Options

As shown in the ordering guide, a variety of options are available for these Series 5500 Reactors. In addition to the

options described here, there are a number of additional accessories such as glass or PTFE liners, special stirrers, gages, gas and liquid feed systems, custom valves, etc., as described in the Options Section of the Ordering Guide.

The Series 5500 Reactors have been designed and packaged to provide the basic functions of a small laboratory reactor and not all of the options available for the more versatile Series 4560 and 4590 Reactors can be incorporated into these units.



Series 5500 Ordering Guide

The Order No. for the Base System is: 55____-T-SS-115-VS-M-3000-4848

A composite identification number to be used when ordering a 5500 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Model			
Model No.	Size	Cylinder, I.D.	
5511	25 mL	1.0-inch	
5512	50 mL	1.3-inch	
5513	100 mL	1.3-inch	
5521	300 mL	2.5-inch	
5522	450 mL	2.5-inch	
5523	600 mL	2.5-inch	
5524	160 mL	2.5-inch	
5525	100 mL	2.0-inch	

B Gasket / N	laximum Temperature
-0V	FKM O-ring, 225 °C
-0K	FFKM O-ring, 300 °C
-T	PTFE Flat Gasket, 350 °C

C Materials	of Construction (MOC)	
-SS	T316 Stainless Steel	
-M0	Alloy 400	
-IN	Alloy 600	
-HB	Alloy B-2 / B-3	
-HC	Alloy C-276	
-CS	Alloy 20	
-TI2	Titanium Grade 2	
-TI4	Titanium Grade 4	
-ZR702	Zirconium Grade 702	
-ZR705	Zirconium Grade 705	
See name 10 or 24 for complete list of available alloys		

See <u>page 10</u> or <u>24</u> for complete list of available alloys.

D Electrical	Supply	
-115	115 VAC	
-230	230 VAC	

	IVIOT
-VS	

```
or
       1/17 HP Variable Speed
```

F Magnetic Stirrer Drive			
-M	General Purpose Magnetic		

Wag. Drive	materials of Construction
-MOC Symbol I	ndicate Material of Construction

-	
H Pressur	e Gage
-3000	3000 psi / 207 bar
-2000	2000 psi / 138 bar
-1000	1000 psi / 69 bar
-600	600 psi / 40 bar
-200	200 psi / 14 bar
-100	100 psi / 7 bar
I Control	ler
-4848 (included in base system)	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software. (RS-485 to USB cable not included) For use with up to three additional display modules.
-4848B	Same as above but for use with up to six additional display modules.
-A2110E	Motor Controller
-4871	Process Controller (for enhanced control options)
See <u>Chapter 6</u>	for a complete list of controllers and options.
J 4848 C	ontroller Options
-TDM	Tachometer Display Module
-MCM	Motor Control Module w/Tachometer
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM*	Motor Torque Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging)

-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E	RS-485 Daisy Chain for Multiple Controller (must be used with A1925E6)
-A3504HC	SpecView Software Package for 4838/4848

* The MTM must be installed in conjunction with the MCM.

K Custom	Options (List All Desired)
-GE	Gas Entrainment Stirrer
-BF	Removeable Baffle Set
-CAD	Internal Catalyst Addition Device
-XCAD	External Catalyst Addition Device
-SCP	Solids Charging Port (Ball Valve)
-RC	Reflux Condenser
-RTC	Reflux/Take-Off Condenser
See Chanter 7	for a complete list of ontional accessories

See <u>Chapter /</u> for a complete list of optional accessories.

L Certifications		
-ASME	ASME Documentation	
-CE	CE Documentation	
-P	Parr Certification	

(M) Spare F	Parts Kit
-5509M	Spare Parts Kit for 5500 Series

Please note that all options and combinations are not compatible with all models.

Multi Reactor Systems Chapter 3



Inside this chapter you will find: 5000 Multiple Reactor System (MRS) 2500 Micro Batch Reactor System Parallel Reactor Systems



Type: Multiple Vessel

Stand: Bench Top

Vessel Mounting: **Moveable**

Vessel Sizes, mL: **45 and 75**

Standard Pressure MAWP Rating, psi (bar): **3000 (200)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 300 w/ PTFE Flat Gasket

5000 Multiple Reactor System (MRS)

The Parr Series 5000 Multiple Reactor System

has been designed to provide an integrated system for running multiple reactions simultaneously and applying the principles of high throughput experimentation to reactions conducted at elevated temperatures and pressures.

The principal features of the 5000 System include:

- Six reactors with internal stirring.
- Operating pressures to 3000 psi.
- Operating temperatures to 300 °C.
- Individual temperature control.
- Continuous individual pressure monitoring.
- Computer control and data logging.
- Manifold system for rapid turn around and to allow two different input gases.
- Volumes and reactor geometry designed for three phase reactions.
- Flexible Control Software.

Stirred Batch Reaction Vessel

This multiple reaction system has been designed around a vessel with 75 mL total volume. This will accommodate between 15 mL and 40 mL of liquid reactants which is close to the minimum volume appropriate for heterogeneous catalytic reactions.

The vessel valves and accessories are designed for maximum operating pressures up to 3000 psi at operating temperatures up to 300 °C. A system with 45 mL vessels is also available.

Stirring System

All six vessels are stirred with a single magnetic stirrer system specifically designed and built for this application. The magnetic drives and fields are focused on the stirrer bars within each vessel. High strength compact magnets are used to provide coupling forces which will operate through the heaters and vessels. The stirring speed of the stirrer bar is variable from 0-1200 rpm. All vessels will have the same stirring speed during a single run of the apparatus. This arrangement ensures that the difference in reaction rates between vessels run in parallel are due to set conditions other than variations in stirrer speed.

Heaters

ADDATE PRACTICE

The external heaters surround the vessel walls for rapid and uniform heating and temperature control. Each vessel is individually temperature controlled. The 250-watt heater used on each vessel produces heating rates up to 15 °C per minute. An optional reactor cooling support rack is available for air-cooling.



Model 5000 Multiple Reactor System, six 75 mL Vessels with Flat Gaskets and Head-Mounted Valves, shown with included 4871 Process Controller.

Shaded bar indicates specifications th Model Number			000	
		5000 45 75		
Approximate Sizes, mL No. of Reaction Vessels	4	-	75	
Maximum Pressure (MAWP)	6 3000 psi (200 bar)			
		3000 psi	(200 bar)	
Maximum Temperature			E 00	
with FKM O-ring	225 °C			
with FFKM O-ring	300 °C 300 °C			
with PTFE Flat Gasket		30	U ⁻ C	
Closure		0	0	
with O-ring	Screw Cap			
with Flat Gasket	Screw Cap (6 Compression Bolts)			
Material of Construction	T316SS			
Process Controller	Model 4871			
Analog Inputs	6 Temperature			
	6 Pressure			
	1 Motor Speed			
Analog Outputs	1 Stirrer Speed (Optional)			
Digital Outputs	6 PID Temperature Control			
Temperature Measurement	6 Thermowells			
Heater Style	6 Band Heaters, Aluminum Block			
Heater Power, Watts	250W Per Station, 1500W Total			
External Thermocouple	Optional			
Stirrer Motor Type	Manual or Computer Controlled			
Stirrer Style	PTFE-	or Glass-Coate	d Magnetic Sti	rrer Bar
Electrical Supply				
Volts, AC	115 / 230			
Maximum Load, amps, 115 / 230		15	/ 7.5	
Vessel Dimensions				
Inside Diameter, inches	1.18 1.50		.50	
Inside Depth, inches		2.69 Flat Gasl	ket, 2.50 O-ring	
Weight of Vessel, pounds	3 6 (w/ Panel-Mounted Valves) (w/ Head-Mounted Valves)			
Dimensions	Width, in.	Depth, in.	Height, in.	Weight, Ibs
Heater	25.75	9.25	2.875	31
Stirrer	28	9.5	7.625	12
4871 Controller	13	11	15	14
Manifold, Remote	26.5	9.0	15	36
Manifold, Head Mount	26.5	9.0	15	18
Spare Parts Kit		50	09M	

Operating Modes

The Series 5000 Multiple Reaction System provides an apparatus for running up to six reactions in parallel to build a database for comparing and optimizing operating conditions. The user can design experiments to:

- Run all reactions at the same temperature and pressure while varying catalyst loading or reactant concentrations to optimize these parameters.
- Run all reactors with identical loads varying pressures at a common temperature to study the effect of pressure on reaction rates.
- Run individual reactors with individual loading and temperature and pressure to screen multiple options for activity.

A comparison of the plots of pressure drop versus time within the reactors running under parallel conditions will usually be the most useful means of measuring reaction rates and comparing operating conditions. The internal thermocouple also provides a means of detecting temperature changes due to exothermic reactions.

Reactor Options

As Parr customers have come to expect with our line of laboratory pressure reactor equipment, these reactors are offered with a number of options which permit the user to configure the system to their reactions and intended operating conditions. These options include:

O-ring or Flat Gasket Seals. Vessels with O-ring seals are closed by simply tightening the screw cap down hand tight. The maximum operating temperature will depend upon the O-ring material. When equipped with FKM (Viton®) O-rings operating temperatures up to 225 °C are permitted. By substituting FFKM (Kalrez[®]) O-rings this limit can be raised to 300 °C. Careful consideration of chemical compatibility must also be given when selecting O-ring materials. PTFE gaskets can be used to temperatures up to 300 °C and offer virtually universal chemical compatibility. Six compression bolts are used to develop the sealing forces on the PTFE gaskets in this design.

5000 Multiple Reactor System (MRS)

Head Configuration. Each reaction vessel is equipped with an inlet valve, exhaust valve, safety rupture disc, and pressure transducer in addition to an internal thermocouple. Vessels can be modified to include a dip tube for liquid sampling or a cold finger for cooling. The user can choose to have the valves, transducer and rupture disc mounted on a multi-port adapter connected directly to the vessel head, or remotely mounted on the back panel.

- Head Mounted. The head mounted design makes it possible to remove the pressurized vessels from the heater/stirrer assembly or to prefill the vessels in a remote location. PEEK flexible tubing with a quick slip connector is provided for each vessel for inlet of gas. Tubing is easily removed after gas fill.
- **Panel Mounted.** The remote panel mounted arrangement connects all the pressure inlets/ outlets to each vessel with a single PTFE lined stainless tubing. Alternate stainless steel tubing is offered if required. In the panel mount valve configuration the gas inlet tubing generally remains attached to the vessels during operation.

Materials of Construction. Type 316 Stainless Steel is the standard material of construction for



al of construction for both the vessel with its wetted parts and the valve and head fittings exposed to vapors. For investigators working with systems containing strong mineral acids or other more corrosive systems these vessels can be made of most of the Parr standard materials of construction.

Stirrer Configuration.

Stirring is accomplished by use of either PTFE coated or glass coated magnetic stirrer bars.

Thermocouple Configuration.

Thermocouples are mounted inside the vessel for the best temperature monitoring and control. The thermocouples

MRS with individually controlled Gas Burettes. are protected by stainless steel sheaths which are inserted into a protective thermowell. These thermowells make it easy to install and remove thermocouples from the vessels, and also provide additional chemical and mechanical protection for the thermocouple.

We also offer an external thermocouple option with the thermocouples positioned in contact with each cylinder wall.

Gas Manifold. The brass gas inlet manifold is designed to handle both a purge gas, usually nitrogen, and a reactant gas, usually hydrogen. This can be set up to automatically fill each vessel to the same initial pressure or to manually fill each vessel to a unique operating pressure. This gas manifold can also be supplied in T316SS to meet more corrosive requirements.

Gas burettes. For applications where gas consumption measurement is important, optional high pressure gas burettes can be provided (see photo on this page). These burettes allow monitoring/ calculation of the number of moles of gas consumed during a reaction. A complete description and list of available volumes and delivery pressure ranges can be found on page 126 of this catalog (or at https://www.parrinst.com/hpb.

4871 Process Controller

The Series 5000 Multiple Reaction System is controlled by a dedicated Parr 4871 Process Controller. A detailed description of this controller is found in <u>Chapter 6 of this catalog, page</u> 109.

For this application the controller is set up to provide:

- Temperature monitoring and PID control of each individual reactor
- Pressure monitoring of each individual reactor
- Data logging of temperature and pressure in each reactor
- Control and logging of the common stirring speed of the reactors.

The controller provides Ramp & Soak programming for individual reaction vessels, digital inputs and outputs for interlocks, alarms or other safety features, and additional analog and digital inputs and outputs to control flow meters or other accessories which might be added at some future date. The user's control station is a PC running any current Windows operating system. A simplified graphical user interface has been designed for the control and monitoring of the Series 5000 Multiple Reaction System. The PC is used strictly as the user interface and data logging module. All control actions are generated in the 4871 Process Controller (not the PC).

Series 5000 Ordering Guide

The Order No. for the Base System is: 5000(45)-T-SS-115-P-MV-2000-MB-CC

A composite identification number to be used when ordering a 5000 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Model						
Model No. Size						
5000 45 or 75 mL						
-						
B Gasket	B Gasket / Maximum Temperature					
-0V	FKM 0-ring, 225 °C					
-OK FFKM O-ring, 300 °C						
-T PTFE Flat Gasket, 300 °C						

C Materials of Construction						
-SS	T316 Stainless Steel					
-M0	Alloy 400					
-IN	Alloy 600					
-HB	Alloy B-2 / B-3					
-HC	Alloy C-276					
-CS	Alloy 20					
-TI2	Titanium Grade 2					
-TI4	Titanium Grade 4					
-ZR702	Zirconium Grade 702					
-ZR705	Zirconium Grade 705					
See page 10 or 24 for complete list of available alloys.						

D Electric	al Supply
-115	115 VAC
-230	230 VAC

E Thermocouple				
-No Symbol	Type J (Standard)			
-К	Туре К			

F	Valve Mounting					
-H		Head Mount w/ PEEK Tubing				
-P		Remote Panel Mount w/ PTFE Lined SS Hoses				

G Inlet Valve					
-AC	Automatic Check Valve				
-MV	Manual Valve				

H Transdu	cer Range
-1000	0-1000 psi
-2000	0-2000 psi
-3000	0-3000 psi

I Stirrer Type					
-MB	Magnetic Bar Stirrer, PTFE				
-GB	Magnetic Bar Stirrer, Glass				

J Stirring Control						
-M	Manual					
-CC	Computer Controlled					
-RPM	Digital RPM Display					
K Certifica	tions					
-No Symbol	No Certification					
-ASME	ASME Certification					
-PED	PED Certification					
-P	Parr Certification					
L Options						
-SV*	Dip Tube with Sampling Valve					
-CF*	Cold Finger					
-MPG	Manifold Pressure Gage					
-FMH	Flexible SS Hoses					
-E-TC	External Thermocouple					
-R-TC	Redundant Thermocouple					

Gas Burettes (see page 126) * Dip Tubes & Cold Fingers cannot be installed at the same time.

Reactor Cooling Support

M Spare Parts					
-5009M	Spare Parts Kit for 5000 Series				
-5050	Replacement Vessel, T316SS, 75 mL, including head, cylinder, screw cap and thermowell				
-5051	Replacement Vessel, T316SS, 45 mL, including head, cylinder, screw cap and thermowell				

Other Available Options

Glass I	_iner
PTFE L	iner

-RCS

-2280

Parr also designs and builds a wide range of multiple reactor systems with overhead magnetic drive stirrers. Many of these multiple reactor systems have been based upon our Series 4590 Micro Reactors, Series 4560 Mini Reactors, and our Series 5500 High Pressure Compact Reactors. We have also supplied multiple reactor systems based on our largest 4555 Stirred Reactor Series. Please see page 74 of this catalog for examples. Contact our Customer Service Department for details and proposals for custom systems.

71



Type: Micro Batch

Stand: Bench Top

Vessel Mounting: Moveable

Vessel Sizes, mL: **5 and 10**

Maximum Operational Pressure, psi (bar): **3000 (200) Standard**

Standard Maximum Operating Temp., °C: 150 w/ EP 0-ring 225 w/ FKM 0-ring 275 w/ FFKM 0-ring 300 w/ Carbon-Filled PTFE Flat Gasket

Series 2500 Micro Batch Reactor System

The 2500 Micro Batch Reactor System is our smallest ever family of stirred reactors, available in a number of configurations. With reactor volumes of 5 mL or 10 mL, independent reactor gas filling, PID temperature control, and a choice of two closure styles, the flexibility of the MBS 2500 is sure to bring value to your small volume, high pressure project.

Based on the popular Parr MRS 5000, the MBS 2500 offers even smaller volume options, a smaller footprint, a stand-alone controller and data logging system, and a lower per reactor price.

Features of the Micro Batch System 2500 include:

- Three reactors with 5 mL or 10 mL volume, available in stainless steel, high-nickel alloys, titanium, zirconium, and other exotic materials
- Small 12" x 12" (30.5 cm x 30.5 cm) footprint, with remotely located 4848-based controller for temperature/stirring control and pressure display.
- Standard 3000 psi (207 bar, 200 bar for CE) at 300 °C rating in T316 Stainless Steel
- Trusted Parr seal geometry employing either a captive PTFE flat gasket and screw cap with compression bolts, or a self-sealing O-ring with tool-free screw cap closure

- Standard head fittings to include safety rupture disc, pressure transducer, inlet/vent valve, and internal thermocouple
- 3-valve inlet manifold for control of gas addition
- Common aluminum block heater with PID control
- Individual reactor temperature monitoring
- Individual reactor pressure monitoring
- Common stirring speed via individual PTFE- or Pyrex-encapsulated stir bars

Parr Model 2500 Micro Batch Reactor alongside the 4848 Micro Batch System Controller

Series 2500 Ordering Guide

The Order No. for the Base System is: 2500(5)-T-SS-115-RCS-MPG-3000-MCM-4848MBSC

Series 2500 Millio Balcii						
Shaded bar indicates specificatio	ns that chang	e within series				
Model Number	2500		A Base Model			
Sizes, mL	5 10			10	Model No.	Size
No. of Reaction Vessels	3			10	2500	5 or 10 mL
				B Gasket	/ Maximum Temperature	
Maximum Pressure		3000 psi (207 bar, 200 bar for CE)			-0V	FKM 0-ring, 225 °C
Maximum Temperature (nternal)				-0K	FFKM 0-ring, 275 °C
with EP O-ring		1	50 °C		-T	Carbon-filled PTFE Flat Gasket,
with FKM O-ring		2	25 °C		-EP	EP O-ring, 150 °C
with FFKM 0-ring		2	75 °C		C Materia	Is of Construction
with Carbon-Filled					-SS	T316 Stainless Steel
PTFE Flat Gasket		3	0° 00		-M0	Alloy 400
Closure					-IN	Alloy 600
	1				-HB	Alloy B-2 / B-3
with O-ring			ew Cap		-HC	Alloy C-276
with Flat Gasket	5	Screw Cap (6 C	ompression B	olts)	-CS	Alloy 20
Material of Construction		T3	16SS		-TI2	Titanium Grade 2
Controller		Model	4848MBSC		-TI4	Titanium Grade 4
		4 Terr	perature		-ZR702	Zirconium Grade 702
Analog Inputs			l Heater Block)	-ZR705	Zirconium Grade 705
		3 Pi	essure		See page to o	$\frac{24}{24}$ for complete list of available a
		-	or Speed		D Electric	al Supply
			•		-115	115 VAC
Analog Outputs			er Speed		-230	230 VAC
Digital Outputs		-	erature Contro		E Thermocouple	
Temperature Measurement	Indiv	idual Thermov	vell inside eacl	n vessel	-No Symbol	Type J (Standard)
Heater Style		Alumir	ium Block		-K	Туре К
Heater Power, Watts		4	50W		F Transducer Range	
Block Mounted Thermocouple		Inc	luded		-1000	0-1000 psi
Stirrer Motor Type		Manual or Co	nputer Control	led	-2000	0-2000 psi
			·		-3000	0-3000 psi
Stirrer Style	PTFE- or Pyrex-encapsulated Magnetic Stirrer Bars				G Stirrer	Тиро
					-MB	Magnetic Bar Stirrer, PTFE
Electrical Supply	1				-GB	Magnetic Bar Stirrer, Glass
Volts, AC		11	5 / 230			
Maximum Load, amps,		4.6 / 2.3			-No Symbol	
115 / 230		4.07 2.3				No Certification PED Certification
Vessel Dimensions					-PED -P	Parr Certification
Inside Diameter, inches	0.	0.62 0.80				
Inside Diameter, inches	1.06 1.22		Options			
Weight of Vessel, pounds	3 (w/ Panel-mounted Valves)		-PTFE	PTFE Liner		
Dimensions	Width in	Width, in. Depth, in. Height, in. Weight, Ibs.		J Spare F	Parts	
	-	•			-2509M	Spare Parts Kit with three Inco
Main Unit	-	<u>8</u> 10.5 22 41		-2009101	rupture discs, 18 gaskets, and	
4848MBSC Controller Spare Parts Kit	13.5	13.5 11 10 14 2509M			-2550	Replacement vessel, T316SS, 1 including head, cylinder, and s
Other options availab	ther options available. See Ordering Guide, visit <u>www.parrinst.com</u> , or call for more information.			-2551	Replacement vessel, T316SS, 5 including head, cylinder, and s	

B Gasket /	Maximum Temperature
-0V	FKM O-ring, 225 °C
-0K	FFKM 0-ring, 275 °C
-T	Carbon-filled PTFE Flat Gasket, 300 °C
-EP	EP O-ring, 150 °C
	s of Construction
-SS	T316 Stainless Steel
-MO	Alloy 400
-IN	Alloy 600
-HB	Alloy B-2 / B-3
-HC	Alloy C-276
-CS	Alloy 20
-TI2 -TI4	Titanium Grade 2 Titanium Grade 4
-ZR702	Zirconium Grade 702 Zirconium Grade 705
-ZR705	
See paye to or .	24 for complete list of available alloys.
D Electrica	al Supply
-115	115 VAC
-230	230 VAC
E Thermoo	couple
-No Symbol	Type J (Standard)
-К	Туре К
F) Transdu	cer Range
-1000 -2000	0-1000 psi 0-2000 psi
-2000 -3000	0-3000 psi
~	0-0000 poi
G Stirrer T	уре
-MB	Magnetic Bar Stirrer, PTFE
-GB	Magnetic Bar Stirrer, Glass
H Certifica	tions
-No Symbol	No Certification
-PED	PED Certification
-P	Parr Certification
I Options	
-PTFE	PTFE Liner
J Spare Pa	arts
-2509M	Spare Parts Kit with three Inconel rupture discs, 18 gaskets, and lube
-2550	Replacement vessel, T316SS, 10 mL, including head, cylinder, and screw cap
-2551	Replacement vessel, T316SS, 5 mL, including head, cylinder, and screw cap

Series 2500 Micro Batch System Specifications

Parallel Reactor Systems

Parr designs and builds custom parallel reactor and multiple reactor systems for use in high-throughput screening, combinatorial chemistry, corrosion testing, and catalyst testing.

These systems are based on our wide range of reactor systems and feature overhead magnetic drive stirrers. Contact our Customer Service Department for details and proposals for custom parallel reactor systems.



Sixteen Reactor Parallel System

This system is a combination of sixteen standard 4560 Mini Reactors with heaters, valves, pressure gages and rupture disc assemblies and two 4871 Process Controllers with sixteen 4875 Power Controllers. It allows the user to run multiple reactions simultaneously, applying the principles of high-throughput experimentation. Individual variables that can be controlled are gas mixtures, liquids, catalysts or other solids, stirring speed, temperature, pressure and time.



Twelve Reactor Parallel System

This system makes use of the lower cost 5500 High Pressure Compact Lab Reactors that feature a modified stand, aluminum block heaters, removable vessels and a standard gage block assembly. A control system (not pictured) automates the process, monitors the parameters and collects the data.



Five Reactor Parallel System

This Parallel Reactor System incorporates five 160 mL. 4560 Series Stirred Mini Reactors rated for use at 350 °C and 3000 psig / 200 bar. Each of the reactors has a dedicated high pressure (5000 psig / 345 bar) 100 mL general purpose vessel associated with it serving as a reactant gas feed reservoir. These feed vessels are located above and immediately behind the main reactors and are each equipped with a dedicated pressure transducer, thermocouple and a constant pressure regulator. The regulator maintains the downstream reactor at a constant pressure while the pressure transducer/thermocouple combination allows the process controller to accurately monitor and record the real-time consumption of reactant gas. This system would normally be used to study various aspects of hydrogenation or carbonylation reactions. The 4871 Process Controller also controls and records the reactor temperature, stirring speed and monitors the reaction pressure. The compact reactor system stand incorporates not only reactant gas feed and vent/ purge manifolds but cooling water feed and drain manifolds to simplify the installation and hook up.



Three Reactor Parallel System

This system incorporates three 1.8 L 4570 Series High Pressure/High Temperature reactors rated for use at 5000 psi (345 bar) at 500 °C. The reactors in this system are manufactured of Alloy C-276 and are equipped with internal coupon holders to facilitate study of corrosion processes in harsh environments. This particular application does not require precisely controlled gas delivery, so Individual gas feed vessels are not present. However, gas supply and gas vent manifolds are provided. A 4871 Process Controller is also included (not shown).



Four Reactor Parallel System on Moveable cart

This system consists of three 160 mL and one 450 mL 4560 Series Reactors rated for use at 3000 psi (200 bar) at 350 °C. This configuration allows both high throughput screening and scale-up, if desired. Gas supply and vent manifolds for the entire system and pressure control for the larger reactor are provided. Reactor temperatures and stirring speed are controlled by four 4848 Reactor controllers (not shown).



Four Reactor Benchtop Parallel System

This system combines four 300 mL, 4560 Series Mini Bench Top Reactors dual rated for 2900 psi (200 bar) at up to 350 °C and 2000 psi (137 bar) at 500°C. This multiple reactor system was designed for hydrogenation testing in the petroleum industry and includes individually controlled gas feed lines, a liquid pump to sequentially fill the reactors, and an automated vent system. Manifolds for gas supply & vent, as well as cooling water supply and drain are also included. Independent temperature and motor control for each reactor is provided by a 4871 Process Controller (not shown).

Parallel systems continued next page >



Parallel Reactor Systems, continued

< Parallel systems continued from previous page



Six Reactor Parallel System with Automatic Sampling System

This system includes six 25 mL, 4590 Series Stirred Micro Reactors rated for use up to 3000 psi (200 bar) at 350 °C. Gas supply and vent manifolds are provided, as well as manifolds to deliver and drain cooling water to/from either the internal cold finger or the external aluminum block heaters with cooling channels. Control is provided by a 4871 Process Controller.

Note, an automated liquid sampling system installed on the far left reactor allows the sequential collection of several ~1 mL liquid samples under full reactor operating pressure. The system automatically clears the lines between samples.

For more information on automated liquid sampling, please see the 4878 Automated Liquid Sampler System on page 130.

Tubular Reactor Systems Chapter 4



Inside this chapter you will find: 5400 Tubular Reactor System Fluidized Bed Reactor Systems

Series 5400 Continuous Flow Tubular Reactor Systems



Model 5403 with a 1" inside dia. x 24" length, 3-zone split tube furnace with gas feed system, cooling condenser, and gas/liquid separator vessel.

Tubular reactors are used in a continuous flow mode with reagents flowing in and products being removed. They can be the simplest of all reactor designs. Tubular reactors are often referred to by a variety of names:

- Pipe reactors
- Packed-bed reactors
- Fixed-bed reactors
- Trickle-bed reactors
- Bubble-column reactors
- Ebullating-bed reactors

Single-phase flow in a tubular reactor can be upward or downward. Two-phase flow can be co-current up-flow, counter-current (liquid down, gas up) or, most commonly, co-current down-flow.

Tubular reactors can have a single wall and be heated with an external electric furnace or they can be jacketed for heating or cooling with a circulating heat transfer fluid. External furnaces are typically rigid, split-tube heaters. Tubular reactors are used in a variety of industries:


- Petroleum
- Petrochemical
- Polymer
- Pharmaceutical
- Waste Treatment
- Specialty Chemical
- Alternative Energy

Tubular reactors are used in a variety of applications, including:

- Carbonylation
- Dehydrogenation
- Hydrogenation
- Hydrocracking
- Hydroformulation
- Oxidative decomposition
- Partial oxidation
- Polymerization
- Reforming

Tubular reactors may be empty for homogenous reactions or packed with catalyst or other solid particles for heterogeneous reactions. Packed reactors require upper and lower supports to hold particles in place. Upper packing often includes inert material to serve as a pre-heat section. Preheating can also be done with an internal spiral channel to keep incoming reagents close to the heated wall during entry, as shown above.

It is often desirable to size a tubular reactor to be large enough to fit 8 to 10 catalyst particles across the diameter and be at least 40-50 particle diameters long. The length to diameter ratio can be varied to study the effect of catalyst bed length by equipping the reactor with "spools" placed into the bottom of the reactor to change this ratio.





Open 3-zone Split Tube Furnace with 1" I.D. Tubular Reactor.

Tubular reactor systems are highly customizable and can be made to various lengths and diameters and engineered for various pressures and temperatures, and materials of construction.

A split-tube furnace is provided for heating these vessels. Insulation is provided at each end to minimize heat losses and prevent the end caps from being heated. The heater length is normally divided into one, two, or three separate heating zones, although it can be split into more zones if required.

A fixed internal thermocouple in each zone can be furnished or a single movable thermocouple in a centerline thermowell can be used to measure the temperature at points along the catalyst bed. External thermocouples are typically provided for control of each zone of the heater, as can be seen in the photo above.

Gas Feed Systems

Various gas feeds can be set up and operated from a **Gas Distribution Rack**. In order to deliver a steady flow of gas to a reactor, it is necessary to provide gas at a constant pressure to an electronic **Mass Flow Controller**. This instrument will compare the actual flow rate delivered to the set point chosen by the user, and automatically adjust an integral control valve to assure a constant flow.

Series 5400 Tubular Reactor System Specifications

Shaded bar indicates specifications that change within series.

Shaueu bar multates specifications that change within series.				
Model Number	5401	5402	5403	5404
Sizes	3/8 in.	1/2 in.	1.0 in.	1.5 in.
O.D. / I.D. (in.)	0.38 / 0.28	0.50 / 0.37	1.50 / 0.99	1.88 / 1.44
O.D. / I.D. (mm)	9.5 / 7.0	13 / 9.5	38 / 25	48 / 36
Heated Length (in.)	6, 12	2, 24	12, 2	4, 36
Maximum Pressure (psi)	30	00	5000	3000
Maximum Temperature	55	50	550	350
Support Spools No Opti		onal		
Spiral Pre-Heat	N	0	Opti	onal
No. Ports in Top Head	lead 1		4	
No. Ports in Bottom Head	1		4	
Internal Thermocouple	Optional (Moveable or multi-point fixed)			

Care must be taken to size these controllers for the specific gas, flow rate range, and maximum pressure of operation. A mass flow controller needs a power supply and read-out device, as well as a means of introducing the desired set point.

When ordering mass flow controllers, you will need to specify:


- 1. Type of gas to be metered (e.g. N₂, H₂, CH₄)
- 2. Maximum operating pressure of the gas (100 or 300 bar)
- Maximum flow rate range in standard cc's per minute (sccm)
- 4. Pressure for calibration of the instrument


Mass flow controllers are available for use to 1500 psi and to 4500 psi. Considerable savings can be obtained if the mass flow controller is to be used only to 1500 psi.

The schematic at right depicts the installation of a mass flow controller for the introduction of gas to a continuous-flow reaction system. Such installations are enhanced with the addition of a by-pass valve for rapid filling or flushing.

A purge line can also be added. It is typically used for feeding nitrogen or helium to remove air before reaction or to remove reactive gases before opening the reactor at the end of a run. The purge line includes a shut-off valve, filter, metering valve, and a reverse-flow check valve.

Shut-off valves can be automated when using a 4871 Control system.





Series 5400 Continuous Flow Tubular Reactor Systems



This 1 liter Tubular Reactor System has two gas feeds, one purge line, and one liquid feed. Custom pressure controls enhance the heated gas/liquid separator.

Liquid Metering Pumps

High pressure piston pumps are most often used to inject liquids into a pressurized reactor operating in a continuous-flow mode. For low flow rates, HPLC pumps, many of which are rated for 5000 psig, are excellent choices. Typical flow rates for pumps of this type range up to 10 or 40 mL per minute. Pumps are available to accommodate manual control from their digital faceplate or computer-control from a 4871 Process Controller.

Chemical feed pumps are our recommendation for continuous feeding of liquids when the desired flow rate is greater than 2 liters per hour. Parr can assist with the feed pump selection. We will need to know the type of liquid; the minimum, typical, and maximum desired feed rate; the maximum operating pressure; and any special operating considerations such as corrosion possibilities.

Cooling Condensers

Cooling condensers are available to cool the products of the reaction. An adaptation of our standard condensers provides an excellent design.

Back Pressure Regulators

The reactor pressure is maintained by a Back Pressure Regulator (BPR) installed downstream of the reactor. This style of regulator will release products only when the reactor pressure exceeds a value preset by the operator.

When a BPR is used in conjunction with mass flow controllers, the user can maintain a constant flow of gas through a reactor held at an elevated constant pressure. This provides for the highest degree of control and reproducibility in a continuous-flow reactor system.

An alternate BPR may be available to permit pressure let-down of a two-phase stream from the reactor. This BPR requires that the operator provide a source of nitrogen or air at a pressure slightly above the desired operating pressure. With this style of BPR, the highpressure gas/liquid separator can be replaced with a low-pressure liquid product receiver allowing collection of near real-time liquid samples.



This continuous flow stirred reactor system is on a cart with our Modular Frame System. This modular frame allows for easy access and flexibility in hook-ups, accessories, and flow, including an interchangeable tubular reactor.

Gas/Liquid Separators

Tubular reactors operating in continuous-flow mode with both gas and liquid products will typically require a Gas/Liquid Separator. The separator is placed downstream of the reactor, often separated from the reactor by a cooling condenser. In the separator vessel, liquids are condensed and collected in the bottom of the vessel. Gases and non-condensed vapors are allowed to leave the top of the vessel and pass to the back pressure regulator. It is important to operate the standard BPR with a single fluid phase to prevent oscillation of the reactor pressure.

The gas/liquid separator can be sized large enough to act as a liquid product receiver that is drained periodically. Many of the non-stirred pressure vessels made by Parr are ideally suited for use as gas/liquid separators. Vessels of 300, 600, 1000, or 2000 mL are commonly chosen. Upon request, the bottom of the separator can be tapered to facilitate draining.

Control and Data Acquisition Systems

A variety of solutions exist to meet the needs of system operators. System accessories such as heaters, mass flow controllers, and pumps can be obtained with individual control packages to create a Local Control System (LCS) based on our 4838 and 4848 Controllers.


As the number of channels to be controlled increases, economics and convenience will often dictate that the system of individual controllers should be replaced with the computer-based Model 4871 Process Controller (PCC).




This tubular reactor system features a single-zone 12-inch long heater to take reactor temperatures to 500 °C. The gas feed system includes automated shut-off valves on the three MFCs and on the purge line. The Gas/Liquid Separator vessel is jacketed to allow for subambient cooling. Operator controls the system via a remote PC using the included Model 4871 Process Controller.

Series 5400 Continuous Flow Tubular Reactor Systems


On this page are schematic representations of typical tubular reactor systems, along with a symbols chart to facilitate understanding. We have provided an ordering number for each of these examples.



Order No. for this system would be: 5402B-SS-115-ST1(6)-1500-LCS-GF(1)-PL-LF(1)-ITW-CCD-GLS(300)-MPC



Order No. for this system would be: 5403F-SS-230-ST3(24)-3000-PCC-GF(2)-PL-LF(1)-ISP-CSS-ITW-GLS(600)-APC-ASV(3)



Series 5400 Ordering Guide

A composite identification number to be used when ordering a 5400 Series Reactor can be developed by combining individual symbols from the separate sections below.

A Base Model		
Model No.	Size (0.D. / I.D.)	
5401	3/8 in. (0.38" / 0.28")	
5402	1/2 in. (0.50" / 0.37")	
5403	1.0 in. (1.9" / 1.0")	
5404	1.5 in. (2.0" / 1.5")	

Add suffix F for Floor Stand mounting Add suffix B for Bench Top mounting

B Materials of Construction		
T316 Stainless Steel		
Alloy 276		
Titanium		
Alloy 600		
Alloy 400		

C Electrical Supply		
-115	115 VAC, 50/60Hz	
-230	230 VAC, 50/60Hz	
D Heater Options		

-ST1(#)	Split Tube, 1-zone
-ST3(#)	Split Tube, 3-zone
-WJ(#)	Welded Jacket
Add suffix (6), (12), (24), (36) for heated length (in.)	

E Maximum Operating Pressure		
-1500	1500 psi / 100 bar	
-3000	3000 psi / 200 bar	
-4500	4500 psi / 300 bar	

F Controll	er
-PCC	PC-based Process Control (4871-style)
-LCS	Local Control System (4838-style)

G Custom	Options
-GF(#)	Number of Gas Feeds
-PL	Purge Gas Feed Line
-LF(#)	Number of Liquid Feeds
-ISP	Internal Pre-heat Spiral (5403/5404 only)
-CSS	Catalyst Support Spools (5403/5404 only)
-ITW	Internal Thermowell, with Movable T/C
-IZT	Internal, 3-PT, Fixed T/C
-CCD	Cooling Condenser
-GLS(#)	Gas/Liquid Separator (300, 600, 1000, 2000 mL)
-SPH	Separator Heater
-MPC	Manual Pressure Control
-APC*	Automated Pressure Control
-ASV(#)*	Automated Shut-off Valves (1-12)
*Available only with 4871 Process Control (PCC)	

H Certifications		
-No Symbol	No Certification Required	
-PARR	Parr Certification	
-ASME	ASME Certification	
-PED	PED Certification	



The system above has three 250 mL tubular reactors operating in parallel and controlled by a 4871 Process Controller with operator interface on a single PC. This system has weighed and/or heated feed tanks and a two-stage pressure let down.



Fluidized Bed Reactors



Fluidized Bed Reactors are used extensively in the chemical process industries.


The distinguishing feature of a fluidized bed reactor is that the bed of solid particles or catalyst is supported by an up flow of gas. This reactor provides easy loading and removal of catalyst. This is advantageous when the solids bed must be removed and replaced frequently. A high conversion with a large throughput is possible with this style of reactor. Such reactors inherently possess excellent heat transfer and mixing characteristics.

Fluidized beds have been significantly utilized in chemical processes in which parameters such as diffusion or heat transfer are the major design parameters. Compared to packed bed, a fluidized bed has notable advantages such as better control of temperature, no hot spot in the bed, uniform catalyst distribution and longer life of the catalyst.

Nearly all significant commercial applications of fluidized bed technology concern gas-solid systems. Applications of fluidized bed reactors include but are not limited to gas-solid reactions, Fisher-Tropsch synthesis, and catalytic cracking of hydrocarbons, and related high molecular weight petroleum fractions. Gasification in a fluidized bed can be utilized to convert coal, biomass and other waste materials into synthesis gas.


Flexible Mantle Heaters are wrapped around this 36-inch-long, 1-inch I.D. Fluidized Bed Reactor (A) and the Cyclonic Separator (B) to maintain temperatures to 350 °C. This system is also equipped with three gas feeds with automated shut-off valves (C), a cooling condenser (D), a product receiver (E), a back-pressure regulator (F), and a Model 4871 Process Controller (not shown). The reactor system pictured at left (page 84) includes the following key components:

- A gas handling and mixing sub-system used to blend and regulate the flow of reactant gas to the bottom of the reactor.
- A reactor roughly one meter long with a 2.5 cm ID. The lower portion of the reactor incorporates an easily replaced porous metal gas diffusion plate and the top of the reactor widens abruptly to form a disengaging zone for the fluidized bed. Separate heaters are provided for both the main reactor and disengaging zone.
- A moveable thermocouple in a thermowell is provided for monitoring the internal reactor temperature distribution.
- A heated cyclone separator or filter is provided immediately downstream of the reactor to capture the fines resulting from particle attrition.
- The reaction products are then cooled by a condenser and collected in a 600 mL product receiver.
- The system pressure is maintained by an automated, dome-loaded, back pressure regulator.
- All system functions and parameters are monitored and maintained by a Parr 4871 Process Controller (not shown, see Chapter 6, page 109).



Parr Instrument Company

Fluidized Bed Reactors



Ordering Guide

A Base Model		
Model No.		
5410		
	rials of Construction	
-SS	T316 Stainless Steel	
-HC	Alloy 276	
-HT	High Temperature Alloy to be determined (for use up to 900 °C)	
See Materia available all	ls of Construction on <u>page 10</u> for other oys	
C Electi	rical Supply	
-115	115 VAC	
-230	230 VAC	
\sim	num Temperature	
-350	350 °C	
-600	600 °C	
-900	900 °C	
E Cyline	der Length	
-(36)	36-inches	
100/	30-11101163	
-(xx)	Other	
-(xx)	Other	
-(xx)	Other num System Pressure	
-(xx) F Maxir -500	Other num System Pressure 500 psi / 35 bar	
-(xx) F Maxir -500 -1500	Other num System Pressure 500 psi / 35 bar 1500 psi / 100 bar	
-(xx) F Maxir -500	Other num System Pressure 500 psi / 35 bar	
-(xx) F Maxir -500 -1500 -3000	Other num System Pressure 500 psi / 35 bar 1500 psi / 100 bar 3000 psi / 200 bar	
-(xx) F Maxir -500 -1500	Other num System Pressure 500 psi / 35 bar 1500 psi / 100 bar 3000 psi / 200 bar col	
-(xx) F Maxir -500 -1500 -3000 G Contr	Other Other Other System Pressure 500 psi / 35 bar 1500 psi / 100 bar 3000 psi / 200 bar Col PC-based Process Control (4871-style)	
-(xx) F Maxim -500 -1500 -3000 G Contr -PCC -LCS	Other Other num System Pressure 500 psi / 35 bar 1500 psi / 100 bar 3000 psi / 200 bar rol PC-based Process Control (4871-style) Local Control System (4838-style)	
-(xx) F Maxin -500 -1500 -3000 G Contr -PCC -LCS H Custo	Other num System Pressure 500 psi / 35 bar 1500 psi / 100 bar 3000 psi / 200 bar rol PC-based Process Control (4871-style) Local Control System (4838-style) om Options (List all desired)	
-(xx) F Maxir -500 -1500 -3000 G Contr -PCC -LCS H Custo -GF(#)	Other Other Other Solution System Pressure 500 psi / 35 bar 1500 psi / 100 bar 3000 psi / 200 bar O PC-based Process Control (4871-style) Local Control System (4838-style) Om Options (List all desired) Number of Gas Feeds (1-3 or TBD)	
-(xx) F Maxir -500 -1500 -3000 G Contr -PCC -LCS H Custo -GF(#) -PL	Other Other Other Other Other Other Other Other Solo psi / 35 bar Solo psi / 100 bar Solo psi / 200 bar Other PC-based Process Control (4871-style) Local Control System (4838-style) Other Options (List all desired) Number of Gas Feeds (1-3 or TBD) Purge Gas Feed Line	
-(xx) F Maxir -500 -1500 -3000 Contr -PCC -LCS H Custo -GF(#) -PL -LF(#)	Other Other Other Solution Solut	
-(xx) F Maxir -500 -1500 -3000 Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW	Other Other Other Solution Solut	
-(xx) F Maxir -500 -1500 -3000 Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW -IZT	Other Other Other Solution Solut	
-(xx) F Maxir -500 -1500 -3000 Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW	Other Other Other Other Solo psi / 35 bar Solo psi / 35 bar Solo psi / 100 bar Solo psi / 200 bar O PC-based Process Control (4871-style) Local Control System (4838-style) O Options (List all desired) Number of Gas Feeds (1-3 or TBD) Purge Gas Feed Line Number of Liquid Feeds (1-3 or TBD) Internal Thermowell, with Movable T/C Internal, 3-PT, Fixed T/C Cooling Condenser	
-(xx) F Maxir -500 -1500 -3000 Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW -IZT	Other Other Other Solution Solut	
-(xx) F Maxir -500 -1500 -3000 G Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW -ITW -IZT -CCD	Other Other Other Solution Solut	
-(xx) F Maxin -500 -1500 -3000 G Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW -IZT -CCD -GLS(#)	Other Other Other Other Other Solo psi / 35 bar Solo psi / 35 bar Solo psi / 100 bar Solo psi / 200 bar O PC-based Process Control (4871-style) Local Control System (4838-style) O Options (List all desired) Number of Gas Feeds (1-3 or TBD) Purge Gas Feed Line Number of Liquid Feeds (1-3 or TBD) Internal Thermowell, with Movable T/C Internal, 3-PT, Fixed T/C Cooling Condenser Gas/Liquid Separator (300, 600, 1000, 2000 mL)	
-(xx) F Maxir -500 -1500 -3000 G Contr -PCC -LCS H Custo -GF(#) -PL -LF(#) -ITW -IZT -CCD -GLS(#) -SPH	Other Other Other Other Solo psi / 35 bar Solo psi / 35 bar Solo psi / 100 bar Solo psi / 200 bar O PC-based Process Control (4871-style) Local Control System (4838-style) O Options (List all desired) Number of Gas Feeds (1-3 or TBD) Purge Gas Feed Line Number of Liquid Feeds (1-3 or TBD) Internal Thermowell, with Movable T/C Internal, 3-PT, Fixed T/C Cooling Condenser Gas/Liquid Separator (300, 600, 1000, 2000 mL) Separator Heater	

I Certifications		
-No Symbol	No Certification Required	
-ASME	ASME Certification	
-PED	PED Certification	
-P	Parr Certification	

Please note that all options and combinations are not compatible with all models.

Specialty & Custom Reactor Systems Chapter 5

Inside this chapter you will find: Biofuels and Alternative Fuels Research Systems Horizontal Reactors Rotating Cage Apparatus A1000G Tri-Axial Core Holder Supercritical Fluids GTL (Gas-To-Liquid) System Apparatus for Vapor Pressure Determination Disbonding Apparatus for ASTM G146

Apparatus for Corrosion Studies

NUM

Introduction to Specialty Custom Reactors

Part Instrument Company designs and builds a wide variety of specialty and custom reactor systems. These reactor systems embody the spirit of innovation that has made Parr the world leader in laboratory pressure apparatus. Our technical sales, engineering and production teams are happy to work directly with our customers' scientists and/or engineers to very quickly, accurately, and economically proceed from concept to manufacturing. The following pages illustrate several examples of specialty and custom reactor systems designed and built using exactly this collaborative process.

To facilitate collaboration between Parr and our customers as well as to support the high quality of our equipment, we have invested heavily in both design and manufacturing technology. All custom reactor systems are designed with 3D modeling software to create a virtual prototype prior to manufacturing. This model supports clear communication between the customer and Parr and shortens the design phase of such projects. On the manufacturing floor, state of the art, five axis, computer controlled machines programmed and operated by Parr's highly skilled machinists allow rapid and accurate production of even the most complex parts.



Parr uses the latest technology in milling and manufacturing techniques to produce your custom system to the highest level of precision possible in the field today.

To match the variety of reactor systems available, our control design and programming team offers multiple types of control systems, ranging from simple, manual temperature control to completely automated and integrated PC-based control, developed to each customer's specification. Control systems are discussed more completely in <u>Chapter 6</u>.

If you have an idea for a reactor system that is not described in this catalog, please contact Parr's world-class technical sales department — there is a very good chance we can build a reactor system to meet your needs.





Biofuels and Alternative Fuels Research Systems



Custom Reactor Systems like the above Biofuels Research System are a product of collaboration between the researchers and the Engineers at Parr Instrument Company.

Parr Instrument Company manufactures non-stirred vessels for the decomposition of biomass in ammonia and steam. Parr stirred reactors, including a new horizontal reactor technology, have been designed for research processes that include hydrogenation, isomerization, and metathesis reactions. In addition, fully customizable continuous-flow tubular reactor systems have been developed with continuous reactant feed and product handling capabilities.

The above photo illustrates a complete pilot scale plant used for hydrogenating feedstock that originates from a proprietary fermentation process. The system is used to develop and optimize the process conditions necessary for a much larger demonstration-scale system, ultimately leading to full-scale production of renewable fuels.

The system is comprised of five major subsystems: from left to right, a gas and liquid feed system, the jacketed tubular reactor module including a reactant pre-heater and circulating bath, product recovery and backpressure control, an auto-sampler and a Parr 4871 Control System (not pictured). The system is completely automated and includes an auto-sampling subsystem that periodically samples the reactor output stream in order to accurately monitor product quality.

Horizontal Stirred Reactors

Stirring biomass is not easy to do in a vertical reactor. A line of reactors has been developed by Parr Instrument Company to function horizontally. Horizontal reactors enable researchers to mix bio feed stock and other cellulosic materials such as:

- straw - plants


_

- corn stover wood pulp
- sugar cane wood chips
- grasses fibers


In the horizontal position, the internal stirrers tumble the material rather than try to stir it with a standard vertical stirrer. These custom reactors are secured vertically for loading and then safely tilted horizontally for stirring. They can also be designed to be operated in both a horizontal and vertical position. One Liter Horizontal Stirred Reactor with heavy-duty stirring motor in its upright position. A 4848 Controller is used to monitor the system pressure and control the temperature and stirring speed.



The 4578 Pressure Reaction System shown above and left is a 1.8 L, fixed head reactor mounted on a floor stand and rated for temperatures up to 500 $^{\circ}$ C with a maximum pressure of 5000 psig. This system is designed to operate in either a vertical or horizontal position.



The reactor shown in the two images above and right has a five gallon capacity, and is used for stirring vertically or horizontally. Pneumatic controls on the left can tilt the system upright to open the bottom drain valve. The product is filtered and collected in the lower heated sample collection vessel. This system is used for producing synthetic fabric from biomass as an alternative to a petroleum-based process.



Rotating Cage Apparatus

For Evaluation of Corrosion Inhibitors per ASTM G170 and Other Methods

valuation of corrosion inhibitors for use in refinery and oilfield applications typically requires apparatus that can simulate the interaction between the corrosive fluid and the carbon/low alloy steel used in the construction of the refinery, pipeline, or drilling equipment. ASTM G170 describes the use of a rotating cage (RC) for such testing.

Parr can provide a complete high pressure, high temperature RC system (pictured) which typically includes the following components:

- Reactor system made of Alloy C276 or other corrosion resistant alloy, with included rotating cage
- Gas feed system to deliver CO₂, H₂S, or N₂
- Liquid feed system to deliver acid or inhibitor
- Reflux condenser and pressure controller
- Automated or manual control and data logging system

Pressure and temperature ratings, reactor volume, material of construction, number of gas and liquid feeds, and control system specifics can be customized to meet individual testing requirements.

Contact Parr's Technical Sales team to discuss details.



5 Liter Vessel lowered to reveal the rotating cage. In this configuration, up to 8 coupons at up to 2,000 rpm, can be tested at one time.



Rotating Cage Apparatus Specifications			
Maximum Operating Pressure	200 bar		
Maximum Operating Temperature	200 °C		
Wetted Materials	C276		
Coupon # / geometry	8 / Variable		
Reactor Volume	5 L		
Max rotation rate	2000 rpm		

A1000G Tri-Axial Core Holder

n the oil and gas industry it is common practice to take one or more cylindrical core samples of hydrocarbon containing or other formation and to subsequently perform one or more tests with this core sample in a laboratory. During these tests the core sample is contained within a polymeric sleeve inside a pressurized sample holder assembly known as a core holder.

The operating temperature and pressure conditions are substantially similar to the hydrocarbon containing or other formation from which the sample has been taken. A treating, completion and/or stimulation fluid is injected into the sample to determine an optimum composition of the treating, completion and/or stimulation fluid and/or an optimum injection regime.

Core holders are used to facilitate a variety of studies,

including gas and liquid permeability, formation damage, porosity, carbon dioxide capture and storage (CCS), relative permeability, secondary water flooding, enhanced oil recovery and steam flooding.

A tri-axial type core holder is one where the axial force on the core sample and radial pressure on the sleeve are independent of each other and can be varied during the test. This arrangement is useful in studying the physical characteristics of the core under various conditions. Radial confining pressure is applied through the body wall, along the outer diameter of the sleeved core sample. The axial ram floats against the rock as axial pressure is applied through the end piece and the rock compresses. A specially designed spider web type groove pattern is utilized on the end faces of the core holder's distribution plugs which come in contact with end faces of the core sample. This is to ensure that the fluids, before entering or on exiting the coresample, are evenly distributed on the whole face of the core-sample.

The A1000G Tri-Axial Core Holder

With this design, the core holder does not need to be completely disassembled in order to interchange core test samples. After the confining pressure has been released the retainer and distribution plugs are removed and the core sample is easily removed from the sleeve. The sleeve and end caps remain in place in the core holder. Spacers are provided to accommodate shorter cores. The distribution plugs on either end of the core holder are provided with a single inlet / outlet.

Internal cutaway drawing of the A1000G

Parr can provide a complete hydrothermal core flooding system which typically includes the following components:

- The core holder
- One or more high pressure, low flow rate pumps
- A means to maintain the sample at elevated temperature and pressure
- Subsystems to maintain radial and axial pressure on the core sample within the core holder
- Means to monitor and record the differential pressure across the core sample
- A control system to coordinate these functional elements

A1000G Tri-Axial Core Holder Specifications		
Maximum Operating Overburden Pressure	7500 psig	
Maximum Operating Temperature	200 °C	
Core Diameter	1-in. Nominal	
Core Lengths	Up to 4-in.	
Wetted Materials 316 SS or other corrosion resistant materials		

Variations of this product are not limited to these typical specifications. If a special material or design characteristic is needed, please contact a technical representative at Parr for further information.




Supercritical Fluids



The batch supercritical water extraction vessel pictured above is constructed of Hastelloy C-276 and is designed for use to 6000 psi (410 bar) at 400 °C. supercritical fluid is any substance at a temperature and pressure above its critical point. Such fluids can diffuse through solids like a gas and dissolve materials like a liquid. Near the critical point, small changes in pressure or temperature result in large changes in density, allowing many properties of a supercritical fluid to be "finetuned". Supercritical fluids are often suitable substitutes for organic solvents in a range of industrial and laboratory processes.

Carbon dioxide is one of the many commonly used supercritical fluids. It is relatively simple to exceed its critical point (31 °C, 1057 psi). Applications that involve supercritical fluids include extractions, nanoparticle and nano-structured film formation, supercritical drying, carbon capture and storage, as well as enhanced oil recovery studies. Parr has provided systems for all the aforementioned applications.

Water is another substance that is often used in its supercritical condition (374 °C, 3185 psi). Its excellent thermal conductivity properties make it the fluid of choice in pressurized nuclear reactors for electricity generation. The extremely aggressive and reactive nature of supercritical water makes it an excellent choice for the oxidative destruction of some hazardous waste materials. This corrosive nature prohibits the use of T316 stainless steel for the vessel and requires the use of a special alloy.



Above is a schematic representation of the batch carbon dioxide extraction system shown at the right.



This batch carbon dioxide extraction system modelled above includes a 1.2 L extraction vessel rated to 5000 psi (345 bar), both CO_2 and co-solvent pumps, automatic pressure control and collection vessel with manual ball valve for extract collection, as well as custom control and data logging system with pre-programmed recipe list.

GTL (Gas-To-Liquid) System

Parr GTL System

This system incorporates three tubular reactors that can be configured as required to operate in a strictly parallel fashion or in a cascade arrangement where the products from one reactor are immediately directed to a second reactor. This type of system can support reaction schemes including but not limited to the Fisher-Tropsch process, methanation reactions, steam reforming and other similar processes.

The Fisher-Tropsch process converts carbon monoxide and hydrogen into oils or fuels that can substitute for petroleum products. The reaction uses a catalyst based on iron or cobalt and is fueled by the partial oxidation of coal or wood-based materials such as ethanol, methanol, or syngas. This reaction scheme offers a promising route to producing economical renewable transportation fuels. By carefully controlling the temperature and oxygen content, resulting products can range from syngas to "green diesel". One of the unique features of this system is a gas blending subsystem capable of mixing up to four reactant gases followed by a controlled delivery of this blended mixture to each of the three reactors via dedicated mass flow controllers.

Downstream components for each reactor include a heat exchanger/condenser, a gas/liquid separator (product receiver) and a fully automated back pressure regulator. The system includes support for introducing liquid reactants via a high pressure metering pump. The system comes completely automated with the addition of the highly versatile 4871 Process Controller (not pictured).

Other Specialty and Custom Reactor Systems




Apparatus for Vapor Pressure Determination


Vapor pressure (also known as equilibrium vapor pressure or saturation vapor pressure), is the pressure of a vapor in equilibrium with its non-vapor phases in a closed container. At a given temperature, the pressure at which a pure liquid boils is called its vapor pressure. Vapor pressures are dependent only on temperature. The vapor pressure of a pure liquid does not depend on the amount of the liquid in the container. Vapor pressures have an exponential relationship with temperature and always increase with temperature.

Investigation on vapor-liquid equilibrium of solutions can be derived from absolute vapor pressure measurements. Vapor-liquid equilibrium (VLE) data are essential for development and design of chemical processes. In the static total pressure method, the total composition in the vessel is determined by loading precise volumes of the components. The measurements are isothermal and the equilibrium pressure is measured. The amount of components in the liquid phase and the vapor phase is calculated iteratively. The static total pressure method is a fast and efficient method for determining VLE for binary mixtures with the additional advantage that the apparatus using the static total pressure method can be easily automated.

Pictured above is a custom dual vessel system for the accurate determination of the saturation pressure of specialty organic materials as a function of temperature. All the associated head fittings, valves and plumbing, including the pressure transducers, are enclosed in a heated, temperature controlled, aluminum block to minimize temperature gradients and to ensure that the average system temperature closely matches the temperature of the vessels, including their screw cap closures, housed in the lower temperature controlled block. The photo shows the exposed valves and fittings which are covered by an aluminum plate during normal operation of the system.



CAD Model for the Vapor Pressure CAD Model for the Vapor Pressure Apparatus





Disbonding Apparatus for ASTM G146

Parr Instrument Company has developed an automated reactor system for the evaluation of disbonding of bimetallic stainless alloy/steel plate. This hydrogen disbonding reactor system is capable of attaining the conditions necessary for treating samples consisting of bimetallic plate typically used in the refining industry for High Pressure/High Temperature hydrogen gas service. The test temperature and hydrogen gas pressure are selected to simulate those conditions found in refinery hydrogencontaining environments. These typically range from 14 to 20 MPa hydrogen gas pressure and temperature from 300 to 500 °C depending on actual refinery service conditions under consideration. The purpose of this test scheme is to allow for comparison of data among test laboratories on the resistance of bimetallic stainless alloy/steels to hydrogen induced disbonding (HID). This test procedure provides an indication of the resistance or susceptibility, or both, to HID of a metallurgically bonded stainless alloy surface layer on a steel substrate due to exposure to hydrogencontaining gaseous environments under HP/HT conditions. This system can be used over a broad range of pressures, temperatures, cooling rates, and gaseous hydrogen environments where HID could be a significant problem. These tests can be used to assess the effects of material composition, processing methods, fabrication techniques, and heat treatment as well as the effects of hydrogen partial pressure, service temperature, and cooling rate. This testing regime is fully described in ASTM Standard Practice G 146.

The apparatus developed by Parr consists of an HP/HT test cell made from Inconel 625. Two gas booster pumps, one for nitrogen and the other for hydrogen, are used for pressurizing the test cell. The system incorporates an automated back pressure regulator and a variety of automated valves for pressure control. The control system for the apparatus automatically purges the vessel with argon, performs a high pressure leak test with nitrogen and then sequences the main portion of the test. The main test consists of filling the test cell with hydrogen, heating the vessel at a user defined temperature set point for a period of typically 48 hours. After this, the vessel is cooled at a predefined rate, typically 150 °C/hr, until the inside temperature reached 200 °C. The test cell is then automatically vented, purged and cooled further so that the test samples can be safely retrieved. All important test parameters are user adjustable.



Apparatus for Corrosion Studies

This custom vessel was designed specifically to be used for long term (up to 4000 hours) corrosion testing of zirconium alloy samples in aqueous solutions near the critical point of water.


The vessel illustrated is a nonstirred, fixed head design with a volume of 3.7 liters with a maximum working temperature and pressure of 450 °C and 275 bar, respectively.

The vessel is heated with three clamp-on style band heaters. The system is capable of maintaining a uniform temperature of ± 1 °C over a working zone of 30 cm. Valves are provided on the head of the vessel for purging with inert gas prior the start the test as well as periodic sampling of the contents of the vessel. The controller is equipped to record temperature and pressure data for the duration of the test.

For more information on any of these custom systems please visit <u>www.parrinst.com</u> or contact the Parr Technical Service Department.



Notes



Reactor Controllers Chapter 6



Inside this chapter you will find: **Controller Overview** 4848 Reactor Controller 4848B Expanded Reactor Controller 4848M Master Controller **4848A** Reactor Controller for AC Motors 4848T Touchscreen Controller **4838** Temperature Controller A3504HC Specview Software **4871** Process Controller 4875 Power Controller 4876 Power Controller 4877 Valve Controller A2110E Motor Controller A2200E Flow Controller Interface

Parr Controller Overview

The Model 4848 Reactor Controller is our general-purpose reactor controller. It can control temperature and stirring speed, and it can be equipped to monitor a redundant temperature and pressure. It can datalog and be operated remotely from a PC.







The Model 4848B Reactor Controller is an expanded reactor controller. It has the same capabilities as the 4848 but with a larger chassis. The 4848B can accommodate up to six meters, plus the Primary Temperature Module.



The Model 4838 Reactor Controller is offered to control temperature in non-stirred pressure vessels. It can be expanded to include an optional pressure or redundant temperature module.

The Model 4871 Process Controller is a full featured process controller, which can handle either a single reactor with a wide variety of inputs and outputs or multiple reactors running independently. It controls the entire process including gas and liquid flows through a PC interface custom tailored to the application.



The Model 4848T Reactor Controller is a touchscreen controller for use with a single reactor system. It has the same control capabilities as a fully-loaded 4848 with Motor Control Module, Pressure Display Module, and secondary temperature module, with touchscreen interface and VNC communication capability.

Parr Instrument Company offers a full line of reactor controllers to monitor, control, datalog, and archive test data.

The Model 4848M is a controller for Master/Slave setups. It is most often used with large reactors with multi-zone heaters where each zone draws enough power to require its own power supply.

The Model 4848A is a reactor controller for use with AC Motors.

The A2110E is a Motor Controller. It can be provided with or without a tachometer display.

The A2200E is a Mass Flow Controller Interface. It can handle up to four Mass Flow Controllers at one time.

The 4875 Power Controller is utilized with a 4871 Process Controller or a 4848M Master Controller.

The 4876 Power Controller is utilized with a 4871 Process Controller or 4848M Master Controller. It can handle up to three heater outputs and is only available in 230V.

The 4877 Valve Controller is utilized with a 4871 Process Controller to actuate valves.

The A3504HC SpecView Software is enhanced software for use with the 4848 or 4838 Controllers. It provides superior datalogging and has an excellent user interface.

	4848	4848B	4848T	4838	4871
Application:	One Stirred Reactor	One Stirred Reactor	One Stirred Reactor	One Non-Stirred Pressure Vessel (No Motor Control or cooling)	Advanced Systems*
Features:					
Temperature Control	One	One	One	One	Multiple
Pressure Monitoring	Optional	Optional	One	Optional	Multiple
Motor Speed Control	One	One	One	None	Multiple
Auxiliary I/O	None	Up to two inputs	None	None	Multiple
Redundant Temperature Monitoring	Optional	Optional	One	Optional	Multiple
PC User Interface with Datalogging	Optional	Optional	VNC	Optional	Included
Digital Communications	RS-485	RS-485	Ethernet	RS-485	Ethernet
Expansion Modules	Up to Three	Up to Six	Three Included	Up to One	Unlimited
Number of Reactors Controlled	One	One	One	One	Multiple

* Contact a Parr Customer Service Representative for your custom application needs.

4848 Reactor Controller



4848 Reactor Controller shown with PTM, MCM, PDM, and HTM Modules installed.

The 4848 Reactor Controller

brings digital communications to all of the functions of this modular reactor controller. The 4848 offers all of the features expected in a Parr general purpose reactor controller, namely:

- PID programming with Autotuning capability for precise temperature control and minimum overshoot
- Ramp and soak programming
- Separate heating and cooling control loops
- Optional Solenoid Valve Module for cooling control
- Motor speed control
- High or low power heater switch
- Lockout relay and reset for over temperature protection
- Optional expansion modules for tachometer, pressure, and high temperature alarm

With the Parr 4848 Reactor Controller, all of the expansion modules as well as the primary temperature control module are equipped with bi-directional digital communications [RS-485] that enable the user to not only log all current readings to a PC, but also to send set points, stirrer speeds, and alarm values from the PC to the 4848 Controller.

Modular Design

A total of seven different modules are offered for the 4848 Controller. A maximum of three expansion modules in addition to the primary temperature control can be installed in the 4848 Controller. The user can select either the Tachometer Display Module or the Motor Control Module. This either/or option also applies to the High Temperature Cut Off and External Temperature Limit Modules. The 4848B has a larger chassis to accommodate up to six expansion modules.

1. Primary Temperature Control Module (PTM)

The temperature control module can accept either thermocouples or RTD temperature sensors. It has three outputs that are used for heating and cooling control and for alarm-actuated heater cut off. The control function is a full proportional, integral and derivative (PID) control with auto-tune capabilities. The controller provides ramp and soak programming with up to 49 segments.

2. Pressure Display Module (PDM)

This pressure monitoring module is set up to accept its input from a pressure transducer mounted on the reactor or attached accessory. It can be set to accept a wide variety of operating ranges. Operating pressures are transmitted continuously to the PC. These modules are available calibrated in either psi, bar, or mPa. The output from the pressure monitoring module is connected to the alarm relay to shut off power to the heater if the high pressure limit set by the operator is reached during operation.

3. Tachometer Display Module (TDM)

In this configuration, the module will display the stirrer speed and will continuously transmit it to the PC for display and logging. The stirrer speed is set manually using a potentiometer on the face of the 4848 Controller.

4. Motor Control Module (MCM)

In this configuration, the module provides true closed loop feedback control of the reactor stirring speed. The primary output of this module is wired to dynamically adjust the motor voltage in response to changes in motor loading. This provides better reactor stirring speed regulation than the standard open loop speed control, especially with reactions that involve changing viscosities. Additionally, the use of this module allows the stirring speed set point to be adjusted remotely from the host PC. A by-product of this closed loop speed control scheme is that the value of the primary controller output directly reflects the degree of loading on the motor in order to maintain a constant stirring speed. While not a direct torque measurement, this is a useful option for those who want to monitor the progress of a reaction where there is a change in viscosity as the reaction proceeds.

5. Motor Torque Module (MTM)

The MTM will display motor output from an MCM. It is particularly useful for applications with changing viscosities.

6. High Temperature Cut Off Module (HTM)

The high temperature cut off module or limit controller augments the operation of the main control module. Its redundant sensor can be mounted either internally or externally to the reactor. The primary output of the module is wired to activate the lockout relay in order to provide safety shutdown should the reactor reach an unsafe temperature.



4848 back panel for 115V model.

7. External Temperature Limit Module (ETLM)

This configuration uses the same aforementioned HTM Module with its sensor mounted in such a way to monitor the reactor's outside wall temperature. The primary output of this module is used to limit the external temperature of the reactor. This is done by interrupting the control signal from the main temperature controller when the external temperature exceeds a predetermined value. The secondary output of this module is used to activate the lockout relay in a non-latching manner if the outside wall temperature exceeds a preset unsafe temperature. The use of this module provides an effective alternative to cascade control, offering improved temperature

regulation in systems with large thermal lags, such as those found in non-stirred reactors or systems that use PTFE liners, as well as systems where the reactants have low heat capacities, such as gas phase reactions.

8. Solenoid Valve Module (SVM)

This package includes a solenoid valve and a flow adjustment valve with all of the parts required to assemble an automatic system to control the flow of coolant through a cooling coil in any reactor. It plugs into the cooling output socket on the 4848. It is designed for use with tap water as the cooling media.

4848 Continued next page >

4848 Reactor Controller



The 4848B Reactor Controller is a larger version of the 4848 with up to six expansion modules

< Continued from previous page

4848B Expanded Reactor Controller

The 4848B Expanded Reactor Controller is a larger version of the 4848 Reactor Controller. It has the same Primary Temperature Control Module as the 4848, and can accommodate up to six expansion modules. It will also communicate with the A3504HC SpecView software.

The 4848B is most useful when the process requires more modules than a single 4848 can handle. A common application would be a stirred reactor which requires a pressure display showing internal pressure and the pressure of an external gas burette. The 4848B has a chassis with enough space to handle both burettes and associated cables and other infrastructure.

The 4848B also has two optional Auxiliary modules. These will read an analog signal, most commonly 4-20 mA or 0-5 VDC. These can be fed into the chassis through a dedicated auxiliary input slot.

4848M Master Controller

The 4848M Master Controller is typically used with multi-zone heaters such as those used on our 4555 Series Stirred Reactors. It will interface with one or more slave boxes, allowing a dedicated power supply to be used for each zone.

4848A Reactor Controller for AC motors

The 4848A can be used with AC Motors. DC motors are more or less ubiquitous in Parr equipment, but some specialty applications require AC motors. The 4848A allows a tachometer display or motor control module to be used with an AC motor.

4848 Ordering Guide

A composite identification number to be used when ordering a 4848 Reactor Controller can be developed by combining individual symbols from the separate sections.

Example: A 4848 Reactor Controller, 115V electrical, with Tachometer Display Module, Solenoid Valve Module, RS-485 to USB Cable, and SpecView Package would be listed as:

No. 4848-EB-TDM-SVM-A1925E4-A3504HC

			\leq	$\leq =$	\leq
Α.	В.	С.	C .	С.	C .
Model	Voltage	Options	Options	Options	Options
4848	-EB	-TDM	-SVM	-A1925E4	-A3504HC

A Base Model

PID, Ramp & soak digital communications with motor speed control and software

Model No.	Description
4848	Reactor Controller for use with up to three additional display modules
4848B	Reactor Controller for use with up to six additional display modules
4848M	Master Controller
4848A	Reactor Controller for AC Motors

B Electrical Supply			
-EB	115 VAC		
-EE	230 VAC		

C Options	
-TDM	Tachometer Display Module
-MCM	Tachometer w/Motor Control Module
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM*	Motor Torque Module
-SVM	Solenoid Valve Module (for Cooling Control)
-AUX	0-5 VDC, 4-20 MA (4848B Only)
-A1925E4	RS-485 to USB Cable for 4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E**	RS-485 Daisy Chain for Multiple Controllers, 10-ft
-A3504HC	SpecView Software Package

*The MTM module must be installed in conjunction with the MCM.

**Must be used with A1925E6.

Dimensions				
Model	Width, in.	Height, in.	Depth, in.	
4848	11.1	9.7	11.3	
4848B	13.6	9.6	11.3	
4848M	11.1	9.7	11.3	
4848A	13.6	9.6	11.3	

4848T Touchscreen Reactor Controller

The Parr Model 4848T

Controller brings a new touchscreen design to our controller line. This full featured controller has the same functionality as a fully populated 4848 Reactor Controller but with a touchscreen interface.

Connectivity

The Model 4848T is designed for superior connectivity via network communications. After connecting to your network, the Model 4848T can communicate with your tablet, phone, or PC via a built-in web or VNC server.

Module Features of 4848T

The Model 4848T Controller is offered as a fully populated controller with Primary Temperature Control, Motor Control Module, Pressure Display Module, and secondary temperature module.

Primary Temperature Control Module (PTM)

The temperature control module can accept either a thermocouple or RTD temperature sensor. It has three outputs that are used for heating and cooling control and for alarm-actuated heater cut off. The control function is a full proportional, integral and derivative (PID) control with auto-tune capabilities. The controller provides ramp and soak programming with up to 64 steps.

Pressure Display Module (PDM)

This pressure monitoring module is set up to accept its input from a pressure transducer mounted on the reactor or attached accessory. It can be set to accept a wide variety of operating ranges. Operating pressures are displayed



4848T Reactor Controller

continuously. These modules are available calibrated in either psi, bar, or mPa. The output from the pressure monitoring module is connected to the alarm relay to shut off power to the heater if the high pressure limit set by the operator is reached during operation.

Motor Control Module (MCM)

In this configuration, the module provides true closed loop feedback control of the reactor stirring speed. The primary output of this module is wired to dynamically adjust the motor voltage in response to changes in motor loading. This provides better reactor stirring speed regulation than the standard open loop speed control, especially with reactions that involve changing viscosities. A by-product of this closed loop speed control scheme is that the value of the primary controller output directly reflects the degree of loading on the motor in order to maintain a constant stirring speed. While not a direct torque measurement, this is a useful option for those who want to monitor the progress of a reaction where there is a change in viscosity as the reaction proceeds.

High Temperature Cut Off Module (HTM)

The high temperature cut off module or limit controller augments the operation of the main control module. Its redundant sensor can be mounted either internally or externally to the reactor. The primary output of the module is wired to activate the lockout relay in order to provide safety shutdown should the reactor reach an unsafe temperature.

4848T Continued next page >



4848T Touchscreen Reactor Controller



4848T back panel for 115V model.

< Continued from previous page

External Temperature Limit Module (ETLM)

This configuration uses the same aforementioned HTM Module with its sensor mounted in such a way to monitor the reactor's outside wall temperature. The primary output of this module is used to limit the external temperature of the reactor. This is done by interrupting the control signal form the main temperature controller when the external temperature exceeds a predetermined value. The secondary output of this module

is used to activate the lockout relay in a non-latching manner if the outside wall temperature exceeds a preset unsafe temperature. The use of this module provides an effective alternative to cascade control, offering improved temperature regulation in systems with large thermal lags, such as those found in non-stirred reactors or systems that use PTFE liners, as well as systems where the reactants have low heat capacities, such as gas phase reactions.

	CONTROL	ATEA
Suprem	FC = 22.7	88-11.4
Nia Spink	Wied.]IP+0
howe	W-S	investor cuty
Aux Temp	P(+28	jasentor cely
Disponen Program Name		1
	ALTOMATED IN	CONTRACK.
		1
Program Steam		Descrive
Cycles Returnining		1
Total Yape	12	3
Colorest Imp		μ
Step Time Result		0.00
Reasoning Step To	ine .	2.00
Couges SP		31.4
No Specia		1
Prop on Name		1
Program Stemm		Esective
Cyrles Remaining	6	1
Deal Steps		μ
Cuerry Nep		0
htep Totes Raman	r	0.00
Researcing Step To	ier .	0.90
Dapar 17		10

4848T Communications via IP Network

Series 4848T Ordering Guide

A composite identification number to be used when ordering a 4848T Reactor Controller can be developed by combining individual symbols from the separate sections.

Example: A 4848T Reactor Controller, 115V electrical, with High Temperature Module, for use with 1/4 HP motor would be listed as:

No.4848T-EB-HTM-VS.25

Α.	В.	С.	D.
Model	Voltage	Options	Motor
4848T	-EB	-HTM	-VS.25

A Base Model

Touchscreen, PID, Ramp & soak digital communications with motor speed control and software

Model No.	Description			
4848T	Touchscreen Controller with Primary Temperature Control, Motor Control Module, Pressure Display Module, and either High Temperature Module or External Temperature Limit Module			
B Electrica	Supply			
-EB	115 VAC			
-EE	230 VAC			
C Seconda	ry Temperat	ure Option		
-HTM	High Tempera thermocouple	High Temperature Module with redundant		
-ETLM	External Temperature Limit Module with external thermocouple mounted through heater VAC			
Dimensions				
Model	Width, in.	Height, in.	Depth, in.	

Model	Width, in.	Height, in.	Depth, in.
4848T	11.5	12.7	13.3

Test of Vener	Va
Utile e [Invelor [at interes cape	1
States -	Sec. 1

4848T Communications via VNC

Series Number:

4838 Reactor Controller

The 4838 Reactor Controller is designed to control the temperature in our line of non-stirred pressure vessels. No provision is made for the current or future expansion to control motor stirring speeds or to actuate cooling water. The controller contains the same Primary Temperature Control Module used in the 4848 Controller and provides identical control and communications capabilities. A high/low heater switch feature, as well as a lockout relay and reset for over temperature protection, are also included. The 4838 Reactor Controllers can be enhanced with the addition of either a Pressure Display Module (PDM) or High Temperature Cut Off Module (HTM). An External Temperature Limit Module (ETLM) can be substituted for the HTM.

Features found on the 4838 Reactor Controller are:


- PID programming with Auto-tuning capability for precise temperature control and minimum overshoot
- Ramp and soak programming
- High or low power heater switch
- Lockout relay and reset for over temperature protection.

Series 4838 Ordering Guide

A composite identification number to be used when ordering a 4838 Reactor Controller can be developed by combining individual symbols from the separate sections.

Example: A 4838 Reactor Controller, 115V electrical, with optional Pressure Display Module and SpecView Software would be listed as:

No.4838-EB-PDM-A3504HC-A1925E4






Parr 4838 Reactor Controller shown with optional HTM Module

A Base Model							
For Non-Stirred	$\label{eq:Vessels} \textit{PID}, \textit{Ramp} \textit{ and } \textit{Soak}, \textit{Digital Communications and Software}$						
Model No.	Description						
4838	Reactor Controller						
B Electrica	al Supply						
-EB	115 VAC						
-EE	230 VAC						
C Options							
-PDM	Pressure Display Module						
-HTM	High Temperature Cut Off Module						
-ETLM	External Temperature Limit Module						
-A1925E4	RS-485 to USB Converter with 30-ft. Cable						
-A1925E6	RS-485 to USB Converter, isolated, with 30-ft. Cable						
-A3504HC	SpecView Software Package for 4848 and 4838						

Dimensions					
Model	Width, in.	Height, in.	Depth, in.		
4838	10.19	7.11	7.92		



SpecView[®] Software for 4848 and 4838 Controllers



SpecView Graphic User Interface for 4848 and 4838 Controllers

SpecView Software

The unified Modbus communication architecture of the 4848 allows us to bring the power and convenience of SpecView to the 4848 and 4838 Series Controllers. This version of SpecView provides support for up to eight individual loop or limit controllers.

SpecView Features and Benefits:

- · Real time charting
- Datalogging
- Alarm and event logging
- Superior alarm annunciation
- Customizable graphical user interface

This powerful program is highly customizable and configurable to show variables beyond the process value such as heater output, motor output or alarm status. Notes are easily inserted into the interface and data, and the entire interface can be rearranged and edited for the intended process.

PC Requirements:

SpecView runs on the following Windows operating systems:

Windows 10, 8, 8.1, & 7; 32 bit & 64 bit
 Note: SpecView does not work on Windows 8
 RT (as found on some Microsoft Surface tablets)

Parr recommends a minimum of:

- Intel Pentium processor or equivalent AMD
- 2GB RAM
- 1GB hard disk space for SpecView and related utilities
- Additional disk space for datalogging
- Display adapter capable of operating at 1024 x 768 pixels and at least 16-bit color
- 2 USB ports
- A CD-ROM drive

Series Number:

4871 Process Controller

The 4871 Process Controller has been developed to provide an integrated stand alone control system for controlling either a single reactor with multiple feed and product controls or multiple reactors operating independently or in parallel. The 4871 is the ideal choice for all of Parr Instrument Company's custom systems, including: tubular reactors, multiple reactor systems (both batch and cascade systems), continuous flow stirred reactors, and fluidized bed reactors.

Control Module

The control module of the Parr 4871 Process Controller is a Honeywell HC900 Hybrid Controller. This controller combines analog and logic control into a versatile, cost-effective controller designed specifically for process applications requiring analog measurement combined with programmable control actions.

Function	
1 unction	

i dilotion	
Analog Input (AI) Maximum Channels	480
Analog Output (AO) Maximum Channels	200
Digital I/O Maximum Channels	1920
Remote I/O	Yes
Control Loops	Per Memory
Logic Scan	25-50 ms
Loop Scan	500 ms
Ethernet Communications	Yes
Peer to Peer Communications	Yes
Modbus Master	Yes
Modbus Slave	Yes
On-line Programming	Yes
Function Blocks	2000
Sequential Functions	Yes
E-mail Alarms	Yes

Input / Output


The controller is adapted to each user's requirements by adding to the control chassis input and output modules. Each module provides for between four and sixteen individual inputs or outputs fully isolated from one another.

These modules include:

Input Modules

The analog inputs are of universal type and are most commonly used for thermocouple or RTD temperature sensors, strain gage type pressure transducers, and similar devices with mV, V, or resistance inputs. Input isolation, cold junction compensation, and burnout protection are incorporated into the circuitry. Each analog input module provides for eight separate inputs.







Process Flow Diagram for stirred reactor system.

The digital inputs can be logic inputs or contact closures. These are typically used for sensing valve positions or conditions of safety devices. Each digital input module provides for sixteen separate inputs.

Continued next page >



4871 Process Controller



4871 Process Controller Input/Output Diagram.

< Continued from previous page

Output Modules

The analog outputs are 0-20 mA. A suitable dropping resistor can be used to convert this to 0-5 or 0-10 VDC. Analog outputs are commonly used to set stirrer motor operating speeds, position control valves, or drive mass flow controllers or pumps. Each analog output module is capable of controlling four separate devices or functions.

The digital outputs are open collector type capable of sinking up to 300 mA. They are commonly used to control heaters, solenoid valves for cooling or other flow control, system safety shut down, visual and/or audible alarms, and similar devices. Each digital output module is capable of controlling sixteen separate devices or functions.

Control Loops

The controller can provide any number of PID or ON/OFF control loops, limited solely by the available CPU memory.

The PID control algorithm includes auto-tuning and fuzzy logic overshoot suppression for each control loop. For heating and cooling control, the PID control loops provide time proportioning of the associated digital output.

Many temperature control applications utilize two separate time proportioning outputs with one PID controller; one for heating and one for cooling.

Control loops can be linked together to provide cascade, feed forward or ratio control for difficult or advanced control applications.

Both high and low limit values can be entered for each control loop to sound alarms or initiate safety control schemes.

The 4871 Process Controller combines the following components into a single control system:

Control Module:


- Honeywell HC900 Hybrid Controller

Flexible, Powerful Software:


- SpecView SCADA software
- System set up with graphical user interface configured to individual requirements

Power Controller:

- 4875 or 4876 Power Controllers for handling heating, cooling, safety, and motor control devices



Process Flow Diagram for a tubular reactor with real-time process renderings.



 $\label{eq:screenshot} \ensuremath{\mathsf{Screenshot}}\xspace{0.5ex} \ensuremath{\mathsf{otstroller}}\xspace{0.5ex} \ensuremath{\mathsf$

4871 Dimensions					
Model	Width, in.	Height, in.	Depth, in.		
4871A	12.86	15.12	11.00		
4871B	16.98	21.56	13.00		
4871C	29.78	21.56	13.00		

Set Point Programming

Ramp/Soak profiles for controlling the entire process of a reactor can be written using the set point profiler incorporated into the control firmware. A single profile may be from 2 to 50 segments in length.

A typical profile might be a ramp and soak of the reactor temperature but, in addition, the analog and digital outputs can be tied to the basic profile to start and stop flows, activate stirrers or accessories, or change alarms. Any of the set points within the profile can be protected with the set point guarantee function that assures that the process variable will be within the entered limits before the profile can proceed.

The number of set point profiles is limited by the amount of available memory in the controller CPU. Typically, at least eight separate profiles can be running simultaneously. For example, eight reactors can be running a unique program at the same time.

While a maximum of 99 profiles can be stored in the controller itself, an unlimited number can be stored in the operator's PC for rapid transfer to the controller.

In addition to the set point profiling capability, the controller is also equipped with a set point scheduling function. This feature can operate up to 8 profiles operating on a common time base.

The 4871 Controller includes Operational Sequence Control

The sequence control function offered by the 4871 Controller greatly expands the capabilities of this control for users who wish to control reactor systems. The operation of valves, pumps and other peripheral devices can be programmed on either a time or an event driven basis. Sequences can be very simple timed events or they can be very complex with multiple nested default sequences programmed to occur only if process feedback indicates a need to take alternative actions.

Communications Channels

Each 4871 Controller is equipped with an RS-485 and Ethernet communication ports. The RS-485 port is used for initial configuration and communications setup. The Ethernet port provides communication with the host PC when using the SpecView GUI program. Multiple controllers, each with a unique address, can be networked on the Ethernet interface with a single connection to the PC.

The principal advantage of the Ethernet interface is that it allows the user to use an existing network infrastructure to connect the controller to the PC. As a result, one can operate the controller over the network from anywhere within your facility.

Continued on next page >



4871 Process Controller

< Continued from previous page

Additionally, internet access from remote locations becomes possible. This type of connectivity offers unique possibilities, for example, related to remote diagnostics and system troubleshooting.

Three Models Available

Parr 4871 Process Controllers are available in three different models to cover a wide range of applications. The 4871A will accept four I/O modules (typically up to 36 inputs and outputs). The 4871B will accept up to eight I/O modules, and the 4871C will accept up to twelve I/O modules. If more than 100 inputs and outputs are required, multiple controllers can be linked.

SpecView SCADA Software

SpecView describes their product as "Software for people with other jobs". That seems to be an excellent description of this software package used with the 4871 Controller to:

- Interface with the control package
- Develop the graphical screen layout
- Establish the data logging profiles
- Prepare custom reports
- Create bar graphs
- Generate time trend graphs
- Monitor alarms
- Create flexible recipes
- Retrieve and replace logged data
- Operate the reactor system(s).

The full software package, not just a run-time version, is supplied with the 4871 Controller so operators can enhance their system as their needs change or expand/change their applications. Download a demo at www.specview.com.


Graphical User Interface

An integral part of the Honeywell Controller is the hybrid control designer software. This is the "Drag and Drop" software that enables Parr to rapidly establish the controller's internal logic and adapt it to individual systems requirements.

The user can employ this same software to change or enhance the fundamental logic of the controller as additional components are added to the system or as functions need to change.

Current Industrial Standards

Modbus, Ethernet, auto-tuning, fuzzy logic, auto-configuring Man Machine Interface (MMI), supervisory control and data acquisition (SCADA), multi-loop control; these are a few of the terms and capabilities designed into the hardware and software incorporated into the Parr 4871 controller.



Screenshot from remote PC showing logged data in Data Logging Mode.

This modern and powerful package enables Parr to configure and build turnkey systems within weeks of order at very attractive prices.

PC Requirements

In most laboratory and pilot plant applications, a PC will be used for the operator interface. A touch screen interface is also available. For plant or production applications, an industrial type user interface box with a color graphic LCD is available.

The PC used with the 4871 serves several functions. It is the operator's user interface for controlling the process. It also logs all of the operating data generated during a process run. It can also store multiple setups for rapid transfer to the 4871 Process CPU.

Any modern PC with current Windows operating system can be used with these controllers.

It is important to note that the control of the process always resides within the 4871 Controller, and not in the PC itself.

4871 Ordering Guide

Each 4871 Controller will be specified and assembled to match the user's intended application. These applications can vary from a single reactor with full process control to eight reactors operating in parallel. The 4871 Process Controller connects to a customer-supplied PC with a current Windows Operating System.

Please contact a Parr Customer Service Representative to discuss configuring a 4871 Process Controller to your specific application.

Please contact a Parr Customer Service Representative to discuss configuring a 4871 Process Controller to your specific application.

4871 Accessories



The 4875 Power Controller module is flexible and can be used in remote locations from the system being controlled.



4875 Power Controller back panel for 230V model.

4875 Power Controller

The power controller(s) will be designed for each individual system, but as an example, a 4875 Power Controller includes:

- 1. A solid-state relay sized to handle the current drawn by the vessel heater. This is commonly a 25-amp relay with its protective fuses.
- 2. A solid-state relay sized to drive a solenoid valve to control the flow of cooling water to the vessel.
- 3. A motor speed controller that converts the analog output signal from the controller to the electrical signal required to drive the specific stirrer motor. A circuit breaker for the motor is also provided.
- 4. A lockout relay to shut down the heater circuit should an alarm condition be detected.
- 5. Status lights for the principal functions.
- 6. Connections of appropriate style for the power input and device outputs.

4875 Dimensions						
Model	Width, in.	Height, in.	Depth, in.			
4875	9.40	5.00	10.50			

Power Controllers

Parr designs, builds and furnishes power controllers to adapt the analog and digital outputs from the 4871 Controllers to the reactors or systems being controlled. These power controllers handle all of the high current power circuits so that the control circuitry is isolated from these loads. This also makes it possible to install the controller in a control room some distance from the system being controlled.



4876 Power Controller



Back panel for 4876 Power Controller

4876 Power Controller

The 4876 Power Controller is equipped with three solid state relays, three lockout relays, and a single motor control. It is used with the 4871 to control up to three separate heaters and one motor. It can only be used in 230V systems.

4876A Power Controller

The 4876A Power Controller is the same as the 4876 described above but without motor control capability.

4876 Dimensions						
Model	Width, in.	Height, in.	Depth, in.			
4876	16.18	4.98	13.00			
4876A	9.40	5.00	10.50			



Additional Controllers

4877 Valve Controller

The 4877 Valve Controller is used with the 4871 Process Controller to remotely actuate valves. It requires a dry air supply.

4877 Dimensions							
Model	Width, in.	Height, in.	Depth, in.				
4877 (1-4 valves)	15.00	5.00	16.00				
4877 (5-8 valves)	19.00	5.00	16.00				
4877 (9-12 valves)	24.00	5.00	16.00				

Series 5000 Multi Reactor System (MRS): a six station multi-reactor system. Each reactor is equipped with it's own constant pressure gas delivery system. The 4871 Controller maintains all of the important system parameters, including temperature and stirring speed, and records the gas consumption of each

of the reactors.

11日



A2200E Flow Controller Interface

The A2200E Flow Controller Interface can be used to send a set point and receive flow rate information from up to four MFCs. It is a basic control box for use with distributed control systems. Users needing more advanced features like PC communication and totalization should consider our 4871 Process Controller.

A2200E Dimensions							
Model	Width, in.	Height, in.	Depth, in.				
A2200E	9.40	5.00	10.50				



A2110E Motor Controller

The A2110E Motor Controller is a stand alone box which can be used for a single motor. It may be equipped with a Motor Control Module expansion for setpoint control of RPM. A version of this controller with 4-20mA I/O capability for integration into existing systems is available.

A2110E Dimensions						
Model	Width, in.	Height, in.	Depth, in.			
A2110E	9.00	4.64	13.00			



This custom order is set up to run sixteen reactors, two 4871 Process Controllers, with sixteen 4875 Power Controllers all through one PC.

Optional Accessories Chapter 7



Heaters

Standard Heaters

Parr has designed standard electrical heaters for all of the reactors in our product line. Different types of heaters are used for individual reactors to best meet the operational needs, heating load, and expected operating temperatures. The standard heater type and power rating for each reactor model is listed in the reactor specification tables.



Ceramic Heaters

These are special purpose heaters with an electric element embedded in a shaped ceramic body which is held within an insulated metal housing. They are used for reactors designed for temperatures to 600 °C and for large multi-zone heaters.



Calrod-Type Sheathed Element Heaters

These are rugged heaters with Calrod-type elements held within a metal shell. They are used for medium to large reactors for operating temperatures to 350 °C. In some cases the heater shell itself forms a part of the reactor support. An advantage of Calrod heaters is that the heating elements are easily replaceable.



Rigid Heating Mantles These are quartz fabric mantles housed in aluminum shells. They are used for moderate sized reactors in designs where the heater can be moved on or off the vessel. They are light weight and easy to handle, but they are not used to support the weight of the vessel and they are generally limited to operating temperatures of 350 °C or less.

Optional and Custom Heaters

Parr offers a variety of heater designs which can be substituted for the standard heater normally furnished with each reactor. Most of these can also be used with Parr non-stirred pressure vessels as well.



Clamp-On Band Heaters These are normally used for very small reactors where maximum watt densities and heat transfer are required due to the limited surface area available on the vessel.



Aluminum Block Heaters Aluminum block heaters are available as an option for vessels of two gallons or less. These heaters are machined from solid blocks of aluminum with heater wells machined into the walls of the block. Optional cooling channels are also available.

Aluminum block heaters have distinct features that make them desirable for certain applications. When compared to a circulating jacket, they offer the convenience of direct electric heat control as well as no plumbing requirements for hot oil. The heating elements, cooling channels, and associated wiring are embedded within the machined aluminum enclosure. This allows for even heat distribution, physical isolation of the electrical connections, and rapid cooling when needed. These attributes, in combination with a surface temperature

limiting device and proper external wiring, have allowed these heaters to be used in some potentially ignitable atmospheres (Hazardous Locations). Having integral machined cooling channels, aluminum block heaters have also been used for control of some exothermic reactions when internal cooling is not available.



Flexible Heating Mantles

These can be furnished for many different applications. These are similar to our rigid type heating mantles except they are not held in an aluminum housing. They have a flexible fabric outer case for electrical and thermal insulation. This type of mantle is particularly useful for heating vessels with irregular shapes, such as those with windows in the cylinder wall, since they are flexible and can be split and laced onto a vessel around any external protrusion. As with rigid mantles, they will produce temperatures up to 350 °C, but they are limited to watt densities of 10 watts per square inch. This type of heater can be made to cover any of the vessels offered by Parr, and they are sometimes preferred when only moderate temperatures are required. Since they are constructed of cloth, an electrical ground wire cannot be provided.





Circulation Jackets

A jacket can be welded to the outer wall of most Parr pressure vessels to provide a means for heating or cooling the vessel with a hot or cold liquid or steam. This type of heating is ideal for users who want to duplicate plant operating conditions, using a jacketed reactor comparable to jacketed equipment used in their plant. Since there are no electrical components in a jacket, and since the maximum temperature can be controlled by controlling the temperature of the heating medium, a jacketed vessel may be a good option for use in hazardous atmospheres.

Rapid and uniform heating can be attained with a jacketed vessel since the heating medium is in direct contact with the vessel. By controlling the temperature of the heating medium, temperature overshoots can be avoided when working with sensitive materials. Standard jackets are designed for operating pressures up to 100 psig (7 bar) within the jacket. Higher pressure jackets can be provided if required.

117

Stirrer Motors and Drives

Torque vs. Stirring Speed

The standard, open-type, variable speed motor installed on each Parr reactor will produce stirring speeds from zero to between 600 and 800 rpm with a torque adequate to drive the installed impellers in average viscosity mixtures. Higher horsepower motors and special stirrers can be provided for higher viscosities. Alternate drive pulleys are available to produce higher stirring speeds, but several basic rules must be considered when changing any of these components.

The highest torque from any motor is obtained at lower stirring speeds. Increasing the stirring speed reduces the torque in inverse proportion to the speed. For operations involving high viscosity mixtures, the motor size, the type of impeller and the stirring speed must be matched to provide an effective mixing system.

As a general rule, the magnetic coupling installed on each Parr reactor will have a torque rating considerably higher than the torque obtainable from the motors offered for use with that apparatus. The goal is to make the motor the weak link so that the magnetic stirrer will be protected. Reference torque rating for applicable magnetic drive.

Explosion-proof Motors

Explosion-proof motors designed for Class I, Groups C and D and Class II, Groups F and G with variable speed control can be furnished for most Parr reactors.

Flameproof "d" Motors

ATEX certified AC Flameproof Motors designed for use with group IIC gases are available for all of Parr's stirred reactors. These inverter duty motors are available with simple variable frequency drives for control of the stirring speed. Please consult Parr for any specific ATEX certification related requirements.

Air Motors

Air-driven motors can be installed on most reactors. The horsepower rating, torque, and available speed are all dependent upon the pressure and available volume of the driving air source. Maximum torques are delivered at relatively slow speeds and maximum horsepower is delivered at high speed.

			Standard Pul	ley	Optional Pull	ey		
Motor Designation	HP (kW) Rating	Hazardous Location*	Variable Speed	Туре	Max Speed, RPM**	Max Torque, in-lb	Max Speed, RPM**	Max Torque, in-lb
-VS.12	1/16 (0.05)	No	Yes	PMDC	600	6.75	1700	2.25
-VS.25	1/4 (0.18)	No	Yes	PMDC	600	27	1700	9
-XP.25	1/4 (0.18)	Class I, Div. 1 & 2 Groups C & D, E & F	Yes	PMDC	600	27	1700	9
-XP.25X	1/4 (0.18)	Ex 2G de IIC T4	Yes	Inverter Duty Constant Torque AC	450	33	1350	11
-AM.25**	1/4 (0.18)***	Ex II 2 GD c T4	Yes	Air	1000***	30	—	—
-VS.50	1/2 (0.37)	No	Yes	PMDC	600	54	1700	18
-XP.50	1/2 (0.37)	Class I, Div. 1 & 2 Groups C & D, E & F	Yes	PMDC	600	54	1700	18
-XP.50X	1/2 (0.37)	Ex 2G Ex de IIC T4	Yes	Inverter Duty Constant Torque AC	450	66	1350	22
-AM.50**	1/2 (0.37)***	Ex II 2 GD c T4	Yes	Air	1000***	66***	—	—
-VS.75	3/4 (0.55)	No	Yes	PMDC	600	81	1700	27
-XP.75	3/4 (0.55)	Class I, Div. 1 & 2 Groups C & D, E & F	Yes	Yes	600	81	1700	27

Values represented are nominal.

Stirrer Drive Motors

*For more information on Hazardous Locations information see Tech Note 230.

**Maximum speed values based on "no load"

***HP, RPM, and torque values for air motors are based on a 40 psi supply capable of 34 cfm for the AM.50 and 10 cfm for the AM.25.

118 Parr Instrument Company

Gear Box Torques							
	3:1 Gear Box 5.				10:1 Gear Box		
Motor HP Rating	Max Speed, RPM	Max Torque, in-Ib	Max Speed, RPM	Max Torque, in-lb	Max Speed, RPM	Max Torque, in-lb	
1/4	600	27	360	45	180	90	
1/2	600	54	360	90	180	180	
3/4	600	81	360	135	Not Recommended		

Contact Parr for ATEX Gear Box.

Geared, Direct Drive Motors


A geared, direct drive motor can be installed on most fixed head floor stand reactors. This is an attractive arrangement for handling heavy stirring loads.

Any 1/4 hp or larger, variable-speed standard or explosion-proof motor can be used. Gear box drives are available with ratios of 3:1, 5:1 and 10:1. The 1700 rpm maximum speed will be reduced in an amount determined by the reduction ratio of the gear box, and the associated torque values from the table will be increased in the same ratio.

Magnetic Drive				
Description	Maximum Torque, in-lb			
General Purpose	16			
Footless General Purpose	16			
Heavy Duty	60			
Footless Heavy Duty	60			
Extra Heavy Duty	120			
Footless Extra Heavy Duty	120			



Stirrer Options



(5) Paddle Type Anchor

Turbine Type Impellers (1)

Parr reactors are usually equipped with turbine type impellers which produce an excellent mixing action over the range of stirring speeds at which these reactors typically operate. These impellers are made in four-blade and six-blade styles, with the smaller four-blade impellers used only on Micro and Mini Reactors. These impellers, for reactors with 300 mL volume or greater, may be positioned anywhere on the stirring shaft, with one impeller usually located near the bottom of the vessel to keep solids up in suspension and a second impeller positioned near the base of the vortex to pull reactant gases down into the liquid phase. These impellers generally provide excellent mixing for systems with effective viscosities up to approximately 25,000 centipoise (cP) with a 16 in-lb magnetic drive or up to 50,000 cP with 60 in-lb magnetic drive.

Spiral Stirrers (2,3)

Spiral stirrers are offered in two styles: stamped and machined. The standard spiral stirrer includes blades constructed of stamped sheet metal and provides a cost-effective option. The more robust machined spiral stirrer is manufactured from a solid piece of bar stock which results in added strength and facilitates cleaning. Either of these spiral stirrers can be installed in any 1 liter, 2 liter or 1 gallon reactor to produce a positive downward thrust or upward thrust action when working with viscous polymers or other high viscosity mixtures. They work best in floor stand reactors with adjustable speed and heavy duty drive systems. Either left-hand (downward thrust) or right hand (upward thrust) spirals are available. The downward thrust spiral is generally preferred for heavy suspensions.

Note: All stirrer options may not be appropriate or available for each reactor size. Additional internal fittings may be required to adapt some stirrer styles to existing reactors in the field. Please contact the Parr Technical Service Department for assistance in selecting a stirrer suitable for the intended operating volume and viscosities.

Anchor Stirrers (4-6)

Anchor stirrers are available in several configurations for use with moderate to high viscosity materials. This type of stirrer usually works best in vessels with an inside depth to diameter ratio of 1.5 to 1 or less. They are intended to operate at relatively slow speeds and generally require a heavy duty drive system capable of generating and delivering sufficient torque to the agitator. Footless magnetic drives work well with anchor or spiral stirrers.

Three basic types are offered:

- A U-shaped, flat bar anchor.
- A flat blade, paddle type anchor.
- A two-arm or three-arm, self centering anchor with PTFE wiper blades.

Gas Entrainment Impellers

Parr offers a series of gas entrainment impellers for users who want to obtain maximum gas dispersion into a liquid system. This is obtained with a unique impeller attached to a hollow stirring shaft through which gases are continuously recirculated from the head space above the liquid through the impeller into the liquid phase. As with all impellers, the speed of the stirrer creates a vacuum at the tip of the impeller. Gas enters openings near the top of the shaft and is pulled through dispersion ports located at the tips of the impellers. In the Parr system with dispersion ports located at the very tips of the impellers, the higher the stirring speed - the higher the vacuum — and the higher the driving force for this very effective gas dispersion system.

When ordered with a new reactor, these impellers are offered as a complete package which includes the impeller, the hollow shaft with coupling, and any required foot bearings and brackets for the intended reaction. The baffles are a separate option which must be ordered individually.

The gas entrainment stirrers may be ordered as an optional stirrer when purchasing a new reactor system or



easily installed in an existing system in the field. With the wide variety of reactor head styles and magnetic stirrers furnished on Parr reactors it is best to contact us with the numbers stamped on the head of your vessel so that we will be able to furnish the correct gas entrainment assembly for a particular reactor system.

Since these gas entrainment impellers operate best in the 1000-1200 rpm range, users will want to ensure that their stirrer drive system is set up to deliver these operating speeds; alternate pullevs and belts are available to convert existing reactor systems.

Baffles

Because it is the relative speed of the tip of the impeller to the liquid phase that governs the mass transfer, baffles, which impede the rotation of the liquid with the impeller, can greatly enhance the operation of these gas entrainment impellers. While some natural baffling is provided by the internal thermowell, dip tube and cooling coils, the removable baffles are recommended for use with these gas entrainment impellers. These baffles may also be beneficial with the more traditional turbine type impellers for certain applications.



Gas Entrainment Impeller with Hollow Shaft

Removable Baffle

Catalyst Baskets

Catalyst baskets can be provided for holding a supported catalyst so that it will not be destroyed or changed by the stirring action of the impeller. These can be installed in reactors with volumes ranging from 300 to 2000 mL. Two interchangeable styles are available. Special heads, internal cooling coils, thermowells and dip tubes are required to provide clear space in the vessel for these baskets.

The Static Design

In the static design the mesh basket holding the catalyst remains stationary while impellers on the stirring shaft and baffles outside of the basket direct the flow of reactants over the surface of the contained catalyst. A unique gas entrainment impeller provides a uniform flow of both gas and liquid over the fixed catalyst bed held within the annular basket. The Parr design for these baskets includes a rigid bottom support which permits high speed stirring without excessive vibration. Cooling coils, internal temperature measurements, and liquid and gas sampling operations can be continued as usual without interference from the installed catalyst basket.

The Dynamic Design

In the dynamic design the catalyst is held in an annular shaped, mesh basket which is attached to the stirrer drive in place of the stirring shaft. The rotating basket then serves as an impeller for stirring the reactants. Fixed baffles and coaxial impellers ensure good circulation over the surface of the contained catalyst. The dynamic baskets are available for reactors with volumes of 1000, 1800, and 2000 mL. Dynamic baskets must be installed in reactors equipped with at least 1/4 hp motors to ensure that sufficient stirrer torque and speeds are available for proper operation. Dynamic baskets are interchangeable with static baskets in 1 liter and larger vessels.



Catalyst Basket Dynamic Design



Catalyst Basket Static Design



Catalyst Basket Static Design with Uniflow Stirrer

Optional Accessorie

Condensers

Parr offers two styles of condensers for attachment to the head of a stirred reactor or pressure vessel. These can be made in various sizes to match the size of the reactor.

Reflux Condenser

The reflux condenser consists of a length of tubing connected directly to the head of a vessel and equipped with a water cooling jacket. Condensed vapors are returned directly to the vessel and any non-condensable gases can be released through a needle valve at the top of the condenser. A helical insert in the condenser ensures maximum effectiveness in a rather short length.

Reflux/Take-Off Condenser

The reflux/take-off condenser consists of the same water jacketed tube described above, assembled with a receiving vessel attached to the lower end of the condenser. Any vapor, such as water from a polymerization reaction, can be condensed and collected in the receiver, from which it can be withdrawn through a bottom valve. Any non-condensable gases can be released through a needle valve at the top of the condenser. If condensate collection is not required, the receiver can be removed and the condenser can be mounted directly above the reactor for direct reflux into the vessel.

Modifications

Many users opt to install a ball valve at the head of the reactor below the condenser to use as a shut-off to the condenser. Alternate collection vessel volumes are available upon request.


The installation of a condenser on any of the Parr reactors requires a larger port in the head of the vessel, the size of which will vary with the volume of the reactor system. Due to the limited space on the 4560 mini reactors, either the gage opening or one of the cooling coil ports can be enlarged to 1/4" NPT for use with a condenser. This modification would then either combine the gage and condenser functions or eliminate the internal cooling loop to accommodate the condenser. Reactors with volumes of 1 liter and greater would be modified with a 3/8" NPT opening or larger depending on the reactor volume. The standard head fittings would be rearranged to accommodate this port.



Condensers					
Reactor	Style	Note	Inner Tube O.D. Diameter, in.	Standard Receiver, mL	Part No.
4560 / 4590	Reflux/Take-off	Mod. Gage Opening 1/4" NPT	3/8	150	A2011HC
4560 / 4590	Reflux	Mod. Gage Opening 1/4" NPT	3/8		A2012HC
4560	Reflux/Take-off	Mod. Cool Coil Opening 1/4" NPT	3/8	150	A2013HC
4560	Reflux	Mod. Cool Coil Opening 1/4" NPT	3/8		A2014HC
4520 / 4530 / 4550 Reflux/Take-off	3/8" NPT	1/2	300	A2001HC	
4520 / 4530 / 4550	Reflux	3/8" NPT	1/2		A2002HC
4530 HD*	Reflux/Take-Off	3/8" NPT	1/2	300	A2003HC
4530 HD*	Reflux	3/8" NPT	1/2		A2004HC
4540 / 4570 / 4580	Reflux/Take-off	3/8" NPT	3/8	300	A2016HC
4540 / 4570 / 4580	Reflux	3/8" NPT	3/8		A2017HC
4555 Reflux/Take-off		1/2" NPT	3/4	1000	A2018HC
4555	Reflux	1/2" NPT	3/4		A2019HC

*with HD Magnetic Drive

Safety Rupture Discs



Parr Pressure Vessels are protected by custom built rupture discs. Examination of these discs will show that each of these discs is domed. This dome was produced at the factory by taking the individual disc to 70% of its burst pressure.

ASME as well as other pressure vessel codes dictate that pressure vessels must be equipped with a rupture disc designed to burst no higher than the design pressure of the vessel. For pressure loads that do not cycle rapidly such as in our vessels, we suggest limiting the actual

Typical Rupture	Typical Rupture Discs for 1/4" Orifice					
Burst Rating, psig	Inconel Disc	Gold-Faced Inconel Disc				
1000	526HCPD	581HCPD				
2000	526HCPF	581HCPF				
3000	526HCPG	581HCPG				
4000	526HCP40CT	581HCP40CT				
5000	526HCPH	581HCPH				
8000	526HCPJ	581HCPJ				
10000	526HCP100CT					

Note: For a complete list of rupture disc part numbers, burst ranges, materials, and temperature ratings, see manual 231M. operating pressure to no more than 90% of the disc burst pressure. This combination will limit operating pressures to no more than 90% of the design pressure of the vessel in most cases.

We have selected Alloy 600 as the standard material for these rupture discs. It provides excellent corrosion resistance while retaining over 90% of its room temperature rating at temperatures up to 450 °C. For added corrosion resistance we can furnish these discs with gold facing or replace them with discs made of Alloy C-276. Discs can be produced to match any operating pressure and temperature above the stated minimums.

Parr reactors and pressure vessels from 25 mL to 2000 mL use the 526HC Series Alloy 600 disc or 581HC Series Alloy 600 with gold facing. The 1 gallon and larger use the 708HC series discs. The 4580 reactor systems use the 1415HC series discs.

In general, the 1000 psi disc in the 526HC/581HC series discs and the 800 psi in the 708HC are the lowest available ranges in the Alloy 600 material. Alternate disc materials are available but they do not offer the same corrosion resistant properties and temperature capabilities.

For applications where users prefer a lower range pressure gage, we would add a spring loaded relief valve set to protect the gage and a 1000 psi rupture disc as the fail safe protection.

Safety Rupture Discs

Most rupture discs furnished by Parr come with CE certification. Upon request, Parr can furnish ½" discs with flat seat (both holder and rupture disc) with ASME certification markings and documentation (with UD certification designator).

Users are invited to contact the Parr Technical Support Staff with requirements for your rupture disc needs.

Typical Rupture Discs for 1/2" Orifice

Burst Rating, psig	Inconel Disc
1000	708HCP10CT
1500	708HCP15CT
2000	708HCP20CT
3000	708HCP30CT
3000	1415HCP30CT
4500	1415HCP45CT

Pressure Relief Valves

Spring-loaded relief valves should be viewed as a supplement to and not a substitute for a safety rupture disc. The rupture disc is the primary means for protection of the operator and the vessel in case of accidental over-pressure. Spring-loaded relief valves can be added to a reactor or vessel to:

- Relieve pressures near the maximum _ operating pressure.
- Reseal once excess pressure has been relieved. _
- Protect low pressure components at pressures below available rupture disc ranges.

The relief valves listed below can be installed on any Parr vessel. The relief valves are constructed of stainless steel and have FKM O-rings. Other valve and O-ring materials are available on special order. Relief valves with CE Certification are also available.





Relief Valve

Pressure Relief Valves for Low Pressures					
Part No.	Preset psi	Adjustable psi	Bar	Discharge Connection	
A140VB2PA	100	50-150	3.45-10.3	1/4" NPT (M)	
A140VB2P1	150	50-150	3.45-10.3	1/4" NPT (M)	
A140VB2P2	145	50-150	3.45-10.3	1/4" NPT (M)	
A140VB2PB	200	150-350	10.3-24.1	1/4" NPT (M)	
A140VB2PC	600	350-600	24.1-41.4	1/4" NPT (M)	
A140VB2PD	300	150-350	10.3-24.1	1/4" NPT (M)	

Pressure Relief Valves for High Pressures					
Part No.	Pressure Relief Range, psi	Bar	Discharge Connection		
A175VB	750-1500	52-103	1/4" NPT (F)		
A175VB2	1500-2250	103-155	1/4" NPT (F)		
A175VB3	2250-3000	155-207	1/4" NPT (F)		
A175VB4	3000-4000	207-276	1/4" NPT (F)		
A175VB5	4000-5000	276-345	1/4" NPT (F)		
A175VB7	350-750	24-52	1/4" NPT (F)		
A175VB8	50-350	3.4-24.1	1/4" NPT (F)		

* Note: When ordering any of the above relief valves, the user may specify a desired set pressure.



Gages



4-1/2" Dia.

Gages for Parr pressure vessels can be furnished with either

3-1/2" or 4-1/2" dials in any of the ranges shown in the table below. All have stainless steel Bourdon tubes and 1/4" NPT male connections.

Alloy 400 gages are available on special order. Accuracy is 1.0 percent of full scale for the 4-1/2" size and 1 percent for the 3-1/2" gages. All standard gages include dual scales, with graduations in both pounds per



2633HCP10AD 3-1/2" Back Mount Gage

square inch (psi) and bar. Gages in Pascal units are available on special order. Compound gages which show vacuum to 30 inches of Mercury and positive pressures to 300 psi (20 bar) are also available.

When ordering a special gage, specify the gage diameter, the desired range and scale units.

The gage on a pressure vessel should be 150 percent of the maximum operating pressure. This allows the gage to operate in the most accurate pressure range and prevents the gage from being stressed repeatedly to its full range, which will effect the calibration.

Gages				
Pressure, psi	Range, bar	4-1/2" Dia. Gage No.	3-1/2" Dia. Gage No.	3-1/2" Dia. Back Mount Gage No.
0-100	0-7	56HCPA	593HCP1AD	2633HCP1AD
0-200	0-14	56HCPB	593HCP2AD	2633HCP2AD
0-600	0-40	56HCPC	593HCP6AD	2633HCP6AD
0-1000	0-69	56HCPD	593HCPD	2633HCP10AD
0-2000	0-138	56HCPF	593HCPF	2633HCP20AD
0-3000	0-207	56HCG	593HCPG	2633HCP30AD
0-4000	0-275	N/A	593HCP40AD	N/A
0-5000	0-345	56HCPH	593HCP50AD	2633HCP50AD
0-7500	0-517	56HCP75AD	N/A	2633HCP75AD
0-10000	0-690	56HCPK	593HCPK	N/A
30" Hg Vac/300 psi		56HCP3YB	593HCP3YB	2633HCP3YB

Parr's standard stainless steel reactor gages are typically manufactured to meet ASME/ANSI B40.1 Grade A accuracy specifications.



Gas Measurement Systems

Parr offers a variety of accessories for its line of pressure reaction vessels to enable the investigator to accurately determine the amount of gas consumed in a reaction conducted at elevated pressures and temperatures. There are essentially two methods used to measure the amount of gas delivered to a reaction vessel. These are:

- The measurement of the pressure drop in an auxiliary supply vessel of known volume.
- 2. The measurement and integration of the flow rates using an electronic mass flow meter.

Each of these methods has its advantages and limitations as discussed below.

Intermediate Supply Tanks

Certainly the simplest method to measure the amount of gas consumed in a reaction is to feed the gas from a vessel of known volume and to measure the pressure drop in this vessel during the course of the reaction. The consideration in this method is to select a supply vessel with a volume matched to the amount of gas that will be consumed in the reaction. It needs to be large enough to contain enough gas to complete the reaction and small enough that the pressure drop will be significant and measurable. This basic technique can be applied in a number of ways:

- The supply tank can be connected directly to the reaction vessel. This is the simplest and least expensive. The principal limitation of this approach is that the reaction pressure will fall as gas is consumed and the reaction will not be conducted at a constant pressure.
- 2. The supply tank can be fitted with a constant pressure regulator. The regulator must be selected to match the planned operating pressure. This regulator will deliver gas to the reaction vessel at constant pressure overcoming the limitation described in (1) above.
- 3. Initial and final pressures in the supply tank can be measured with analog gages, or continuous pressure readings can be made and recorded using pressure transducers. While the transducers add cost, they also add increased resolution and the opportunity to follow the rate of the pressure drop and hence the rate of reaction.
- Enhanced precision can be achieved by measuring the temperature in the supply tank and applying corrections as appropriate.

Parr offers a series of high pressure burettes in complete packages for direct connection to our reactors. The basic ones are listed on the following page.

These burettes can also be equipped with digital pressure transducers, internal thermocouples and data acquisition and reduction support. Please contact our customer support group for information on these possibilities.

Mass Flow

Parr Instrument Company can provide mass flow meters or controllers for quantitative mass flow based analysis. Mass flow controllers are mass flow meters that incorporate an integral control valve, external valve, or feed pump to control the fluid flow. Mass flow controllers are typically used in automated or semi-automated systems. Due to many application and calibration specific requirements, please contact Parr Instrument Company for technical assistance with mass flow meters or controllers.

Parr offers multiple, price driven, electronic interface devices for mass flow meters and controllers. The Parr A2200E Mass Flow Meter/Controller interface system offers a manually operated readout and/or set point module for up to four mass flow meters/controllers. When the A2200E is used with a mass flow controller, a manually operated back pressure regulator is required. The Parr 4871 Process Controller offers remote set point, readout, data logging, totalizing, gas mixing, and process related interfaces with these mass flow devices. Other intermediate interfaces can be provided.

High Pressure Gas Burettes

Parr offers a series of high pressure burettes intended to introduce gas (commonly hydrogen) to a reactor at a constant pressure. The burettes consist of a high pressure reservoir equipped with an inlet valve, a pressure gage and a relief valve. A constant pressure regulator with a check valve, a connecting hose and a support stand are included with each pipette.

The amount of gas consumed in a reaction can be determined by knowing the volume of the high pressure reservoir and observing the pressure drop in the reservoir during a reaction.



Parr high pressure burettes can be furnished in various sizes as shown in the adjoining table, each with a regulator to deliver gas to the reactor over the designated pressure range. The moles of gas shown in the table represent the amount of hydrogen that will be held in the burette at the maximum pressure. The deliverable volume will be a function of the difference in pressure between the pipette and the reactor. The size of the burette selected should be large enough to provide sufficient gas to complete the reaction while still maintaining sufficient pressure in the burette to force gas into the reactor.

Reservoirs with larger volumes are available as are regulators with different delivery ranges. Modifications can be made to these basic systems to add an internal thermocouple to the reservoir and/or a pressure transducer for digital readout and/or recording.

Liquid Charging Systems

Liquid Metering Pumps

Liquid metering pumps are commonly used to introduce liquids into a reactor or vessel at elevated pressures on a continuous basis. A wide variety of pumps are available to meet various pressure, flow, and control requirements. The pumps listed here cover some of the more common pressure and flow requirements associated with Parr reactors and pressure vessels. The pumps described under these catalog numbers include an inlet filter, a reverse-flow check valve and the outlet tubing to the reactor. Special pumps can be furnished to meet requirements outside the range of these pumps. All pumps can be operated from their faceplates and all except the 2312E can also be remotely controlled with a 0-10VDC analog signal, such as from the Model 4871 Process Controller.

Liquid Metering Pumps					
Part No.	Flow Rate, mL/min	Pressure, Max. psi	Wetted Material	Digital Pressure Alarm / Shut-off	
2312E	0.01-10	2000	PEEK	No	
2313E	0.01-10	5000	Stainless	No	
2314E	0.04-40	1500	Stainless	No	
2315E	0.01-10	5000	Stainless	Yes	
2316E	0.04-40	1500	Stainless	Yes	

Gas Burettes Sizing				
Sample Cylinder Volume, mL	Maximum Pressure, psi	Moles of H ₂	Maximum Pressure (bar) For CE (TPED) Applications	Moles of H ₂ For CE (TPED) Applications
150		0.6	100	0.5
300		1.3		1.1
500	1800	2.3		1.8
1000		4.6		3.7
2250		10.4		8.5
150		1.7		1.5
300	5000	3.4	300	3.0
500		5.7		5.0

Liquid Charging Pipettes

To introduce liquids into reactors or vessels at elevated pressures, the most economical way is to use a pressure pipette as a secondary vessel. These are often used for liquid addition to a batch process. Liquid is forced into the reactor from the pipette by applying gas pressure to the pipette greater than the pressure within the vessel. If the passages in the connecting line are large enough, slurries or catalyst suspensions can also be charged into the reactor in this manner.

The pipettes listed below offer a choice of volumes and are rated for pressures to 1800 psi. They include a nitrogen filling connection for attachment to a nitrogen tank. More elaborate pipette systems can be assembled to special order to include additional fittings, such as a pressure gage for the pipette, a pressure relief valve or a large opening ball valve. Special pipettes can also be furnished for higher pressures to 5000 psi.

Liquid Charging Pipettes			
Part No.	Pipette Volume, mL	Pressure Rating, psi	
A2113HC3	50	1800	
A2113HC4	150	1800	
A2113HC	300	1800	
A2113HC2	1000	1800	

A2113HC Liquid Charging Pipette

127

Solids/Slurry Addition Devices

One of the modifications most frequently requested is a port or other means to feed liquids, solids, or slurries into the vessel without removing the head. This can be done in various ways.

Solids Charging Port with Ball Valve

A ball valve with a 3/8" diameter opening can be installed on any one liter or larger vessel and used in conjunction with a high pressure pipette for injecting slurries under pressure. These are opened or closed with a quarter turn of the handle. Larger diameter valves are available for 1 gallon and larger vessels. These ball



A143VB Ball Valve

valves will withstand the full pressure developed in a reactor at moderate temperatures, but their pressure rating falls off rapidly at temperatures above 100 °C.

Solids Charging Ports				
Part No.	Nominal Size	Orifice Diameter, in.		
A143VB	1/4" NPT (F)	0.250		
A132VB	3/8" NPT (F)	0.375		
396VBAD	1/2" NPT (F)	0.406		

Solids Charging Device

A solids charging device in the head of a reactor can serve as a convenient solids charging port at atmospheric pressure. The body of this device is machined with an internal taper to aid in the delivery of the solids into the vessel. It has a convenient screw cap closure with an FKM O-ring seal for use up to 225 °C. Other O-ring materials are available upon request, either for higher temperature operation or material compatibility.



Various views of Solids Charging Device with internal taper

Internal Catalyst Addition Device (ICAD)

Parr has developed a unique device for adding small amounts of solids (or liquids) from a sealed container

held within a reactor. The ICAD is of particular interest to users performing kinetic studies of catalytic reactions. This device consists of a small cylindrical chamber with a cap that is sealed to the body with an O-ring. It attaches to the underside of the vessel head with a 1/8" NPT connection. To discharge the contents of the holder, gas pressure is applied through a valve installed on the top of the head. When the applied pressure is greater than the pressure within the reactor, the cap is forced open and the catalyst or other contents of the holder will be released into the



A550HC Catalyst Addition Device

reactor. This device works best in the taller, 450 mL and 600 mL Mini Reactors, and in the 1 liter and larger Parr Reactors.

Internal Catalyst Addition Devices					
Complete Reactor	Mounting Size, cc	Assembly No.	Thread		
Mini	6	A550HC3	1/8" NPT		
One Liter	8	A550HC	1/8" NPT		
Larger	20	A550HC2	1/8" NPT		



Reactor	Available Fitting Sizes
Mini	1/4" NPT (M)
1 & 2 Liter	3/8" NPT (M)
Gallon & up	1/2", 3/4", & 1" NPT (M)

External Catalyst Addition Device (XCAD)

The addition of a specialized 5 mL or 10 mL chamber with a tools-free closure and two valves to the above solids charging device adds two important capabilities:

- 1. The chamber/valves assembly can be removed from the reactor without tools and filled in a glove box, for example, if an atmosphere-sensitive catalyst is being used.
- 2. Upon reattachment of the chamber/ valves assembly to the reactor, a high pressure gas source can be connected to the top of the device to allow addition of catalyst to the reactor at elevated temperature and pressure.

Similar to the Internal Catalyst Addition Device (ICAD) shown above, the XCAD is ideal for kinetic studies, and is often paired with a gas burette (page 126) for constant pressure hydrogenation catalyst screening/characterization studies.



XCAD

Cooling Coils

A.

Serpentine Cooling Coil 1000 mL



Spiral Cooling Coil 1000 mL

Internal cooling coils are available for all but the smallest Parr reactors. These coils provide an extremely effective means of removing heat from the vessel to control an exothermic reaction or for cooling the reactor at the end of a test. Since heat is transferred through the relatively thin wall of the coil instead of the thick wall of the vessel, cooling rates are generally much faster than heating rates; particularly at temperatures above 80 °C. Water is normally used as the cooling medium although compressed air can be used for modest cooling loads. Cooling coils are offered in three standard configurations:

Single Loop - Single loop coils consist of a vertical run of tubing formed into a "hairpin" shape. These are normally installed on small reactors where there is minimum space available.

Serpentine Coils - Serpentine coils consist of six to eight vertical runs of tubing uniformly spaced around the circumference of the vessel. These coils provide reasonable surface area, minimum interference with stirring patterns, a reasonable amount of baffling, and ease of cleaning and maintenance.

Spiral Coils - Spiral coils consist of multiple loops wound just inside the inside diameter of the vessel. They are normally available only for the 4" and 6" ID vessels although other sizes have been built on special order. They do maximize the cooling area available, but sometimes at the expense of uniform stirring and ease of cleaning. The individual reactor specifications will dictate the style of coil or coils available for each reactor.

Cooling coils are available in the same choice of materials as the reactor bodies themselves. All cooling coils are removable. Plugs are available to close the openings in the head and in most cases these openings can be converted to alternate inlets/outlets if cooling is not required.



Glass Liners 2000 and 1000 mL Sizes Temperature Limit: 565 °C



PTFE Liners 2000 and 1000 mL Sizes Temperature Limit: 225–250 °C

Liners

Removable, open top,

cylindrical liners made either of borosilicate glass or PTFE can be furnished to fit most Parr reactors and general purpose vessels. These liners slide into the cylinder and require no additional fittings, but they may not coordinate with some alternate accessories and stirrers. Although they will not keep corrosive vapors from reaching the surfaces of the cylinder and head, they make it much easier to add and remove liquid reactants, and they give some protection to the cylinder when working with corrosive solutions. It must be noted, however, that adding a PTFE liner will slow the heat transfer rate into and out of the vessel. and it may be necessary to adjust the temperature control method to prevent overheating.

Liners			
Fits ID, in.	Fits Cylinder Size, mL	Glass Liner Part No.	PTFE Liner Part No.
1.3	50	1431HC	1431HCHA
1.3	100	1431HC2	1431HC2HA
1.5	125	2920HC2	2920HC4HA
1.5	200	2920HC3	2920HC3HA
1.5	75	2920HC	2920HC2HA
2-1/2	250	762HC10	N/A
2-1/2	500	762HC2	762HCHHA
2-1/2	300	762HC	762HC4HA
2-1/2	450	762HC2	762HC5HA
2-1/2	600	762HC3	762HC6HA
2	100	762HC7	762HC7HA
2-1/2	160	762HC8	762HC8HA
3-1/4	600	2312HC	2312HC3HA
3-1/4	1200	2312HC2	2312HC4HA
3-3/4	1000	1441HC	1441HCHA
3-3/4	1800	1442HC	1442HCHA
4	1000	398HC	398HCHA
4	2000	399HC	399HCHA
6	1 Gallon	894HC	894HC4HA
6	2 Gallon	894HC2	894HC5HA



Liquid Sampling

Sample Collection Vessel

A sample collection vessel can be added to most reactor systems. Designed to efficiently and safely allow for the withdrawal of liquid or vapor samples at elevated temperatures and pressures, this quick close, O-ring seal vessel has a volume of 5 mL or 10 mL and is designed for operating pressures to 3000 psi (200 bar).

The typical arrangement for this sample vessel includes a cooling sleeve, isolation and vent valves. A drain valve may also be added to the vessel.

The isolation valve is mounted at the head of this vessel and is used to seal the vessel once the sample is transferred. The vent valve is installed in a tee and is used to release any residual pressure in the line between the sample valve and the sample vessel. Samples can be removed either by opening the collection vessel and pouring it out or by use of the drain valve.

Standard material of construction is T316 Stainless Steel but it can be provided in any of the other alloys if required. A high pressure 25 mL or 75 mL sample collection vessel without a cooling sleeve for pressures to 5000 psi is available upon request.

Sample Collection VesselsPart
No.Description4351Sample Collection vessel, 10 mL, with cooling sleeve,
isolation & vent valves for connection to 1/8" NPT valves4352Sample Collection vessel, 10 mL, with cooling sleeve, isolation &
vent valves for connection to 1/4" NPT valves4353Sample Collection vessel, 10 mL, with cooling sleeve, isolation &
vent valves for connection to 3/8" NPT valves4353Sample Collection vessel, 10 mL, with cooling sleeve, isolation &
vent valves for connection to 3/8" NPT valves

-D Optional Drain Valve



4878 Automated Liquid Sampler

safe and reliable method for collecting multiple liquid samples from heated and pressurized reactors has been sought for many vears. Parr Instrument Company is pleased to introduce the Parr 4878 Automated Liquid Sampler which can extract up to six liquid samples at a user-defined interval without the need for the continued presence of an operator. Controlled through a user-friendly touch screen, a series of precision switching valves allows collection and deposition of samples with consistent volumes into individual vials. The 4878 can operate up to the maximum working pressure of all standard Parr reactors and pressure vessels. A multi-step loop sequence ensures clearing of the reactor dip tube between samples to yield samples representative of the bulk reactor fluid.

Features include:

- Touch screen controller with easy-to-navigate graphical displays
- Compact footprint 14-in wide x 16-in deep (35cm x 40cm)
- Stand alone design compatible with any new or existing Parr reactors and pressure vessels
- User-definable parameters including time between samples, number of samples, and number of loop sequences
- Various sample loop volumes available to accommodate a wide variety of sampling scenarios
- Capability to connect with mobile devices on both iOS and Android platforms
- Possible customizations include: design for use in hazardous locations, increased number of samples, special alloy components. Contact Parr Technical Service for additional options to fit your requirements.



Filtered Dip Tubes

To assure you collect only liquid from a vessel, the solids must be removed. Parr recommends press-on filters, or in some instances screw-on filters, at the end of the dip tube to remove solids from the liquid sample. These sintered metal filters are available in a variety of materials, with a variety of nominal pore sizes, often from 0.5 microns to 40 microns. Most of these filters have a relatively small filtration area, so the user should consider how long it will take for a sample of desired volume, viscosity, and solids content to be collected in a sampling device. If clogging of a press-on filter is a possibility, Parr recommends purchasing a second easy-to-install complete filtered dip tube as the filters are factory-installed and not easily replaced in the field. Screw-on filters are field-replaceable and have a larger surface area, but due to their larger size usually only supplied for use with 1 Liter and larger vessels.

Catalyst Testing System

2280 Burette, 4566 Reactor with 4848 Controller, and 4878 Automated Liquid Sampler

he recent introduction of Parr's 4878 Automated Liquid Sampler has allowed the construction of a complete system for activity/kinetic testing of catalysts in a high pressure, high temperature environment. In the system shown, a Parr 2280 Gas Burette with 50 mL reservoir delivers H₂ at a user-settable, constant pressure to a 4566 300 mL stirred reactor. Operation at up to 180 bar at up to 350 °C is possible. A Parr XCAD (eXternal Catalyst Addition Device) allows introduction of solid/powder catalyst to the liquid contents of the reactor at reaction pressure and **temperature**, setting t=0 for kinetic experiments. An included gas entrainment impeller provides for high speed three phase mixing. A Parr 4878 Automated

Liquid Sampler allows unattended collection of filtered liquid samples from the reactor at a user-defined interval for subsequent analysis.



Bottom Drain Valves



Bottom drain valves can be added to most Parr reactors. These valves are particularly useful for those working with polymers or other material that must be discharged from the reactor while they are still hot and before they can solidify. These valves are also quite useful for the 1 gallon and larger vessels which are too large to conveniently lift from the heater for product recovery. Bottom valves are rarely installed on the micro and mini reactors with their small volumes and light vessel weights.

The standard bottom drain valve has a rising stem, that is flush with the inside cylinder bottom so that there is no dead space between the bottom of the vessel and the shut off point of the valve. In the fully open position the stem is retracted completely to open a clear passage for draining the vessel. When the valve is reclosed, any material in this

A465VB Bottom Drain Valve

Bottom Drain Valves					
Part No.	Opening Dia., in.	Outlet Connection	Max. Press., psi	Max. Temp, °C	Seal
A485VB	0.20	1/4" NPT (F)	3000	225	PTFE
A485VB2	0.20	1/4" NPT (F)	3000	350	Silver
A485VB3	0.20	1/4" NPT (F)	3000	350	Silver
A465VB	0.34	3/8" NPT (F)	2000	350	Grafoil
A465VB2*	0.34	3/8" NPT (F)	2000	350	Grafoil
A465VB3	0.34	3/8" NPT (F)	2000	350	Silver
A177VB2	0.31	3/8" NPT (F)	5000	500	Silver
A296VB2	0.69	1" NPT (F)	1900	350	Silver

* Set up for a Band Heater.

passage will be pushed back into the reactor by the rising stem. Valves with 3/8" diameter clear passage are recommended for vessels with volumes from 1000 mL to 2

gallons. A 1/4"

valve is available

for 600 mL and smaller vessels. High pressure and larger diameter valves are available where required.

These valves will withstand the full operating pressures and temperatures of the vessels in which they are installed in the closed position. They are available in all of the current Parr materials of construction. Users can also specify that a reactor ordered with a bottom valve shall have a tapered bottom so that it will drain easily through the valve opening.

Not all Parr reactors will accept a bottom drain valve. Since the valve extends approximately 8 inches below the bottom of the vessel, the entire vessel must be raised by this amount to accommodate the valve. This makes some models too tall for convenient bench top operation. The specification tables for each model will identify those reactors in which a bottom drain can be readily installed, and those which will not accept a bottom drain, or those which will require custom modification of the heater and support stand to accommodate a bottom valve.

Needle Valves and Ball Valves

Needle valves and ball valves can also be installed as bottom outlet valves. Needle valves are generally used on the smaller reactors. While ball valves can be used for large discharge passages, they are generally limited in their operating temperature/ pressure capabilities and they leave a fairly large dead space between the bottom of the vessel and the seat of the valve.

131

Valves and Fittings

Parr stocks and can install a wide variety of valves and fittings for use with reactors and pressure vessels. These include:

- Needle Valves with NPT or tube connection.
- Regulating Valves with NPT or tube connection.
- Ball Valves with NPT or tube connection.
- High Pressure Valves
- Severe Service Valves
- Remote Operating Valves
- Tube Connectors
- Pipe Connectors
- Plugs
- Union Coupling Adapters Please contact our customer service department for details.

Manual Control Valves for Compressed Gas Tanks

ank valves with couplings to fit standard compressed gas cylinders are available in stainless steel for corrosive cases and in nickel plated brass for noncorrosive gases. The brass valves have a 2-1/2" diameter pressure gage which shows the tank pressure. Both styles have a 1/4" NPT female outlet which will accept any pressure hose or gas tube assembly. These valves do not regulate the delivery pressure of the gas. Pressure regulators are available on special order.

T303 Stainless Steel Valves-No Gage			
Fits CGA Tank Valve No.	Outlet No.	Typical Usage	
A120VBPN	510	Propane, butane, ethylene oxide	
A120VBPP	660	Chlorine, sulfur dioxide, nitric oxide	
Nickel-Plated Brass Valves with Cylinder Pressure Gage			
Fits CGA Tank Valve No. Outlet No. Typical Usage			
		л с	
A120VBPQ	320	Carbon dioxide, methyl bromide	
A120VBPR	350	Hydrogen, carbon monoxide, ethylene	
A120VBPS	540	Oxygen	
A120VBPT	580	Nitrogen, argon, helium	
A120VBPU	590	Air	

Note: Can be furnished with DIN/BSP connections on special order

Safety Check Valves



364VB Check Valve

Whenever gases or liquids are introduced into a vessel under

pressure, the supply pressure must be greater than the pressure in the vessel to prevent reverse flow back into the supply system. Protection against reverse flow can be obtained by installing a check valve in the supply line. With a check valve in the line, the valve will snap shut if the supply pressure is lower than the pressure in the vessel. This protection is particularly important on stirred reactors

where gas enters through a dip tube. With liquids in the vessel, any back pressure will force liquid back into the gas tank or into the gas supply system.

Parr stocks the poppet check valves for incorporation into the user's supply lines. These valves are typically constructed of T316SS with FKM O-rings with a 10 psi cracking pressure. Alternate O-ring materials and cracking pressures are available upon request.

Pressure Hoses

Three different pressure hose assembles are

available for high pressure gas connections to both stirred and non-stirred vessels. The standard hose is a 6 foot length with a male "A" socket connector on one side and a 1/8" NPT (M) nipple with a 1/4" NPT (M) bushing on the other end. The "A" socket side of the hose attaches to couplings installed on the inlet valve of all stirred reactors and non-stirred vessels, as well as to a side port of the gage block assembly. The choice of either 1/8" NPT or 1/4" NPT on the opposite end of the hose allows for attachment to most gas tank valves, pressure regulators or other gas supply sources.

The A495HC Hose Assembly is made of nylon. It is rated for 2500 psi and is very flexible and easy

to use with dry, non-corrosive gases (nitrogen, hydrogen and oxygen). Care must be taken to ensure that the nylon hose does not come in direct contact with any hot surfaces on the vessel or heater. One of these hoses is included with each complete Parr Series 4500, 5100, and 5500 Stirred Pressure Reaction Apparatus.

The A490HC Hose Assembly is a braided, stainless steel hose with a PTFE lining, rated for 2500 psi. It is reasonably flexible and recommended for use with corrosive gases and liquids, and for applications



A495HC Pressure Hose

requiring additional abrasion resistance, but it is not intended for high temperature liquids or gases.

The A506HC Tube Assembly is a 6-foot length of 1/8" OD stainless steel tubing, rated for 7500 psi. This small diameter tubing is "bendable", but it is not as flexible as the other hoses. It is recommended for corrosive gases, high temperature transfers and other high pressure applications. Special versions of this assembly can be made of other corrosion resistant materials. Larger tubing can be used, but it is rigid rather than flexible.

Special hoses with different lengths or end fittings can be assembled for special orders.

Pressure	Hose
A495HC	Pressure Hose Assembly, 6-ft, reinforced Nylon
A495HC5	Pressure Hose Assembly, 6-ft, reinforced Nylon, with check valve
A495HC7	Pressure Hose Assembly, 10-ft, reinforced Nylon
A495HC8	Pressure Hose Assembly, 10-ft, reinforced Nylon, with check valve
A490HC	Pressure Hose Assembly, 6-ft, PTFE-lined, braided stainless steel
A490HC5	Pressure Hose Assembly, 6-ft, PTFE-lined, braided stainless steel, with check valve
A506HC	Pressure Tube Assembly, 6-ft, 1/8-in OD, T316SS
A506HC2	Pressure Tube Assembly, 6-ft, 1/8-in OD, T316SS, with check valve

Alternate lengths available upon request.

Thermocouples

Parr offers a variety of thermocouples for use in our reactors and pressure vessels. The standard thermocouple is a Type J (iron-constantan) which is compatible with the operating temperature range of these vessels.

The thermocouples are furnished with a sealed 1/8" OD stainless steel sheath and include a standard plug connection at the end of the probe. Our standard thermocouples are manufactured in accordance to ASTM E230.

Alternate thermocouple materials including Alloys C276 and 600 are readily available. Platinum resistance elements (3-wire RTD) are available as special orders as well as multiple point thermocouples.

Most commonly, in small volume vessels the thermocouple probe is installed directly into the vessel with a compression fitting and in larger vessels the probe sits inside a thermowell. The thermowell arrangement offers protection to the thermocouple from physical damage. We also furnish thermowells in vessels manufactured in

materials other than stainless steel so the thermowell will be the alternate alloy and the thermocouple probe can be stainless steel.

Additionally, dual element thermocouples with two separate thermocouples in a single sheath are furnished in smaller volume vessels for use with accessory temperature meters. We also offer spring loaded thermocouples which are designed to be installed through the heater wall to the outside wall of the pressure vessel.

An extension wire is furnished to connect the thermocouple to the control device. The standard length is 5 feet but longer lengths are available if the control is to be mounted away from the reactor.

Type J Thermocouples with 1/8" Diameter			
Part Number	Stem Length, in.	Sheath Material	
A472E	7.5	T316 Stainless Steel	
A472E2	9.5	T316 Stainless Steel	
A472E3	11.5	T316 Stainless Steel	
A472E6	15.5	T316 Stainless Steel	
A472E5	21.5	T316 Stainless Steel	
A472E4	5.5	T316 Stainless Steel	
A472E8	2.5	T316 Stainless Steel	

Most of the above listed thermocouples are also available as Type K (Chromel-Alumel) or Type T (Copper-Constantan).



Equipment for Use in Potentially Ignitable Atmospheres

Parr reactors are typically equipped with totally enclosed variable speed motors, electric heaters, and controllers intended for use in non-hazardous environments. These standard units can be used in most laboratories without undue hazard, but there will be situations where the installed equipment must be considered for use in ignitable atmospheres. Parr offers various optional stirrer drives and heating solutions to meet these strict requirements.

USA and Canadian Codes (HAZLOC – Hazardous Locations)

Designing electrical equipment to be operated in hazardous locations is a complex subject, which is governed by extensive national electrical codes and supplemented by local regulations. These codes require all electrical equipment that is installed in a governed location must be approved for use with the specific gas, vapor, or dust that can be present in the defined location. USA and Canadian electrical codes classify hazardous locations according to the nature and concentration of specific hazardous or flammable materials. These are divided into three classes:

- Class I Flammable liquids, gases or vapors.
- Class II Combustible or electrically conductive dusts.
- Class III Easily ignitable fibers/flyings.

There are two divisions within each of these classes.

- **Division 1 –** Where the flammable material exists in the atmosphere under normal operating conditions.
- Division 2 Where the hazardous material is confined within a closed system from which it may be released only under abnormal conditions, such as a gas leak in the system.

Class I locations are further subdivided into four groups, A, B, C and D which identify specific explosive gases and vapors. Explosive dusts and fibers in Class II are subdivided into Groups E, F and G. Most hazardous applications for Parr apparatus will occur in atmospheres identified by Class I, Group B for hydrogen and Groups C and D for most other combustible gases and vapors. Class II, Group F covers coal dust. Most other combustible dusts, such as flour and grain, are in Group G. Minimum ignition temperatures and energy levels are established for specific materials in each group.

The European Community has corresponding classifications for "Explosive Atmospheres" referred to as ATEX (ATmospheriques EXplosives). Parr will work with all users to provide equipment compatible with their own local codes.

The components in Parr reactor systems that may be considered hazardous and the steps that can be taken to reduce or eliminate the hazards they represent are described below.

Motors

Because of sparking from brush contacts, permanent magnet DC electric motors clearly represent the principal ignition source introduced by a stirred reactor. Electric motors approved for Class I (Divisions 1 & 2), Groups C



Model 4524 Reactor, 2000 mL, Fixed Head Style with Aluminum Block Heater

and D, and Class II (Divisions 1 & 2), Groups F and G atmospheres are readily available in most sizes and voltages. These totally enclosed motors are suitable for many hazardous applications, and they are sometimes used with hydrogen, though they are not approved for Group B atmospheres. Currently, there are no Division 1 motors available for Group A or B atmospheres. A special air purging system can be used to reduce the classification inside the motor. The motor is pressurized by building up a positive pressure of air, or inert gas, within the motor to prevent explosive gases or vapors from entering the motor housing. Division 2 requirements are not as stringent. Other motor options are available. Please contact Parr for additional information for Class I, Division 2 requirements.

Parr can provide ATEX rated, IEC framed, AC motors when required. These constant torque motors have lower speed and torque characteristics than the DC motors Parr has used in the last several years.

An alternate method of dealing with the explosion hazard is to use an air driven motor. These are powered by compressed air and offer a convenient and satisfactory drive system for use in flammable atmospheres, including hydrogen. They are available in sizes suitable for most Parr reactors.

Heaters

The advisable way to heat a Parr reactor in a potentially ignitable atmosphere is to use a hot oil jacket and ensure that the highest temperature of the heat transfer media is below the minimum ignition temperature for the classified area in which it will be installed.

Purging Parr's aluminum block heaters with air, as well as limiting surfaces temperatures below the auto-ignition temperature, can make them suitable for use in classified areas, but it is likely that the air consumption in such heaters will be quite large. Please see our <u>Heaters section</u> on page 116 for additional information.

Wiring

Parr will provide ordinary location wiring for motors and heaters, which can be used for initial testing, setup, etc. Due to national and local requirements for installation of such equipment, it is the user's responsibility to install mains and motor wiring per code requirements.

Parr will provide wiring for intrinsic safety related sensors, transducers, etc. For application specific information on sensors for use in hazardous locations, please contact Parr Instrument Company.

Controllers

The most commonly used method for dealing with the ignition hazard introduced by a temperature or process controller is simply to locate the controller outside of the hazardous atmosphere. Another choice is to install the controller in an explosion-proof cabinet or a cabinet that can be purged using the purging system described above.



Windows









Custom window

Oblong window

Round, screw-in window

Round, integral window

Windows can be installed in Parr stirred reactors and pressure vessels for visual observations, light transmission and other purposes. They usually are installed in pairs so that light can be introduced through one window while the other is used for viewing. Our standard material for these windows is fused silica. Sapphire is also available for small diameter windows. Alternative window materials (with coatings, if requested) are available for specific transmission requirements. Windows can be mounted in several different ways.

Screw-in Circular Windows

The simplest window is a screw-in type with a ¹/₂-inch diameter viewing area. The windows in these assemblies are sealed in a fitting which screws into the vessel using a standard ½ inch NPT male pipe thread. Obviously, the vessel wall must be thick enough to provide full engagement for this thread. O-ring seals restrict the maximum operating temperature to 225 °C or less, depending upon the O-ring material. Alternatively, a Grafoil® gasket seal option is available to increase temperature rating up to 350 °C. Pressure ratings range from 1900 to 5000 psi, depending upon the window material and its thickness. Although these windows are rather small for straight optical viewing, they work well for small video systems and for laser and other analytical beams. A limitation of this design is that there is a dead space approximately 1.25 inches long between the inner face of the window and the inside wall of the vessel.

Integral Windows

Parr has developed designs for installing windows in the wall of the vessel so that the inside face of the window is very close to the inside wall of the vessel. This eliminates the large dead space associated with screw-in windows. These windows are offered in the two styles described below. The maximum size of the window will depend on the size of the cylinder in which it will be installed.

Custom window

Round Windows with a 1/2-inch diameter viewing area are the standard. Round windows are available in a variety of materials including sapphire for very high pressures (up to 5000 psi). Both O-ring and Grafoil[®] seals are available in this design. This type of window is generally used for visual, photographic or optical sensor observations.



Flange-mounted larger window

Oblong Windows with a viewing area 3.5-inch long and 0.62-inch wide are the standard size and can be installed on 4590/4790 vessels of 100 mL volume or 4560/4760 vessels of 450 mL volume and larger. Only O-ring seals are available in this design. Both Fused Silica and sapphire windows are available – maximum rating for these windows are limited to 1900 psi at 225 °C with FKM or FFKM seals. These windows are commonly used for visual observations of both the vapor and liquid phases or for observing the liquid level in the vessel. Multiple windows can be stacked or staggered on larger vessels.

The windows described above as standard are maintained in our inventory for readily available replacements. Custom windows in both the round and oblong styles can be furnished in larger sizes upon request. All reactors and pressure vessels equipped with windows require custom designed heaters and supports. Flexible heating mantles, integral cartridge heaters, and attached circulating jackets (600 mL mini and larger) are the most commonly used heaters for window vessels.

Externally Welded/Flange Mounted Larger Windows

Large round windows that are either externally welded or flange mounted to the vessel can be installed. These will provide a viewing area of 1-3/4-inch to 7-1/4-inch depending on the size of the window ordered and the size of the vessel. Externally welded windows greatly reduce the maximum working pressure of the vessel to 600 psi or less, although custom configurations rated higher are available.

Certification

The windows described above are considered "proprietary fittings" and as such are excluded from the scope of ASME Section VIII Division 1.

Most of the standard windows furnished by Parr can be installed on vessels bearing the CE mark. The windows furnished on these vessels are subjected to various tests on a per lot basis in order to satisfy requirements of 2014/68/EU (Pressure Equipment Directive). Please contact Parr Instrument Company for further information.

Insulated Electrical Glands

A variety of insulated electrical leads can be installed in most Parr reactors or pressure vessels for electrical connections required inside the vessel. Four commonly used examples are discussed below.

Transducer Glands

Transducer glands are available to allow up to 16 individual wires to be sealed and insulated through a single gland. These types of glands are typically used with low voltages and current.

Electrode Glands

Applications requiring a single electrical conductor having high current or very high voltage capacities can be handled with an electrode gland. These glands provide a seal for a single conductor or electrode with the ends of the conductor threaded so that internal and external lead wires can easily be attached.

Power Leads

Power leads can be provided with either single or multiple flexible wires. Current ratings range from 5 to 20 amperes at up to 600 volts. PTFE or ceramic insulating glands are commonly available. Ceramic glands can be used to the full temperature rating of most Parr vessels. Pressure ratings will vary depending upon the design of the gland, its size, and the type of insulation used.

Miscellaneous Sensors

Parr has installed a number of different sensors in its various reactors and pressure vessels, including single point and continuous liquid level sensors, pH electrodes, and dissolved oxygen electrodes. Each of these installations must be carefully developed in consultation with the user, the electrode or probe supplier, and Parr's Engineering Department. Glass electrodes with O-ring seals will obviously carry rather strict temperature and pressure restrictions. There are also space restrictions which generally dictate that accessories of this type can only be installed in 1000 mL or larger vessels.

Spare Parts Kits

Each stirred reactor is furnished with a set of spare parts and fittings including a 6-foot gas supply hose, head gaskets, rupture discs, and a set of replacement parts for the stirrer drive.

A reserve supply kit of spare parts can be ordered from Parr Technical Service to provide sufficient parts and tools to handle most normal replacements and emergency repairs during the first year of heavy usage. These kits include replacement gaskets, O-rings, rupture discs, drive belts, and seals. These kits are a convenient package of the small perishable items required for normal maintenance of the reactor.

When ordering any kit for an existing reactor please provide the serial number, specify the preferred gasket/seal material, the burst pressure of the rupture disc, material of the reactor, and the length of the drive belt.



Temperature Limits

External Valves and Fittings

There are a number of factors that determine the maximum temperature rating of a pressure vessel. For most applications it is the gasket material. Vessels with O-ring seals are limited to 225 °C and those with FKM or FFKM are limited to 300 °C. Parr's design for contained PTFE gaskets extends the operating temperature range to 350 °C. Flexible Graphite (FG) material essentially removes the gasket as the limiting factor. Maximum temperature limits for the metals used in these vessels are established by ASME code and other standards. Most metals have maximum temperature limits between 400 and 800 °C. The allowable strength for these metals falls off rapidly as they reach maximum operating temperature. Finally, the difficulties encountered with screw threads and other closure components operating at high temperatures establish a practical temperature limit for externally heated vessels. We have found 600 °C to be a reasonable limit.

Internally Heated Vessels

Exposed Heaters. Another approach that has proven useful in extending the maximum temperature limit is to place the heater inside the pressure vessel. The heater is surrounded by a layer of insulation. This creates a hot zone in the center of the vessel and prevents the walls from exceeding their allowable limit. As this system is very energy efficient, internal heaters can be less powerful than external heaters. Internally heated vessels are equipped with insulated electrical feed-throughs to power the heater. Multiple thermocouples are used to control or monitor the temperatures in the hot zone and on the vessel inner wall.

The reactions or studies carried out in internally heated vessels must be limited to those which will not destroy the exposed internal heaters and insulation. These are normally gas-solid reactions or controlled atmosphere heat treatment studies. The heating elements are normally ceramic. Some users have developed induction style heaters and insulators and have extended their investigations to above 2500 °C.

Although internal heaters can be installed in almost any non-stirred Parr pressure vessel, vessels with larger I.D.s are required in most applications.

Protected Heaters. Internally heated vessels have also been manufactured with cartridge type heating elements inserted in specially designed "heater wells". These wells protect the heater from the reactants and expand the applications that can be studied. Cartridge type heaters have a maximum temperature of 760 °C.

Materials of Construction

In the standard configuration, the valves, gage, magnetic drive, and other external parts on Parr reactors are furnished in stainless steel, even when a different material is specified for the head, cylinder, and internal wetted parts. The external stainless components are typically only exposed to the vapor of the reactants and are at much lower temperature than the cylinder and internal fittings. These conditions allow stainless steel external fittings to perform satisfactorily in most cases. If external parts made of a material other than stainless steel are required for safety or other reasons, Parr can accommodate this in most cases. Any request for external parts made of a specific material must be stated clearly when ordering.

Valves

Most reactor valves are also available in Alloy 400 at a reasonable cost premium. Valves made of Alloy C-276 are also available, but generally only on special designs and at a considerable cost premium. Soft materials such as titanium and zirconium generally make poor performing valves.

Gages

Pressure gages are available in stainless steel and Alloy 400. Other materials of construction are not available. The standard method for protecting the gage in a corrosive environment is to install a diaphragm gage protector. These have a flexible diaphragm which isolates the gage from the reactants and a sealed hydraulic connection for pressure transfer to the gage. These assemblies are too large to install on all but the largest Parr reactors. The diaphragm assemblies are large and may not be practical on small reactors.

As an alternative, Parr has designed an oil filled piston isolator gage protector to isolate the gage (and transducer, if required) on reactors and pressure vessels where space is limited. These isolators can be furnished in any of the current Parr materials of construction. A piston style isolator can limit the precision the pressure measuring device.

Pressure Transducers

Pressure transducers are commonly available in stainless steel and Alloy C-276 though other materials may be available. Parr provides a mounting adapter with a water cooling jacket on pressure transducers to protect them from excessive temperatures. These can be augmented with isolators similar to gage protectors when corrosion resistance is required. When a gage and a pressure transducer are installed, a single isolator can protect both.

Magnetic Drives

Magnetic drives can be furnished in all of the current Parr materials of construction except nickel, which is magnetic.

Rupture Discs

The standard materials of construction for rupture discs is Alloy 600, gold-faced Alloy 600, and Alloy C-276. Tantalum and other premium materials are available on special order.

Please see the Safety Rupture Disc Assemblies manual 231M for available rupture discs.

Non-Stirred Reactor Systems Chapter 8



Inside This Chapter You Will Find:

4700 22 & 45 mL

4740 25 & 75 mL High Pressure

4750 125 & 200 mL

4790 25-100 mL

4790 HPHT 25-100 mL High Pressure/High Temperature

4760 100-600 mL

4760 HT 300-600 mL High Temperature

4600 1000 & 2000 mL

4625 600 & 1200 mL High Pressure/High Temperature

4650 250-1000 mL High Pressure/High Temperature

4660 1 & 2 Gallon

4680 1000 & 1800 mL High Pressure/High Temperature

4670 1 & 1.5 Gallon High Pressure/High Temperature

4675 2.6 & 5 Gallon

Gage Block Assemblies

Coned Pressure Fittings

Non-Stirred Pressure Vessels

Parr offers a broad selection of non-stirred pressure vessels in convenient styles, sizes and pressure ranges for many laboratory uses. **Recent additions** to this line have included new sizes, new temperature and pressure combinations and new self-sealing closures, all designed to provide convenient vessels for laboratory procedures that must be performed under pressure at elevated temperatures.



Pressure Vessel Selection Procedure

The steps to be followed in selecting a non-stirred pressure vessel are similar to those used for a stirred reactor. Selection begins with the Guide to Non-Stirred Vessels, from which the user can identify the Model Number for vessels in a range of sizes with temperature and pressure ratings suitable for the intended application. Detailed specifications and catalog numbers can then be found on the referenced page. After the specifications have been established, the options and accessories required to customize a complete pressure vessel system can be selected.

Series Number	Vessel Size	Max. Pressure psi (bar)	Maximum Temperature °C	See Page Number
<u>4700</u>	22 & 45 mL	1700 (115)	300	<u>142</u>
4740	25 & 75 mL	8500 (586)	538	<u>144</u>
<u>4750</u>	125 & 200 mL	3000 (207)*	350	<u>146</u>
<u>4790</u>	25-100 mL	3000 (207)*	350	<u>148</u>
<u>4790HP</u>	25-100 mL	5000 (345)	350	<u>148</u>
<u>4790HP/HT</u>	25-100 mL	5000 (345)	500	<u>148</u>
<u>4760</u>	100-600 mL	3000 (207)*	350	<u>152</u>
<u>4760HT</u>	300-600 mL	2000 (138)	500	<u>152</u>
<u>4600</u>	1000-2000 mL	1900 (131)	350	<u>156</u>
<u>4600HP</u>	970-1900 mL	2900 (200)	350	<u>156</u>
<u>4625</u>	600-1200 mL	5000 (345)	350	<u>160</u>
<u>4650</u>	250-1000 mL	6000 (413)	600	<u>162</u>
<u>4660</u>	1 & 2 gallon	1900 (131)	350	<u>164</u>
<u>4680</u>	1000 & 1800 mL	6000 (413)	600	<u>166</u>
<u>4670</u>	1 & 1.5 gallon	3000 (207)*	500	<u>168</u>
<u>4675</u>	2.6 & 5 gallon	1900 (131)	350	<u>170</u>

Custom options available. Contact Parr's Customer Service Department.



Head Design

Many of the non-stirred vessels offer a basic head with an NPT plug. These should be considered as only a starting point for custom built vessels as they *do not* include a safety relief device. Safety codes, good practice, and common sense dictate that a safety relief device should be installed on all vessels. The NPT plug is provided for the user to install their own safety relief device.

Users who take delivery of vessels without safety relief devices installed must take responsibility for installing adequate protective devices before the vessel is placed in service.

Model 4600-2L, with 1/4-in plug

8



Parr Non-Stirred Vessels with Volumes from 22 to 2000 mL.



Available Accessories Include Heaters, Controllers, Valves, Gage Block Assemblies, and Gas Filling Hoses

Index to Non-Stirred Vessel Design Features, Options and Accessories

Primary Accessories	Reference Page
<u>Heaters</u>	<u>116</u>
<u>Controllers</u>	<u>99</u>
Gage Block Assemblies	<u>173</u>
Valves and Fittings	<u>132</u>

Options	Reference Page
Gaskets and Seals	<u>18</u>
Materials of Construction	<u>10</u>
Pressure Gages	<u>125</u>
Rupture Disc	<u>124</u>
Certification	<u>8</u>

Accessories*	Reference Page
<u>Liners</u>	<u>129</u>
Pressure Hoses	<u>133</u>
<u>Check Valves</u>	<u>132</u>
Liquid Pipettes	<u>127</u>
Gas Filling Systems	<u>126</u>

* May require modifications of vessel to install

141



Type: General Purpose

Mounting Style: **Moveable**

Vessel Sizes, mL: 22 & 45

Standard Pressure MAWP Rating, psi (bar): **1700 (115)**


Maximum Operating Temp., °C: **300**

Series 4700 Pressure Vessels, 22 & 45 mL

These are the smallest of the Parr General Purpose Pressure Vessels. They have a 1-inch inside diameter and are offered in two different lengths with

volumes of 22 and 45 mL. These vessels have a flat PTFE or flexible graphite gasket that is sealed with a screw cap. Special wrenches and bench sockets are offered for tightening these closures. The screw caps and bodies are made of alloy steel for use at temperatures to 300 °C. The maximum working pressure is 1700 psi (115 bar). The standard material of construction for these small, inexpensive vessels is Type 316 Stainless Steel, but they are also available in all of the current materials of construction Parr provides.

There is room for a single opening on the head of these vessels. This is generally an "A" socket which will accept a 4310 Gage Block Assembly.



4700-45mL-A Cross Section

These heads can also be furnished with 1/8" NPT, 1/4" NPT, or a rupture disc assembly. These vessels are normally

heated in ovens, baths, or similar general purpose heating devices. Special heaters for these vessels are not available from Parr.



Model 4700-22mL-A Vessel shown with needle valve.

Shac	led bar indicates specifications that change wi	thin series.	
Model Number	4700		
Approximate Volume, mL	22	45	
Maximum Pressure (MAWP)	1700 psi	(115 bar)	
Maximum Temperature	30	0 °C	
Vessel Details			
Mounting Style	Movea	ble Only	
Closure	Screw C	ap - Steel	
Head Seal	Flat Gasket		
Head Opening Options	One head opening (1/8" NPT, 1/4" NPT, "A" socket, or rupture disc assembly port) Custom port option (Contact our technical sales department to discuss options)		
Maximum Head Openings	1		
Temperature Measurement	N/A		
Stand and Heater Type	N/A		
Weight and Dimensions	•		
Cylinder I.D. x Depth, inches	1.0 x 1.6	1.0 x 3.8	
Vessel Assembly Weight, pounds*	0.8	1.10	
Cylinder Weight, pounds	0.2	0.5	
	ht is based on a moveable head vessel with a ht will vary based on head configuration choice		
Other options availabl	e. See Ordering Guide, visit <u>www.parrinst.com</u> ,	, or call for more information.	

Series 4700 Ordering Guide


An example order number for a vessel this series is: 4700-45mL-T-SS-A-21AC4-A22AC3

A composite identification number to be used when ordering a 4700 Series Pressure Vessel can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Volume		
-22mL	22 mL, 1" ID	
-45mL	45 mL, 1" ID	
B Gasket / N	laximum Temperature	
-T	PTFE Gasket, 300 °C	
-FG	Flexible Graphite Flat Gasket, 300 °C	
C Materials	of Construction	
-SS	T316 Stainless Steel	
-C20	Alloy 20	
-M0	Alloy 400	
-IN600	Alloy 600	
-IN625	Alloy 625	
-HB	Alloy B-2 / B-3	
-HC	Alloy C-276	
-NI	Nickel 200	
-TI2	Titanium Grade 2	
-TI4	Titanium Grade 4	
-ZR702	Zirconium Grade 702	
-ZR705	Zirconium Grade 705	
See <u>page 10</u> or <u>24</u> for complete list of available alloys.		
D Head Cont	figuration	
-RD	Rupture disc port with rupture disc assembly & disc	

-RD	Rupture disc port with rupture disc assembly & disc
-1/8	One head opening, 1/8" NPT port, plugged
-1/4	One head opening, 1/4" NPT port, plugged
-A	One head opening, "A" socket
-NS	Non-standard, custom head opening
E Optional I	Fittings & Custom Options (Check all needed)
-V	Needle Valve
-4310A	Gage Block Assembly (for head with "A" socket opening)
F Required	Accessories
Part No.	Description
-21AC4	Wrench
-A22AC3	Bench Socket
G Certifications	
-ASME	ASME Documentation
-CE	CE Documentation
-P	Parr Certification

Please note that all options and combinations are not compatible with all models.



800 1000

1600

1800

500



Type: **High Pressure / High Temperature**

Vessel Mounting: Moveable

Vessel Sizes, mL: **25 & 75**


Maximum Pressure MAWP Rating, psi (bar): **8500 (586)**

Maximum Operating Temp., °C: **538**

Series 4740 HP/HT Pressure Vessels, 25 & 75 mL

These high pressure/high temperature pressure vessels have a 1-inch. They have a 1-inch inside diameter and offer volumes of 25 or 75 mL. These vessels are closed with an alloy steel screw cap which includes six compression bolts to develop the sealing force on a flat, flexible graphite gasket. Interchangeable PTFE gaskets are also available for users who need the high pressure capabilities provided by these vessels, but who do not need to operate above 350 °C. Maximum pressure drops quickly at temperatures above 350 °C to a rating of 1850 psi (125 bar) at 538 °C. These vessels are available in all of the standard materials currently offered by Parr, although some of these materials will limit the maximum pressures and temperatures available.

The head of these vessels can only accommodate one opening. It is designed to accept the 4310A Gage Block Assembly which can be furnished with or without a thermocouple. Parr offers the 4921 Bench Top Ceramic Fiber heater for use with the 75 mL, 4740 Vessel only. The smaller 25 mL Vessels are generally heated in an oven, bath, furnace or similar general purpose heater (not available from Parr)



4740-75 mL Cross Section

Shaded	bar indicates specifications that change	ge within series.
Approximate Volume, mL	25	75
Maximum Pressure		
Head & Cylinder of T316SS With Standard Screw Cap	4100 psi	i (586 bar) @ 350 °C i (282 bar) @ 500 °C i (127 bar) @ 538 °C
Head & Cylinder of T316SS With Alloy C-276 Screw Cap	8000 psi	i (586 bar) @ 350 °C i (551 bar) @ 500 °C i (544 bar) @ 538 °C
Maximum Temperature		
With PTFE Flat Gasket		350 °C
With Flexible Graphite Flat Gasket		500-538 °C
Vessel Details		
Mounting Style	Moveable	
Closure	Screw Cap (6 Compression Bolts)	
Head Seal	Flat Gasket	
Head Opening Options	One head opening ("A" socket, or rupture disc assembly port) Custom port option (Contact our technical sales department to discuss options)	
Maximum Head Opening	1	
Temperature Measurement (Optional)	Thermocouple (in cylinder bottom or in optional gage block assembly)	
Stand and Heater (optional)		
Stand type		Bench Top Stand with Heater
Heater Style	N/A	Ceramic Fiber
Heater Power, Watts	IV/A	700
Maximum Load, amps, 115V / 230V		6 / 3
Weight and Dimensions		
Cylinder I.D. x Depth, inches	1.0 x 1.7	1.0 x 5.6
Vessel Assembly Weight, pounds*	4	6
Cylinder Weight, pounds	1.4	3.2
*Vessel weight Vessel weight	is based on a moveable head vessel wi will vary based on head configuration c	ith a plugged "A" socket. hoice and fitting options.
Other options available.	See Ordering Guide, visit <u>www.parrinst</u>	.com, or call for more information.

G Stand with I

Series 4740 Ordering Guide

An example order number for a vessel in this series is: 4740-75mL-FG-SS-A

A composite identification number to be used when ordering a 4740 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Vo	lumo	
-25mL	25 mL, 1" ID	
-25mL	75 mL, 1" ID	
	· · · ·	
B Gasket / N	laximum Temperature	
-T	PTFE Gasket, 350 °C	
-FG	Flexible Graphite Flat Gasket, 500-538 °C	
C Materials	of Construction	
-SS	T316 Stainless Steel	
-C20	Alloy 20	
-M0	Alloy 400	
-IN600	Alloy 600	
-IN625	Alloy 625	
-HB	Alloy B-2 / B-3	
-HC	Alloy C-276	
-NI	Nickel 200*	
-TI2	Titanium Grade 2*	
-TI4	Titanium Grade 4*	
-ZR702	Zirconium Grade 702*	
-ZR705	Zirconium Grade 705*	
* Maximum temperature and/or pressure limited See <u>page 10</u> or <u>24</u> for complete list of available alloys. D Head Configuration		
\sim	Rupture disc port with rupture disc	
-RD	assembly & disc	
-A	One head opening, "A" socket	
-NS	Non-standard, custom head opening	
Optional F (Check all	ittings & Custom Options	
-HCSC -TC	Alloy C-276 Screw Cap for HPHT	
-16	Thermocouple (in cylinder bottom)	
-4310A	Gage block assembly (for head with "A" socket opening)	
-4310ATC	Gage block assembly with thermocouple (for head with "A" socket opening)	
F Pressure Gage Option		
-No Symbol	None	
-10,000	10,000 psi / 690 bar	
-7500	7500 psi / 517 bar	
-5000	5000 psi / 345 bar	
-3000	3000 psi / 207 bar	
-3000		

Please note that all options and combinations are not compatible with all models.

2000 psi / 138 bar

Model 4740-75mL Vessel, with bench top stand and Ceramic Fiber Heater (-BTS).

-2000

G Stand	with Heater Options	
-No Symbol	None	
-BTS	Ceramic Fiber Heater w/Stand, 700W	
H Contro		
-No Symbol	None	
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one addi- tional display module.	
-4848	PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.	
See <u>Chapter (</u>	<u>6</u> for complete list of Controller options.	
Contro	ller Options (Check all needed)	
-PDM	Pressure Display Module	
-HTM	High Temperature Cut Off Module	
-ETLM	External Temperature Limit Module	
-SVM	Solenoid Valve Module (for cooling control)	
-A1925E4	RS-485 to USB Cable for 4838/4848 Controller (required for data logging)	
-A1925E6	RS-485 to USB Converter, isolated, 30-ft	
-A2208E	RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)	
-A3504HC	SpecView Software Package for 4838/4848 Controller	
Electric	cal Supply for Heater	
	Controller	
-No Symbol	None	
-115	115V Power	
-230	230V Power	
K Certific	ations	
-ASME	ASME Documentation	
-CE	CE Documentation	
-P	Parr Certification	

Please note that all options and combinations are not compatible with all models.

Model 4740-75mL-A Vessel, with 4310 Gage Block Assembly.




1 - 8 0 0 - 8 7 2 - 7 7 2 0 | 1 - 3 0 9 - 7 6 2 - 7 7 1 6



Series 4750 General Purpose Pressure Vessels, 125 & 200 mL

General Purpose 4750 Series Vessels are available in volumes of 125 and 200 mL, with temperatures up to 350 °C and maximum pressure of 3000 psi (200 bar). They come standard with a Split Ring closure with 6 compression bolts, and a choice of three different head styles: 1/8" NPT w/plug, an "A" Socket, or an "A" Socket with Thermocouple. The maximum amount of openings for this head style is two and would be dependent on the size of the ports and the type of fittings required. For the 4750 Series General Purpose Vessels, Parr offers the 4921 Ceramic Fiber Heater.



4750 Cross Section

Model 4750-200mL-A Vessel, with 4310 Gage **Block Assembly.**

Series 4750 Pressure Vessel Specifications			
Shade	ed bar indicates specifications that change wit	hin series.	
Model Number	4750		
Approximate Volume, mL	125	200	
Maximum Pressure (MAWP)	3000 psi (207 bar, 20	0 bar for CE orders)	
Maximum Temperature	350	°C	
Vessel Style			
Head Style	Move	eable	
Closure	Split-Ring (6 Cor	npression Bolts)	
Head Seal	Flat Gask	et - PTFE	
Head Opening Options	One head opening (1/8" NPT, 1/4" NPT, "A" socket, or rupture disc assembly port) Two head openings ("A" socket and 1/8" NPT or two 1/8" NPT) Custom port option (Contact our technical sales department to discuss options)		
Maximum Head Openings	2		
Temperature Measurement (optional)	Fixed Thermocouple (Thermowell for special alloys)		
Stand and Heater (optional)			
Stand Type	Bench Top Sta	nd with Heater	
Heater Style	Cerami	c Fiber	
Head Style	70	00	
Head Style	d Style 6 / 3		
Weight and Dimensions	Weight and Dimensions		
Cylinder I.D. x Depth, inches	1.5 x 4.4	1.5 x 6.9	
Vessel Assembly Weight, pounds*	4.4	5	
Cylinder Weight, pounds	1.7	2.4	
*Vessel weight is based on a moveable head vessel with a "A" socket and 1/8" NPT port. Vessel weight will vary based on head configuration choice and fitting options.			
Other options available. See Ordering Guide, visit <u>www.parrinst.com</u> , or call for more information.			

Type: **General Purpose**

Mounting Style: Moveable

Vessel Sizes, mL: 125 & 200

Standard Pressure MAWP Rating, psi (bar): 3000 (200)

Maximum Operating Temp., °C: 350

Series 4750 Ordering Guide



Model 4750-125mL-A Vessel, with 4310 Gage Block Assembly.



Available stand assembly with heater

An example order number for a vessel in this series is: 4750-125mL-T-SS-A-21AC4-A22AC3

A composite identification number to be used when ordering a 4750 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vesse	Volume
-125mL	125 mL, 1" ID
-200mL	200 mL, 1" ID
B Gaske	t / Maximum Temperature
-T	PTFE Gasket, 350 °C
-FG	Flexible Graphite Flat Gasket, 500-538 °C
C Mater	ials of Construction
-SS	T316 Stainless Steel
-C20	Alloy 20
-M0	Alloy 400
-IN600	Alloy 600
-IN625	Alloy 625
-HB	Alloy B-2 / B-3
-HC	Alloy C-276
-NI	Nickel 200
-TI2	Titanium Grade 2
-TI4	Titanium Grade 4
-ZR702	Zirconium Grade 702
-ZR705	Zirconium Grade 705
See <u>page 10</u> or <u>24</u> for complete list of available alloys.	

D Head Configuration		
-RD	Rupture disc port with rupture disc assembly & disc	
-1/8	One head opening, 1/8" NPT port, plugged	
-1/4	One head opening, 1/4" NPT port, plugged	
-A	One head opening, "A" socket	
-A &1/8	Two head openings, "A" socket & 1/8" NPT w/ thermocouple	
-(2)1/8	Two head openings, (2) 1/8" NPTF(s) w/ plug and w/ thermocouple	
-NS	Non-standard, custom head opening	

Optional Fittings & Custom Options (Check all needed) E Gage block assembly (for head with "A" -4310A socket opening)

Pressure Gage Option		
-TC	Thermocouple (requires 1/8" NPT opening)	

Gage Block Assembly Required)	
None	
3000 psi / 207 bar	
2000 psi / 138 bar	
1000 psi / 69 bar	
600 psi / 40 bar	
200 psi / 14 bar	
100 psi / 7 bar	

G Stand w	ith Heater Options
-No Symbol	None
-BTS	Bench Top Stand w/ Heater
H Controll	er
-No Symbol	None
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module.
-4848	PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.
0 01 1 01	

See <u>Chapter 6</u> for complete list of Controller options.

Control	er Options (Check all needed)
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E4	RS-485 to USB Cable for 4838/4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E	RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)
-A3504HC	SpecView Software Package for 4838/4848 Controller

J Electrical Supply for Heater and/or Controller				
-No Symbol	None			
-115	115V Power			
-230	230V Power			
K Certifications				
ACNAE				

-ASME	ASME Documentation	
-CE	CE Documentation	
-P	Parr Certification	

Please note that all options and combinations are not compatible with all models.

Series Number:

Type: Micro

Stand: Bench Top

Mounting Style: Moveable or Fixed Head

Sizes, mL: **25, 50, 100**

Standard Pressure MAWP Rating, psi (bar): **3000 (200)**

High Pressure (HP) MAWP Rating, psi (bar): **5000 (345)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

High Temperature (HT) / High Pressure (HP) Maximum Operating Temperature, °C: 500 w/ FG Flat Gasket (Fixed Head Only)

Series 4790 Pressure Vessels, 25-100 mL

These small pressure vessels


are a good choice for chemists working with very expensive materials or materials only available in small amounts. They will also appeal to users who wish to minimize the risks associated with hazardous materials or reactions by restricting the reactants or products to a minimum. Use of the Series 4790 Vessels also helps minimize the quantities of waste products which may require special disposal procedures. The 4790 Vessels are now available in high pressure (5000 psi) or high pressure/high temperature (500 °C) versions.

These micro vessels have been designed to provide as many of the features of the larger vessels as possible in the limited space available. The standard head fittings include a gage, rupture disc, thermocouple, and valves. Optional cooling is available with aluminum block heaters with cooling capability, welded jackets, or cold fingers.


The vessels are offered in both fixed head and moveable head vessel styles with choices of a self sealing O-ring for temperatures up to 225 °C or with FFKM O-ring for temperatures to 300 °C or with a flat PTFE flat gasket for temperatures to 350 °C, or with an FG Flexible Graphite gasket for temperatures to 500 °C with the HP/HT option (for fixed head versions only). A split ring closure is standard.

All three volumes use the same bench top ceramic fiber heater.

These micro vessels can be easily converted from one size to another by simply substituting a larger or smaller cylinder and the corresponding internal fittings. The support system can also be readily adapted to accept any of the vessels from the 4790 Series Vessels. The opportunity to modify these small reactors is restricted because of the limited head space available.



4790-25mL Cross Section with Flat Gasket Seal



4790-25mL Cross Section with O-ring Seal



4790-100mL-VGR-DVD Moveable Head Vessel with PTFE Flat Gasket

4790

Shaded bar indic	ates specifications that change	within series.		
Model Number		4790		
Approximate Volume, mL	25	50	100	
Maximum Pressure (MAWP)				
Standard Rating	30	00 psi (207 bar, 200 bar for CE or	ders)	
High Pressure (HP)		5000 psi (345 bar)		
High Pressure/High Temperature (HP/HT)		5000 psi (345 bar)		
Maximum Temperature				
with FKM 0-ring		225 °C		
with FFKM O-ring		300 °C		
with PTFE Flat Gasket		350 °C		
with Flexible Graphite Flat Gasket (HP/HT only)		500 °C		
Vessel Details				
Mounting Style	Moveable or	Fixed Head (HP/HT option requi	res Fixed Head)	
Closure	Split-Ring (6 Compressi	Split-Ring (6 Compression Bolts for Flat Gasket, no Compression Bolts for O-ring) or Screw Cap (O-ring seal only)		
Head Opening Options	Ports for valve	One head opening (1/8" NPT*, 1/4" NPT*, "A" socket, or rupture disc assembly port) Ports for valve, gage, rupture disc, and thermocouple (VGR**) Custom port option (Contact our technical sales department to discuss options)		
Maximum Head Openings	5 (Depen	dent on opening size, and requir	ed fittings)	
Valve Connections	1/8" NPT Male (standa	1/8" NPT Male (standard rating), 1/4" NPT Male (HP), 1/4" NPT Female (HP/HT)		
Pressure Gage, Size		3.5 inches		
Range, Standard Temperature	0-3000 psi (207 bar)			
Range, HP & HP/HT		0-5000 psi (345 bar)		
Temperature Measurement	Fixed Thermocouple (Thermowell for special alloys)			
Cooling (optional)		Cold Finger		
* NPT ports are not available for vess	els with a High Temperature (H1) rating. ** 4790HPHT requires \	/GR.	
Stand and Heater (optional)				
Stand Type		Only (standard rated and HP ves op Stand with Heater (required		
Heater Style (Standard rating & HP)		Ceramic Fiber Heater		
Heater Power, Watts	700	700	350	
Maximum Load, amps, 115V / 230V	6 / 3	6 / 3	3/2	
Heater Style (HP/HT)		Ceramic Fiber Heater		
Heater Power, Watts		700		
Maximum Load, amps, 115V / 230V	6/3			
Weight and Dimensions				
Cylinder I.D. x Depth, inches	1.0 x 2.0	1.3 x 2.3	1.3 x 4.6	
Vessel Assembly w/ Moveable Head, Weight, pounds**	16	15	16	
Vessel Assembly w/ Fixed Head, Weight, pounds**	17	17	18	
Cylinder Weight, pounds 3.5 3.1 4.2		4.2		

Weights shown are based on vessels with valve, gage, rupture disc, thermocouple, and double valve assembly.

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information.

8

Series Number:

Index to Options

Option	Page
Certification	8
Materials of Construction	n 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

150

Series 4790 Pressure Vessel Systems



Series 4790 Ordering Guide

An example order number for a vessel in this series is: **4790-25mL-T-SS-3000**

A composite identification number to be used when ordering a 4790 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Vo	olume		
-No Symbol	3000 psi / 207 bar, 200 bar for CE orders @ 350 °C		
-HP	5000 psi / 345 bar @ 350 °C		
-HPHT	5000 psi / 345 bar @ 350 °C		
B Vessel Vo	alumo		
-25ml			
	25 mL, 1" ID		
-50mL	50 mL, 1.3" ID		
-100mL	100 mL, 1.3" ID		
C Head Mo	unting Style		
-No Symbol	Moveable Head (Standard configuration)		
-FH	Fixed Head (Required for HP/HT rating)		
D Gasket /	Maximum Temperature		
-T	PTFE Flat Gasket, 350 °C		
-0V	FKM 0-ring, 225 °C		
-0K	FFKM 0-ring, 300 °C		
-FG	Flexible Graphite Flat Gasket, 500 °C (HPHT option only)		
	of Construction		
-SS	T316 Stainless Steel		
-C20	Alloy 20*		
-M0	Alloy 400*		
-IN600	Alloy 600		
-IN625	Alloy 625		
-HB	Alloy B-2 / B-3*		
-HC	Alloy C-276		
-NI	Nickel 200		
-TI2	Titanium Grade 2*		
-TI4	Titanium Grade 4*		
-ZR702	Zirconium Grade 702*		
-ZR705	Zirconium Grade 705*		
	ure option not available. 2 <u>4</u> for complete list of available alloys.		
F Head Cor	nfiguration		
-RD	Rupture disc port with rupture disc assembly & disc		
-1/8	One head opening, 1/8" NPT port, plugged		
-1/4	One head opening, 1/4" NPT port, plugged		
-A	One head opening, "A" socket		
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple (includes thermowell for special alloys)		
-NS	Non-standard, custom head opening		
G Ontional	Fittings & Custom Ontions (Check all peeded)		
-VD	Fittings & Custom Options (Check all needed) Single Valve on head with dip tube (VGR Required)		
-DVD	Double Valve Assembly on head with dip tube (VGR Required)		
-CF	Cold Finger		
-4310A	Gage Block Assembly (for head with "A" socket opening)		
-4310A -CAD	Internal Catalyst Addition Device		
-XCAD	External Catalyst Addition Device		
-SCP	Solids Charging Port (Ball Valve)		
-RC	Reflux Condenser		
-RTC	Reflux/Take-Off Condenser		
-WJ	Welded Jacket		

H Pressu	
-No Symbol	Single Valve on head with dip tube (VGR Required)
-5000	5000 psi / 345 bar
-3000 -2000	3000 psi / 207 bar 2000 psi / 129 bar
-2000	2000 psi / 138 bar 1000 psi / 69 bar
-600	1000 psi / 69 bar 600 psi / 40 bar
-200	200 psi / 14 bar
-100	100 psi / 7 bar
	and Heater Options
-No Symbol	None
-NS	Heater only, no stand (Temperatures up to 350 °C)
-BTS	Bench Top Stand with Heater
	'G" for Welded Jacket option
	·
J Contro	ller
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module.
-4848	PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.
See <u>Chapter</u>	<u>6</u> for complete list of Controller options.
K Contro	ller Options (Check all needed)
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E4	RS-485 to USB Cable for 4838/4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E	RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)
-A3504HC	SpecView Software Package for 4838/4848 Controller
	cal Supply for Heater Controller
-No Symbol	None
-115	115V Power
-230	230V Power
M Certifi	cations
-No Symbol	No Certification
-ASME	ASME Documentation
-CE	CE Documentation
-P	Parr Certification

Please note that all options and combinations are not compatible with all models.

X



Type: Mini

Stand: Bench Top

Mounting Style: Moveable or Fixed Head

Sizes, mL: **100-600**

Standard Pressure MAWP Rating, psi (bar): **3000 (200)**

High Pressure (HP) MAWP Rating, psi (bar): **2000 (138)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

High Temperature (HT) / Maximum Operating Temperature, °C: 500 w/ FG Flat Gasket

Series 4760 Pressure Vessels, 100-600 mL

These are the most popular of our small general purpose

vessels. Parr offers a range of sizes large enough to work with significant sample sizes, yet small enough to be handled with ease by all operators.

They are made in both fixed head and moveable vessel styles and are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

Choosing the high temperature option (HT) boosts the maximum temperature to 500 °C, but reduces the MAWP to 2000 psi (138 bar). Standard temperature vessels can be converted to high temperature vessels (500 °C max temperature and 2000 psi MAWP) by changing the head assembly (contains cone connections, high temperature valves, and Grafoil® gasket) replacing the heater with a ceramic fiber heater and replacing the split ring. Contact Parr for details. Although internal and external

space is limited on these small vessels, custom ports and fittings are available. All vessels in this series can also be easily


 $\rm 4760\text{-}160mL\text{-}VGR\text{-}DVD$ Vessel, BTS Support Stand w/ Heater, and 4838 Controller.

converted from one size to another by simply substituting a longer or shorter cylinder with the corresponding internal fittings and heater.


Clamp-on or mantle heaters are available for the Series 4760 moveable vessels. Stands with heaters are offered for the fixed head vessels.

The heaters furnished for nonstirred vessels are the same as those furnished for the stirred versions of the same size vessels. A full description of the different types of heaters is found in our Heater Options section. Automatic temperature controllers are available for most Parr non-stirred vessels.

Custom heaters are also available for these vessels when it is necessary to meet explosion proof requirements, to accommodate installed windows or fittings, or to convert to a welded jacket for steam or oil heating.



47600 Cross Section with O-ring Seal



4760 Cross Section with Flat Gasket Seal

4760-300mL with plug and rupture disc assembly installed.

4760

8

Shaded bar indicates specifications that change v	vitnin series.				
Model Number			4760		
Approximate Volume, mL	100	160	300	450	600
Maximum Pressure (MAWP)	3000 psi (207 bar, 200 bar for CE orders)				
Maximum Pressure for HT (MAWP)	N/	Ά		2000 psi (138 bar)	
Maximum Temperature					
with FKM 0-ring	225 °C				
with FFKM 0-ring			300 °C		
with PTFE Flat Gasket			350 °C		
with Grafoil Gasket (HT)	N/	Ά		500 °C	
Vessel Details					
Mounting Style		Ν	Noveable or Fixed Hea	ad	
Closure	Split-Rir	ng (6 Compression Bo	lts for Flat Gasket, no	Compression Bolts for	O-ring)
Head Opening Mounts	One head opening (1/8" NPT*, 1/4" NPT*, "A" socket, or rupture disc assembly port) Ports for valve, gage, rupture disc, and thermocouple (VGR) Custom port option (Contact our technical sales department to discuss options)				
Maximum Head Openings		7 (Dependent o	on opening size, and r	equired fittings)	
Valve Connections	1/8" NPT Male, 1/8" NPT Female for HT				
Pressure Gage, Size			3.5 inches		
Range, Standard Temperature			0-3000 psi (207 bar)		
Range, High Temperature	N/	'A		0 - 2000 psi (138 bar)	
Temperature Measurement		Fixed Thermoc	ouple (Thermowell fo	r special alloys)	
Cooling Coil (optional)	N/	Ά		Single Loop	
Bottom Drain Valve (BDV) (optional)	1/4" NPT (Not available on HT Models)				
* NPT por	ts are not available f	or vessels with a High	n Temperature (HT) ra	ting.	
Stand and Heater (optional)					
Cooling Coil (optional)	Bench Top S		(standard rated vess ndard rated vessels a	els only) and nd high temperature ra	ated vessels)
Heater Style (Standard Temperature)	Clamp-c	on Band		Rigid Mantle	
Heater Power, Watts	52	25	510	590	780
Maximum Load, amps (115V / 230V), Nominal	5 /	3	5/3	5/3	7 / 4
Heater Style (High Temperature)				Ceramic Fiber	
Heater Power, Watts	N/	A	800	800	1100
Maximum Load, amps (115V / 230V), Nominal			7 / 4	7 / 4	10 / 5
Weight & Dimensions					
Cylinder I.D. x Depth, inches	2.0 x 2.0	2.5 x 2.0	2.5 x 4.0	2.5 x 6.0	2.5 x 8.0
Vessel Assembly w/ Moveable Head Weight, pounds**	14	13	14	15	17
Vessel Assembly w/ Moveable Head Weight, pounds**	15	14	15	17	18
Cylinder Weight, pounds	3.3	2.4	3.7	4.9	6.2
Weights shown are based	on vessels with valve	sed on head configur a, gage, rupture disc, ide, visit <u>www.parrin</u>	· · ·	ouble valve assembly.	

Series Number:

Index to Options

=	
Option	Page
Certification	8
Materials of Construction	n 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fittings	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4760 Pressure Vessel Systems



4760-300mL, Fixed Head, Quick-release Split Ring, and a 4838 Temperature Controller.

Series 4760 Ordering Guide

An example order number for a vessel in this series is: 4760-100mL-T-SS-VGR-3000

A composite identification number to be used when ordering a 4760 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Rati	na
-No Symbol	3000 psi, 200 bar @ 350 °C
-HT	2000 psi, 138 bar @ 500 °C (300, 450, or 600 mL only)
	s are 970 mL & 1900 mL
B Vessel Volu	
-100mL	100 mL, 2.0" ID
-160 mL	160 mL, 2.5" ID
-300 mL	300 mL, 2.5" ID
-450 mL	450 mL, 2.5" ID
-600 mL	600 mL, 2.5" ID
C Head Moun	ting Style
-No Symbol	Moveable Head (Standard configuration)
-FH	Fixed Head
D Gasket / Ma	
-T	aximum Temperature
-0V	PTFE Flat Gasket, 350 °C
-	FKM 0-ring, 225 °C
-0K -FG	FFKM 0-ring, 300 °C
-FU	Flexible Graphite, Flat Gasket, 500 °C (HT option only)
E Vessel Mat	erial of Construction
-SS	T316 Stainless Steel
-C20	Alloy 20*
-M0	Alloy 400*
-IN600	Alloy 600
-IN625	Alloy 625
-HB	Alloy B-2 / B-3*
-HC	Alloy C-276
-NI	Nickel 200
-TI2	Titanium Grade 2*
-TI4	Titanium Grade 4*
-ZR702	Zirconium Grade 702*
-ZR705	Zirconium Grade 705*
*High temperature	option not available.
See <u>page 10</u> or <u>24</u> f	for complete list of available alloys.
F Head Confi	
-RD	Rupture disc port with rupture disc assembly & disc
-1/8	One head opening, 1/8" NPT port, plugged
-1/8	One head opening, 1/3 NPT port, plugged
-1/4 -A	One head opening, "A" socket
-A	
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple (includes thermowell for special alloys)
-NS	Non-standard, custom head opening
G Optional Fi	ttings & Custom Options (List all needed)
-VD	Single valve on head w/ dip tube (VGR Required)
-DVD	Double valve assembly on head w/ dip tube (VGR Required)
-CL	Cooling Loop
-BDV	Bottom drain valve
-4310A	Gage block assembly (for head with "A" socket opening)
-4310A -CAD	Internal Catalyst Addition Device
-XCAD	External Catalyst Addition Device
-SCP	Solids Charging Port (Ball Valve)
-SCF -RC	Reflux Condenser
-nc -RTC	Reflux/Take-Off Condenser
-WJ	Welded Jacket

	e Gage Option nd Gage Block only)
-No Symbol	None
-3000	3000 psi / 207 bar
-2000	2000 psi / 138 bar
-1000	1000 psi / 69 bar
-600	600 psi / 40 bar
-200	200 psi / 14 bar
-100	100 psi / 7 bar
I Stand a	nd Heater Options
-No Symbol	None
-NS	Heater only, no stand (Temperatures up to 350 °C)
-BTS	Bench top stand with heater
See Section "	G″ for Welded Jacket option
J Control	ler
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module.
-4848	PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.
See <u>Chapter 6</u>	for complete list of Controller options.
\sim	ler Options
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-SVM	Solenoid Valve Module (for cooling control)
-A1925E4	RS-485 to USB Cable for 4838/4848 Controller (required for data logging)
-A1925E6	RS-485 to USB Converter, isolated, 30-ft
-A2208E	RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)
-A3504HC	SpecView Software Package for 4838/4848 Controller
L Electric	al Supply for Heater and/or Controller
-No Symbol	None
-115	115V
-230	230V
M Certific	ations
-ASME	ASME Documentation
-CE	CE Documentation
-6E	



Туре:

General Purpose

Stand: Bench Top

Mounting Style: Moveable or Fixed Head


Sizes, mL: 1000 and 2000 970 and 1900 HP

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

High Pressure (HP) MAWP Rating, psi (bar): **2900 (200)**

Maximum Operating Temperature, °C at High Pressure (HP): **350** @ **2900 psi**



4600-1L Cross Section

The Series 4600, 1.0 and 2.0 L Vessels are the largest of the Parr Vessels that can be handled on a bench top.

Series 4600 Pressure Vessel Systems, 1 & 2 Liter

These reactors are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to

300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum. An optional HP (High Pressure) flat gasket version has been added for maximum allowable working pressure of 2900 psi (200 bar) at 350 °C. Both fixed head and moveable vessel designs are available.

An optional pneumatic lift is available for the heavier 2 liter cylinder and heater. It should be noted, however, that the 2 liter, fixed head model is tall and may not be convenient to operate on a standard height bench top. We recommend using a floor stand.

With their larger diameter, these mid-size reactors have sufficient space for special modifications, such as an internal cooling coil, bottom drain valve (not on 2 liter bench top), ball valve for a solids charging port, catalyst addition devices, condensers, electrical feed-throughs and more.



4600-1L with Rupture Disc Assembly

4600-1L Pressure Vessel with Valve, Gage, Rupture Disc, and Thermowell (VGR) and single inlet valve (-VD)

4600

8

Shaded bar indicates specifications that change within	series.			
Model Number		46	00	
Approximate Volume, mL	1000	2000	970	1900
Maximum Pressure (MAWP)	1900 psi	(131 bar)	2900 ps	i (200 bar)
Maximum Temperature				
with FKM 0-ring	225	o°C		1/A
with FFKM 0-ring	300	٥°C	ľ	J/A
with PTFE Flat Gasket		350	0° (
with Flexible Graphite Flat Gasket		350	0° (
Vessel Details				
Mounting Style		Moveable o	r Fixed Head	
Closure	Split-Ring (6 Compression No Compression	on Bolts for Flat Gasket, Bolts for O-ring)	Split-Ring with 12	Compression Bolts
Head Opening Options	Ports for valv	ing (1/8" NPT, 1/4" NPT, "I ve, gage, rupture disc, an ption (Contact our technic	d thermowell with thermo	ocouple (VGR)
Maximum Head Openings		7 (Dependent on opening	size, and required fitting	s)
Valve Connections		1/8" Ma	ale NPT	
Pressure Gage, Size		4.5" Di	ameter	
Range	0-2000 ps	i (138 bar)	0-3000 p	si (207 bar)
Temperature Measurement		Thern	nowell	
Cooling Coil (Optional)		Serpe	entine	
Bottom Drain Valve (optional, requires stand)	3/8"	NPT	1/4	I" NP
Stand and Heater (optional)				
Stand Type*	(Also Available: Be	Bench Top Heate Inch Top Stand, Floor Star	er with Base Plate nd, Floor Stand w/ Pneun	natic Lift, Floor Cart)
Heater Style, Bench Top Heater		Cal	rod	
Heater Power, Watts (115V / 230V)	1000	1350 / 1500	1000	1350 / 1500
Maximum Load, amps, (115V / 230V)	8 / 4	11 / 6	8 / 4	11/6
*Heater watts and amps vary between d	ifferent stand with heater	styles, contact Parr Tech	nical Support for more in	nformation.
Weights & Dimensions				
Cylinder I.D. x Depth, inches	4.0 x 5.4	4.0 x 10.5	3.75 x 5.4	3.75 x 10.5
Vessel Assembly w/Moveable Head Weight, pounds**	22	27	26	34
Vessel Assembly w/Fixed Head Weight, pounds**	25	30	28	35
Cylinder Weight, pounds	8.9	14	11	18

Other options available. See Ordering Guide, visit <u>www.parrinst.com</u>, or call for more information.

Series Number:

Index to Options

Option	Page
Certification	8
Materials of Construction	n 10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	118
Rupture Disc	124
Pressure Gage	125
Cooling Coil	129
Bottom Drain Valve	131
4848 Reactor Controller	102
4871 Process Controller	109
A2110E Motor Controller	114

Custom Option	Page
Stirrers	120
Heaters	116
External Valves & Fitting	s 138
Solids Charging Ports	128
Condensers	123
Electrical Glands	137
Catalyst Addition Device	128
Catalyst Baskets	122
Relief Valves	125
Windows	136
Gear Drives	119
Explosion Proof Options	134
Liquid Sampling	130

Accessories	Page
Liners	129
Spare Parts Kits	137
Pressure Hoses	133
Check Valves	125
Liquid Pipettes	127
Gas Burettes	126

Series 4600 Pressure Vessel Systems



4600-2L-VGR Vessel with bench top heater (-BTH) and base plate

8

Series 4600 Ordering Guide

An example order number for a vessel in this series is: 4600-1L-T-SS-VGR-2000

A composite identification number to be used when ordering a 4600 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Rati	na
- No Symbol	Standard rating, 1900 psi at 350 °C (1 & 2 Liter Vessels)
No bymbol	High pressure rating, 2900 psi at 350 °C
-HP	(970 & 1900 mL Vessels)
B Vessel Volu	
-1000mL	1000 mL, 4" ID (standard rating only)
-2000mL	2000 mL, 4" ID (standard rating only)
-970mL	970 mL, 3.75" ID (high pressure rating only)
-1900mL	1900 mL, 3.75" ID (high pressure rating only)
C Head Moun	ting Style
-No Symbol	Moveable Head (Standard configuration)
-FH	Fixed Head
D Gasket / Ma	aximum Temperature
-T	PTFE Flat Gasket, 350 °C
-0V	FKM O-ring, 225 °C
-0K	FFKM O-ring, 300 °C
-FG	Flexible Graphite Flat Gasket, 350 °C
	evial of Construction
-SS	erial of Construction T316 Stainless Steel
-SS -C20	
-020 -M0	Alloy 20
	Alloy 400
-IN600	Alloy 600
-IN625 -HB	Alloy 625 Alloy B-2 / B-3
-нь -нс	, .
-NI	Alloy C-276 Nickel 200
-INI -TI2	Titanium Grade 2
-TI2 -TI4	Titanium Grade 4
-ZR702	Zirconium Grade 702
-ZR705	Zirconium Grade 705
	for complete list of available alloys.
~	
F Head Confi	
-RD	Rupture disc port with rupture disc assembly & disc
-1/8	One head opening, 1/8" NPT port, plugged
-1/4	One head opening, 1/4" NPT port, plugged
-В	One head opening, "B" socket
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple w/thermowell
-NS	Non-standard, custom head opening
G Optional Fi	ttings & Custom Options (List all needed)
-VD	Single valve on head w/ dip tube (VGR Required)
-DVD	Double valve assembly on head w/ dip tube (VGR Required)
-SC	Serpentine cooling coil
-BDV	Bottom drain valve
-4310B	Gage block assembly (for head with "A" socket opening)
-CAD	Internal Catalyst Addition Device
-XCAD	External Catalyst Addition Device
-SCP	Solids Charging Port (Ball Valve)
-RC	Reflux Condenser
-RTC	Reflux/Take-Off Condenser
-WJ	Welded Jacket

-No SymbolNone-30003000 psi / 207 bar-20002000 psi / 138 bar-10001000 psi / 69 bar-600600 psi / 40 bar-200200 psi / 14 bar-100100 psi / 7 bar 1 Stand and Heater Options-No SymbolNone-BTHBench top heater with base plate (moveable head only)-FSBench top stand (moveable or fixed head)-FCFloor cart (moveable head only)-FSPFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)-FSPFloor stand with pneumatic lift (fixed head only)-FSPFloor stand with pneumatic lift (fixed head only)-A838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modulePDMPressure Display Module-HTMHigh Temperature Cut Off Module-HTMHigh Temperature Cut Off Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter, isolated-A3504HCSpecView Software Package-No SymbolNone-115115V Heater-230230V Heat		e Gage Option nd Gage Block only)
-30003000 psi / 207 bar-20002000 psi / 138 bar-10001000 psi / 69 bar-600600 psi / 40 bar-200200 psi / 14 bar-100100 psi / 7 bar1Stand and Heater Options-No SymbolNone-BTHBench top heater with base plate (moveable head only)-BTSBench top stand (moveable or fixed head)-FCFloor cart (moveable head only)-FSFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)-SSPFloor stand with pneumatic lift (fixed head only)-A838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display 	-No Symbol	None
-10001000 psi / 69 bar-600600 psi / 40 bar-200200 psi / 14 bar-100100 psi / 7 barI Stand and Heater Options-No SymbolNone-BTHBench top heater with base plate (moveable head only)-BTSBench top stand (moveable or fixed head)-FCFloor cart (moveable or fixed head)-FSFloor stand (moveable or fixed head)-FSFloor stand (moveable or fixed head)-FSPFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)-SSPFloor stand with pneumatic lift (fixed head only)-A838Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter & for complete list of Controller options.•KController Options -PDM-PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software PackageLetectrical Supply for Heater and/or Controller-No SymbolNone -115-115115V Heater	,	3000 psi / 207 bar
-600600 psi / 40 bar-200200 psi / 14 bar-100100 psi / 7 barI Stand and Heater Options-No SymbolNone-BTHBench top heater with base plate (moveable head only)-BTSBench top stand (moveable or fixed head)-FCFloor cart (moveable head only)-FSFloor stand (moveable or fixed head)-FSPFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)-SSFloor stand with pneumatic lift (fixed head only)-Vo SymbolNone-No SymbolNone-No SymbolNone-R838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options -PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software Package-No SymbolNone-115115V Heater	-2000	2000 psi / 138 bar
200200 psi / 14 bar-100100 psi / 7 barI Stand and Heater Options-No SymbolNone-BTHBench top heater with base plate (moveable head only)-BTSBench top stand (moveable or fixed head)-FCFloor cart (moveable or fixed head)-FSFloor stand (moveable or fixed head)-FSPFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)J Controller-No SymbolNone-No SymbolNonePID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module4848Sec Chapter 6 for complete list of Controller options.KController Options-PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter, isolated-A3504HCSpecView Software Package-No Symb	-1000	1000 psi / 69 bar
-100 100 psi / 7 bar I Stand and Heater Options -No Symbol None -BTH Bench top heater with base plate (moveable head only) -BTS Bench top stand (moveable or fixed head) -FC Floor cart (moveable head only) -FS Floor stand (moveable or fixed head) -FSP Floor stand with pneumatic lift (fixed head only) J Controller -No Symbol -No Symbol None -4838 PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module. See Chapter 6 for complete list of Controller options. K Controller Options -PDM -PDM Pressure Display Module -HTM High Temperature Cut Off Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package -No Symbol None	-600	600 psi / 40 bar
Stand and Heater Options -No Symbol None -BTH Bench top heater with base plate (moveable head only) -BTS Bench top stand (moveable or fixed head) -FC Floor cart (moveable head only) -FS Floor stand (moveable or fixed head) -FSP Floor stand with pneumatic lift (fixed head only) J Controller -No Symbol None PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. -4838 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module. See Chapter 6 for complete list of Controller options. K Controller Options -PDM -PTM Pressure Display Module -HTM High Temperature Cut Off Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecVi	-200	200 psi / 14 bar
-No SymbolNone-BTHBench top heater with base plate (moveable head only)-BTSBench top stand (moveable or fixed head)-FCFloor cart (moveable head only)-FSFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)JController-No SymbolNone-4838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter & for complete list of Controller options.KController Options -PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software Package-No SymbolNone-115115V Heater	-100	100 psi / 7 bar
Bench top heater with base plate (moveable head only)-BTKBench top stand (moveable or fixed head)-FCFloor cart (moveable head only)-FSFloor stand (moveable or fixed head)-FSFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)JController-No SymbolNone-4838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options -PDM-PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software PackageLectrical Supply for Heater and/or Controller-No SymbolNone-115115V Heater		nd Heater Options
-BTR (moveable head only) -BTS Bench top stand (moveable or fixed head) -FC Floor cart (moveable or fixed head) -FS Floor stand (moveable or fixed head) -FSP Floor stand with pneumatic lift (fixed head only) J Controller -No Symbol None -4838 PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules. See Chapter 6 for complete list of Controller options. K Controller Options -PDM -PTM Pressure Display Module -HTM High Temperature Cut Off Module -ETLM External Temperature Limit Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package • Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-No Symbol	None
-FCFloor cart (moveable head only)-FSFloor stand (moveable or fixed head)-FSPFloor stand with pneumatic lift (fixed head only)JController-No SymbolNone-4838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options -PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter, isolated-A3504HCSpecView Software Package• Deterical Supply for Heater and/or Controller-No SymbolNone-115115V Heater	-BTH	
-FS Floor stand (moveable or fixed head) -FSP Floor stand with pneumatic lift (fixed head only) J Controller -No Symbol None -4838 PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules. See Chapter 6 for complete list of Controller options. K Controller Options -PDM -PDM Pressure Display Module -HTM High Temperature Cut Off Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package • Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-BTS	Bench top stand (moveable or fixed head)
-FSP Floor stand with pneumatic lift (fixed head only) • FSP Floor stand with pneumatic lift (fixed head only) • Controller -No Symbol -No Symbol None -4838 PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules. See Chapter 6 for complete list of Controller options. (K) Controller Options -PDM Pressure Display Module -HTM High Temperature Cut Off Module -ETLM External Temperature Limit Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package • Lectrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-FC	Floor cart (moveable head only)
Image: Second Symbol None -No Symbol None -4838 PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with one additional display module. -4848 PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules. See Chapter 6 for complete list of Controller options. K Controller Options -PDM -Pressure Display Module -HTM -High Temperature Cut Off Module -SVM -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package • Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-FS	Floor stand (moveable or fixed head)
-No SymbolNone-4838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options -PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software Package• Electrical Supply for Heater and/or Controller -No SymbolNone-115115V Heater	-FSP	Floor stand with pneumatic lift (fixed head only)
-No SymbolNone-4838PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options -PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software Package• Electrical Supply for Heater and/or Controller -No SymbolNone-115115V Heater	J Control	ler
-4838Data logging with Software (RS-485 to USB cable not included). For use with one additional display module4848PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options-PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software PackageLElectrical Supply for Heater and/or Controller-No SymbolNone-115115V Heater	-No Symbol	None
-4848logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.See Chapter 6 for complete list of Controller options.KController Options-PDMPressure Display Module-HTMHigh Temperature Cut Off Module-ETLMExternal Temperature Limit Module-SVMSolenoid Valve Module (for cooling control)-A1925E4RS-485 to USB Converter (required for data logging)-A1925E6RS-485 to USB Converter, isolated-A3504HCSpecView Software PackageLElectrical Supply for Heater and/or Controller-No SymbolNone-115115V Heater	-4838	Data logging with Software (RS-485 to USB cable not included). For use with one additional display
K Controller Options -PDM Pressure Display Module -HTM High Temperature Cut Off Module -ETLM External Temperature Limit Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package L Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-4848	logging with Software. (RS-485 to USB cable not included). For use with up to three additional display
-PDM Pressure Display Module -HTM High Temperature Cut Off Module -ETLM External Temperature Limit Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package • Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	See <u>Chapter 6</u>	for complete list of Controller options.
-HTM High Temperature Cut Off Module -ETLM External Temperature Limit Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package L Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	K Control	ler Options
-ETLM External Temperature Limit Module -SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package • Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-PDM	Pressure Display Module
SVM Solenoid Valve Module (for cooling control) -A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package L Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-HTM	High Temperature Cut Off Module
-A1925E4 RS-485 to USB Converter (required for data logging) -A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package L Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-ETLM	External Temperature Limit Module
-A1925E6 RS-485 to USB Converter, isolated -A3504HC SpecView Software Package L Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-SVM	Solenoid Valve Module (for cooling control)
-A3504HC SpecView Software Package Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-A1925E4	RS-485 to USB Converter (required for data logging)
Electrical Supply for Heater and/or Controller -No Symbol None -115 115V Heater	-A1925E6	RS-485 to USB Converter, isolated
-No Symbol None -115 115V Heater	-A3504HC	SpecView Software Package
-115 115V Heater	L Electric	al Supply for Heater and/or Controller
	-No Symbol	None
-230 230V Heater	-115	115V Heater
	-230	230V Heater
M Certifications	M Certifica	ations
-ASME ASME Documentation		
-CE CE Documentation	-CE	CE Documentation
-P Parr Certification	-P	Parr Certification



Type: High Pressure

Stand: Floor Stand or Bench Top

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: 600 and 1200

High Pressure MAWP Rating, psi (bar): **5000 (345)**

Standard Maximum Operating Temperature, °C: **350** w/ PTFE Flat Gasket



Model 4625-1200mL-VGR-DVD Vessel



4625 Cross Section

This series of vessels has

been designed for users who need higher operating pressures than the 2000 to 3000 psi offered by the General Purpose Reactors, but do not require the high operating temperatures provided by the

Series 4625 High Pressure Vessel Systems, 600 & 1200 mL

Series 4650 High Pressure, High Temperature Reactors. These reactors offer working pressures to 5000 psi (345 bar) at temperatures to 350 °C.

These vessels have been designed with outside dimensions comparable to the Series 4600 so that they can use the same support system and heater as these popular general purpose models. This provides not only an attractively priced high pressure/moderate temperature system, but also reactors that can be interchanged with the 1 and 2 liter sizes. The thicker walls required for higher operating pressures reduce the volumes of these reactors to 600 and 1200 mL.

These vessels can be used in either the bench top, floor stand, or moveable cart mountings. While the 1200 mL reactor is offered as a fixed head bench top model, it is too tall and too heavy to be handled comfortably on a standard height bench top. It is recommended that a floor stand support option should be selected unless the user has an adjustable bench top which will accommodate the overall height of the 1200 mL systems.

Series 4625 Pressure Vessel Specifications		
Shaded bar indi	cates specifications that change within s	eries.
Model Number	46	25
Approximate Volume, mL	600	1200
Maximum Pressure (MAWP)	5000 psi	(345 bar)
Maximum Temperature		
with PTFE Flat Gasket	350	0° (
with Flexible Graphite Flat Gasket	350	℃ (
Vessel Details		
Mounting Style	Moveable o	r Fixed Head
Closure	Split-Ring with 8	Compression Bolts
Head Opening Options	asseml Ports for valve, gage, rupture disc, an	"NPT, "B" socket, or rupture disc bly port) d thermowell with thermocouple (VGR) cal sales department to discuss options)
Maximum Head Openings	7 (Dependent on opening	size, and required fittings)
Valve Connections	1/4" Ma	ale NPT
Pressure Gage, Size	4.5" Di	ameter
Range	0-5000 ps	i (345 bar)
Temperature Measurement	Thern	nowell
Cooling Coil (optional)	Serpentir	ne or Loop
Bottom Drain Valve (BDV) (optional)	1" NPS (Floor Stand o	r Moveable Cart Only)
Stand and Heater (optional)		
Stand Type*	Bench Top Heater with Base Plate Floor Stand, Floor Stand w	e (Also Available: Bench Top Stand, / Pneumatic Lift, Floor Cart)
Heater Style, Bench Top Heater	Са	Irod
Heater Power, Watts	1000	1350 / 1500
Maximum Load, amps, 115V / 230V	8 / 4	11 / 6
*Heater watts and amps vary between different	stand with heater styles, contact Parr Te	echnical Support for more information.
Weight & Dimensions		
Cylinder I.D. x Depth, inches	3.25 x 4.7	3.25 x 9.8
Vessel Assembly w/ Moveable Head, Weight, pounds**	47	58
Vessel Assembly w/ Fixed Head, Weight, pounds**	50	62
Cylinder Weight, pounds	18	29
	ht is based on a vessel with VGR head fitt based on head configuration choice and	
Other options available. See Orde	ering Guide, visit <u>www.parrinst.com</u> , or ca	all for more information.

Series 4625 Ordering Guide

An example order number for a vessel in this series is: 4625-600mL-T-SS-1/8-21AC4-A22AC3

A composite identification number to be used when ordering a 4625 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Vol	ume
-600mL	600 mL, 3.25" ID
-1200mL	1200 mL. 3.25"
B Head Mou	
-No Symbol	Moveable Head (Standard configuration)
-FH	Fixed Head
C Gasket / M	laximum Temperature
-T	PTFE Flat Gasket, 350 °C
-FG	Flexible Graphite Flat Gasket, 350 °C
D Material of	f Construction
-SS	T316 Stainless Steel
-C20	Alloy 20
-020 -M0	Alloy 400
-IN600	Alloy 600
	,
-IN625	Alloy 625
-HB	Alloy B-2 / B-3
-HC	Alloy C-276
-NI	Nickel 200
-TI2	Titanium Grade 2
-TI4	Titanium Grade 4
-ZR702	Zirconium Grade 702
-ZR705	Zirconium Grade 705
See <u>page 10</u> or <u>24</u>	for complete list of available alloys.
E Head Conf	iguration
-RD	Rupture disc port with rupture disc assembly & disc
-1/8	One head opening, 1/8" NPT port, plugged
-1/4	One head opening, 1/4" NPT port, plugged
-В	One head opening, "B" socket
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple and thermowell
-NS	Non-standard, custom head opening
	Non-standard, custom nead opening
F Optional F	ittings & Custom Options (List all needed)
F Optional F	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required)
F Optional F -VD -DVD	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required)
F Optional F -VD -DVD -SC	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil
Optional F -VD -DVD -SC -BDV	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve
Optional F -VD -DVD -SC -BDV -4310B	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve Gage Block Assembly (for head with "A" socket opening)
Optional F -VD -DVD -SC -BDV -4310B -CAD	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve Gage Block Assembly (for head with "A" socket opening) Internal Catalyst Addition Device
Optional F -VD -DVD -SC -BDV -4310B -CAD -XCAD	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve Gage Block Assembly (for head with "A" socket opening) Internal Catalyst Addition Device External Catalyst Addition Device
Optional F -VD -DVD -SC -BDV -4310B -CAD -XCAD -SCP	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve Gage Block Assembly (for head with "A" socket opening) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve)
Coptional F -VD -DVD -SC -BDV -4310B -CAD -XCAD -SCP -RC	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve Gage Block Assembly (for head with "A" socket opening) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser
Optional F -VD -DVD -SC -BDV -4310B -CAD -XCAD -SCP	ittings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve Gage Block Assembly (for head with "A" socket opening) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve)

(VGR re	quired)
-No Symbol	None
-10000	10000 psi / 690 bar
-7500	7500 psi / 517 bar
-5000	5000 psi / 345 bar
-3000	3000 psi / 207 bar

H Stand a	nd Heater Options
-No Symbol	None
-BTH	Bench top heater with base plate (moveable head only)
-BTS	Bench top stand (moveable or fixed head)
-FC	Floor cart (moveable head only)
-FS	Floor stand (moveable or fixed head)
-FSP	Floor stand with pneumatic lift (fixed head only)
Control No Symbol	None
	PID Control, Ramp & Soak Programming, and
-4838	Data logging with Software (RS-485 to USB cable not included). For use with one additional display module
-4838	

See <u>Chapter 6</u> for complete list of Controller options.

J Controller Options		
-PDM	Pressure Display Module	
-HTM	High Temperature Cut Off Module	
-ETLM	External Temperature Limit Module	
-A1925E4	RS-485 to USB Converter (required for data logging)	
-A1925E4	RS-485 to USB Converter, isolated	
-A3504HC	SpecView Software Package	
K Electrical Supply for Heater and/or Controller		

-No Symbol	None	
-115	115V Heater	
-230	230V Heater	

L Certifications		
-ASME	ASME Documentation	
-CE	CE Documentation	
-P	Parr Certification	

Model 4625-1200mL Vessel w/ Bench Top Heater (BTH)

Series Number:

Type: High Pressure High Temperature

Stand: Bench Top

Mounting Style: Moveable or Fixed Head

Vessel Sizes, mL: **250-1000**

High Pressure MAWP Rating, psi (bar): **4200 (290)**

Standard Maximum Operating Temperature, °C:

500 w/ Grafoil Gasket


High Temperature MaximumOperating Temperature, °C:

600 w/ Grafoil Gasket

Series 4650 High Pressure / High Temperature Vessel Systems

4650 Series High Pressure and High Temperature Vessels are available in volumes of 250, 500, and 1000 mL. These vessels can obtain temperatures up to 600 °C and a maximum pressure of 6000 psi (410 bar). The maximum pressure at maximum temperature is 4200 psi (290 bar) at 600 °C. They come standard with a Split Ring Closure with 8 Compression bolts, and a VGR head design that includes Valve, Gage, Rupture Disc, and Thermowell. An additional valve with dip tube can be added.

The maximum allowable head openings are seven but this may be limited by the types of fittings required. These vessels are available as either a moveable or fixed head design.



4650-1000mL Cross Section

Shaded bar indicates sp	ecifications that change	within series.	
Model Number	4650		
Approximate Volume, mL	250	500	1000
Maximum Pressure (MAWP)	5000 psi (345 bar) @ 500 °C 6000 psi (413 bar) @ 350 °C 4200 psi (289 bar) @ 600 °C		
Maximum Temperature			
with PTFE Flat Gasket		350 °C	
with Flexible Graphite Flat Gasket		500-600 °C	
Vessel Details			
Mounting Style		Moveable or Fixed Head	
Closure	Spli	t-Ring (8 Compression Bo	olts)
Head Opening Options	Ports for valve, gage, rupture disc, and thermocouple (VGR) Custom port option (Contact our technical sales department to discuss options)		
Maximum Head Openings	7 (Dependent	on opening size, and req	uired fittings)
Valve Connections	1/4" NPT Female		
Pressure Gage, Size	4.5 inches		
Range	0-5000 psi (345 psi) or 0-7500 psi (517 bar)		
Temperature Measurement	Thermowell		
Cooling Coil (optional)	Cooling Coil (optional) Single Loop		
Stand and Heater (optional)			
Stand Type	Bench Top Heater with base plate (moveable head only) and Bench Top Stand with Heater (fixed head)		
Heater Style	Ceramic Fiber Heater		
Heater Power, Watts		1500	
Maximum Load, amps, 115V / 230V		15 / 9	
Weights & Dimensions			
Cylinder I.D. x Depth, inches	2.5 x 3.25	2.5 x 6.63	2.5 x 13.13
Vessel Assembly w/ Moveable Head Weight, pounds**	38	45	57
Vessel Assembly w/ Fixed Head Weight, pounds**	33	40	52
Cylinder Weight, pounds	8.9	15	26
** Vessel Assembly weight varies bas Weights shown are based on vessels with valve			
Other options available. See Ordering Gui	de, visit <u>www.parrinst.co</u>	om, or call for more inform	nation.

I Controller -No Symbol

-4838

-4848

None

module.

Series 4650 Ordering Guide

PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable

not included). For use with one additional display

PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not

included). For use with up to three additional display

An example order number for a vessel in this series is: 4650-250mL-FG-SS-VGR-VD-5000

A composite identification number to be used when ordering a 4650 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Volume				
-250mL 250 mL, 2.5" ID				
-500mL	500 mL, 2.5" ID			
-1000mL	1000 mL, 2.5" ID			
B Head Mou	unting Style			
-No Symbol	Moveable Head (Standard configuration)			
-FH	Fixed Head			
C Gasket / I	Naximum Temperature			
-T	PTFE Flat Gasket, 350 °C			
-FG	Flexible Graphite Flat Gasket, 350 °C			
D Material o	of Construction			
-SS	T316 Stainless Steel			
-C20	Alloy 20*			
-M0	Alloy 400*			
-IN600	Alloy 600			
-IN625	Alloy 625			
-HB	Alloy B-2 / B-3*			
-HC	· · · · ·			
-HC Alloy C-276 -NI Nickel 200 -TI2 Titanium Grade 2* -TI4 Titanium Grade 4* -ZR702 Zirconium Grade 702*				
		-ZR705	Zirconium Grade 702*	
		-ZK/05 Zirconium Grade 705* *High temperature option not available.		
			<u>4</u> for complete list of available alloys.	
E Head Con	figuration			
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple and thermowell			
-NS	Non-standard, custom head opening			
F Optional	Fittings & Custom Options (List all needed)			
-VD	Single Valve on head w/ dip tube (VGR Required)			
-DVD	Double Valve Assembly on head w/ dip tube (VGR Required)			
-CL	Cooling Loop			
-BDV	Bottom Drain Valve			
-CAD Internal Catalyst Addition Device				
-XCAD	External Catalyst Addition Device			
000				

	modules.			
See <u>Chapter 6</u>	er 6 for complete list of Controller options.			
J Control	ontroller Options (Check all needed)			
-PDM	Pressure Display Module			
-HTM	High Temperature Cut Off Module			
-ETLM	External Temperature Limit Module			
-SVM	Solenoid Valve Module (for cooling control)			
-A1925E4	RS-485 to USB Converter (required for data logging)			
-A1925E6	RS-485 to USB Converter, isolated, 30-ft			
-A2208E	RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)			
-A3504HC	SpecView Software Package			
K Electric	al Supply for Heater and/or Controller			
-No Symbol	None			
-115	115V Heater			
-230	230V Heater			
L Certific	No Certification			
-No Symbol -ASME	ASME Documentation			
-ASIVIE -CE	CE Documentation			
-0E -P	Parr Certification			
-r				

-RTC	Reflux/Take-Off Condenser Welded Jacket	
-WJ		
G Pressure Gage		
-No Symbol	None	
-10000	10000 psi / 690 bar	
-7500	7500 psi / 517 bar	
-5000	5000 psi / 345 bar	
-3000	3000 psi / 207 bar	
-2000	2000 psi / 138 bar	
H Stand and Heater Options		
-No Symbol	None	
DTU		

Solids Charging Port (Ball Valve)

Reflux Condenser

-ino Symbol	None
-BTH	Bench Top heater with base plate
-BTS	Bench Top Stand with Heater

Model 4650-250mL Vessel with VGR, Thermocouple, and Single Inlet Valve (-VD)

Model 4650-500mL Removeable Head Vessel with Bench Top Heater and Base Plate

-SCP

-RC



Type: General Purpose

Stand: Floor Stand

Mounting Style: Moveable or Fixed Head

Vessel Sizes, Gallons (Liters) 1 (3.75) and 2 (7.99)

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket



4660 Cross Section with Flat Gasket Seal



46600 Cross Section with O-ring Seal

Series 4660 Pressure Vessel Systems, 1 & 2 Gallon

The 4660 Vessels extend the range of Parr Pressure Vessels to 1 and 2 gallon (3.75 and 7.99 liter) sizes. These vessels are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

Both fixed head and moveable head designs are available. The moveable head vessels are designed so that they can be opened or closed conveniently without removing the cylinder from the heater and without auxiliary handling equipment. The split-ring cover clamp sections can be moved into place from the sides, and the compression bolts can be tightened with the vessel in place in its heater. The fixed head versions allow the head (and the connections made to it) to remain in place while the cylinder and heater are lowered with the aid of the included pneumatic lift system.

An optional bottom drain valve may be added for convenient product recovery. As with the smaller floor stand models, these larger, self-contained systems can be equipped with a variety of attachments, such as condensers, solids charging ports, bottom drains, special heaters, jacketed vessels and automatic valves and regulators. Because of the higher wattage heaters for these reactors, all models in this series require a 230 volt power supply.

Series 4660 Pressure Vessel Specifications				
Shaded bar indicate	s specifications that change within se	ries.		
Model Number	4660			
Approximate Volume, Gallon (Liter)	1 (3.75)	2 (7.99)		
Maximum Pressure (MAWP)	1900 ps	ii (131 bar)		
Maximum Temperature				
with FKM 0-Ring	22	25 °C		
with FFKM 0-Ring	30	0° 00		
with PTFE Flat Gasket	3	50 °C		
with Flexible Graphite Flat Gasket	3	50 °C		
Vessel Details				
Mounting Style	Moveable	or Fixed Head		
Closure	Split-Ring with 10 Compression Bo	olts, no Compression Bolts for O-ring)		
Head Opening Options	One head opening (1/4" NPT or rupture disc assembly port) Ports for valve, gage, rupture disc, and thermowell with thermocouple (VGR) Custom port option (Contact our technical sales department to discuss options)			
Maximum Head Openings	10 (Dependent on opening size, and required fittings)			
Valve Connections	1/4" NPT Female			
Pressure Gage, Size	4.5 inches			
Range	0-2000 psi (138 bar)			
Temperature Measurement Thermowell		mowell		
Cooling Coil (optional)	Ser	pentine		
Bottom Drain Valve	3/8	" NPT		
Stand and Heater (optional)				
Stand Type	Floor Cart with heater (moveable head only) Floor Stand with Pneumatic Lift and heater (fixed head)			
Heater Style	C	alrod		
Heater Power, Watts	2250	2750		
Maximum Load, amps, 230V	10	12		
Weights & Dimensions				
Cylinder I.D. x Depth, inches	6.0 x 8.5	6.0 x 17.1		
Vessel Assembly w/ Moveable Head Weight, pounds*	84	103		
Vessel Assembly w/ Fixed Head Weight, pounds*	98	117		
Cylinder Weight, pounds	33	52		
* Vessel Assembly weight varies based on head configuration choice and fitting options. Weights shown are based on vessels with valve, gage, rupture disc, thermocouple, and double valve assembly.				
Other options available. See Ordering	Guide, visit <u>www.parrinst.com</u> , or cal	l for more information.		

Series 4660 Ordering Guide

R

An example order number for a vessel in this series is: 4660-1G-T-SS-VGR-2000

A composite identification number to be used when ordering a 4660 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Volume				
-1G 1 Gallon, 6" ID				
-2G 2 Gallon, 6" ID				
B Head Mounting Style				
-No Symbol	Moveable Head (Standard configuration)			
-FH	Fixed Head			
\sim	aximum Temperature			
-T	PTFE Flat Gasket, 350 °C			
-FG	Flexible Graphite Flat Gasket, 350 °C			
	Construction			
-SS	T316 Stainless Steel			
-C20	Alloy 20			
-M0	Alloy 400			
-IN600	Alloy 600			
-IN625	Alloy 625			
-HB	Alloy B-2 / B-3			
-HC	Alloy C-276			
-NI	Nickel 200			
-TI2	Titanium Grade 2			
-TI4	Titanium Grade 4			
-ZR702	Zirconium Grade 702			
-ZR705	Zirconium Grade 705			
See <u>page 10</u> or <u>24</u>	for complete list of available alloys.			
E Head Confi	guration			
-1/4	One head opening, 1/4" NPT port, plugged			
-RD	Rupture disc port with rupture disc assembly & disc			
-VGR	Valve on gage adapter, pressure gage, rupture disc assem- bly, & thermocouple and thermowell			
-NS	Non-standard, custom head opening			
E Optional Fi	ttings & Custom Options (List all needed)			
-VD	Single Valve on head w/ dip tube (VGR Required)			
-DVD	Double Valve Assembly on head w/ dip tube (VGR Required)			
-CL	Cooling Loop			
-BDV	Bottom Drain Valve			
-CAD	Internal Catalyst Addition Device			
-XCAD	External Catalyst Addition Device			
-SCP	Solids Charging Port (Ball Valve)			
-RC	Reflux Condenser			
-RTC	Reflux/Take-Off Condenser			
-WJ	Welded Jacket			
G Pressure G				
-No Symbol	None			
-2000	2000 psi / 138 bar			
-1000	1000 psi / 69 bar			
-600	600 psi / 40 bar			
-200	200 psi / 14 bar			
-100	100 psi / 7 bar			
H Stand and	Heater Options			
-No Symbol	None			
-FC	Floor Cart with Heater			
-FSP	Floor Stand with Pneumatic Lift and heater			

I Contro			
-No Symbol	None PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module. PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.		
-4838			
-4848			
See <u>Chapter (</u>	§ for complete list of Controller options.		
J Contro	ller Options (Check all needed)		
-PDM	Pressure Display Module		
-HTM	High Temperature Cut Off Module		
-ETLM	External Temperature Limit Module		
-SVM	Solenoid Valve Module (for cooling control)		
-A1925E4	RS-485 to USB Converter (required for data logging)		
-A1925E6 RS-485 to USB Converter, isolated, 30-ft -A2208E RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6) -A3504HC SpecView Software Package			
		K Electric	cal Supply for Heater and/or Controller
		-No Symbol	None
-230	230V Heater		
L Certific	ations		
-No Symbol	No Certification		
-ASME	ASME Documentation		
	CE De sum entetien		
-CE	CE Documentation		



Model 4660-2G Vessel, Fixed Head Floor Stand, with **Pneumatic Lift and Heater**





^{Type:} High Pressure High Temperature

Stand: Cart or Floor Stand

Mounting Style: Moveable or Fixed Head

Vessel Sizes, Liters: 1 and 1.8

High Pressure MAWP Rating, psi (bar): **6000 (410)**


Standard Maximum Operating Temperature, °C: 600 w/ Flat FG Gasket



4680-1800mL Fixed Head Vessel shown with Split Rings open and vessel down using a pneumatic lift.

Series 4680 High Pressure / High Temperature Vessel Systems, 1 & 1.8 Liters

The Series 4680 High Pressure/High Temperature Vessels are available with either a 1.0 L or 1.8 L capacity. These vessels can obtain temperatures up to 600 °C and a maximum pressure of 6000 psi (410 bar). The maximum pressure at maximum temperature is 4200 psi (290 bar) at 600 °C. A Split Ring closure with 12 Compression Bolts secures the head which allows a maximum of 8 openings. The Series 4680 Vessels can be configured with either a fixed head on a floor stand or as a moveable head vessel on a moveable cart. The fixed head versions allow the head (and the connections made to it) to remain in place while the cylinder and heater are lowered with the aid of the included pneumatic lift system.



4680 Cross Section

Series 4680 Pressure Vessel Specifications				
Shaded bar indicates specifications that change within series.				
Model Number	4680			
Approximate Volume, mL	1000	1800		
Maximum Pressure (MAWP)	5000 psi (345 bar) @ 500 °C 6000 psi (413 bar) @ 350 °C 4200 psi (289 bar) @ 600 °C			
Maximum Temperature				
with PTFE Flat Gasket	35	50 °C		
with Flexible Graphite Flat Gasket	500-	-600 °C		
Vessel Details				
Mounting Style	Moveable	or Fixed Head		
Closure	Split-Ring with 12	2 Compression Bolts		
Head Opening Options	One head opening ("B" socket, or rupture disc assembly port) Ports for valve, gage, rupture disc, and thermowell with thermocouple (VGR) Custom port option (Contact our technical sales department to discuss options)			
Maximum Head Openings	8 (Dependent on opening	g size, and required fittings)		
Valve Connections	1/4" NPT Female			
Pressure Gage, Size	4.5 inches			
Range	0-5000 psi (345 bar) or 0-7500 psi (517 bar)			
Temperature Measurement	Ther	mowell		
Cooling Coil (optional)	Single Loop			
Bottom Drain Valve	1" NPS			
Stand and Heater (optional)				
Stand Type		er (moveable head only) ic Lift and heater (fixed head)		
Heater Style	Ceramic Fiber			
Heater Power, Watts	2300	2500		
Maximum Load, amps, 230V	10	11		
Weights & Dimensions				
Cylinder I.D. x Depth, inches	3.75 x 6.12	3.75 x 10.5		
Vessel Assembly w/ Moveable Head Weight, pounds*	74	94		
Vessel Assembly w/ Fixed Head Weight, pounds*	88	104		
Cylinder Weight, pounds	36	52		
*Vessel weight is based on a vessel with VGR head fittings. Vessel weight will vary based on head configuration choice and fitting options.				
Other options available. See Ordering	Guide, visit <u>www.parrinst.com</u> , or cal	l for more information.		

Series 4680 Ordering Guide

8

An example order number for a vessel in this series is: 4680-1000mL-FG-SS-VGR-5000


A composite identification number to be used when ordering a 4680 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Vol	uma			
A Vessel Volume -1000mL 1000 mL, 3.75" ID				
-1800mL 1800 mL, 3.75 ID				
B Head Mounting Style				
-No Symbol Moveable Head (Standard configuration)				
-FH	Fixed Head			
C Gasket / Maximum Temperature				
-T PTFE Flat Gasket, 350 °C				
-FG	Flexible Graphite Flat Gasket, 350 °C			
• Material of	f Construction T316 Stainless Steel			
-C20	Alloy 20			
-M0	Alloy 400			
-IN600	Alloy 600			
-IN625	Alloy 625			
-HB	Alloy B-2 / B-3			
-HC	Alloy C-276			
-NI	Nickel 200			
-TI2	Titanium Grade 2			
-TI4	Titanium Grade 4			
-ZR702	Zirconium Grade 702			
-ZR705	Zirconium Grade 705			
See <u>page 10</u> or <u>24</u>	for complete list of available alloys.			
E Head Conf	iguration			
-В	One head opening, "B" socket			
-RD	Rupture disc port with rupture disc assembly & disc			
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple and thermowell			
-NS	Non-standard, custom head opening			
F Optional F	ittings & Custom Options (List all needed)			
-VD	Single Valve on head w/ dip tube (VGR Required)			
-DVD	Double Valve Assembly on head w/ dip tube (VGR Required)			
-CL	Cooling Loop			
-BDV	Bottom Drain Valve			
-4310B Gage Block Assembly (for head with "B" socket open				
-CAD	Internal Catalyst Addition Device			
-XCAD	External Catalyst Addition Device			
-SCP	Solids Charging Port (Ball Valve)			
-RC	Reflux Condenser			
-RTC	Reflux/Take-Off Condenser			
-WJ	Welded Jacket			
G Pressure G				
-No Symbol	None			
-10000	10000 psi / 690 bar			
-7500	7500 psi / 517 bar			
-5000	5000 psi / 345 bar			
-3000	3000 psi / 207 bar			
-2000	-2000 2000 psi / 138 bar			
H Stand and Heater Options				
-No Symbol	None			
1				

-ivo Symbol	None
-FC	Floor Cart with heater (moveable head only)
-FSP	Floor Stand w/ Pneumatic Lift and heater (fixed head)

Controller			
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module.		
-4848	PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.		
See <u>Chapter (</u>	§ for complete list of Controller options.		
J Control	ller Options (Check all needed)		
-PDM	Pressure Display Module		
-HTM	High Temperature Cut Off Module		
-ETLM	External Temperature Limit Module Solenoid Valve Module (for cooling control) RS-485 to USB Converter (required for data logging) RS-485 to USB Converter, isolated, 30-ft RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6) SpecView Software Package for 4838/4848 Controller		
-SVM			
-A1925E4			
-A1925E6			
-A2208E			
-A3504HC			
K Electric	al Supply for Heater and/or Controller		
-No Symbol	None		
-230	230V Heater		
L Certific	ations		
-No Symbol	No Certification		
-ASME	ASME Documentation		
-CE	CE Documentation		
-P	Parr Certification		







4680-1800mL Moveable Head High Pressure Vessel in Floor Cart with Heater

Series Number:

Type: **High Pressure High Temperature**

Stand: Cart or Floor Stand

Mounting Style: Moveable or Fixed Head


Vessel Sizes, Gallons (Liters): 1 (3.75) and 1.5 (5.8)

Standard Pressure MAWP Rating, psi (bar): **3000 (200)**

Standard Maximum Operating Temperature, °C: 600 w/ Flat FG Gasket Series 4670 HP-HT Vessel Systems, 1 and 1.5 Gallon

The Series 4670 High Pressure/High Temperature Vessels are available in either a 1.0 or 1.5 Gallon capacity. These vessels can obtain temperatures up to 600 °C and a maximum pressure of 3000 psi (207 bar). The maximum pressure at maximum temperature is 2200 psi (151 bar) at 600 °C. A Split Ring closure

with 16 Compression Bolts secures the head which allows a maximum of 12 openings. The Series 4670 Vessels can be configured with either a fixed head on a floor stand or as a moveable head vessel on a moveable cart. The fixed head versions allow the head (and the connections made to it) to remain in place while the cylinder and heater are lowered with the aid of the included pneumatic lift system.



4670 Cross Section

	s specifications that change within se	
Model Number	Vodel Number 4670	
Approximate Volume, Gallon (Liter)	1 (3.75)	1.5 (5.8)
Maximum Pressure (MAWP)	3000 psi (207 bar, 200 bar for CE orders) @ 500 °C 2200 psi (151 bar) @ 600 °C	
Maximum Temperature		
with PTFE Flat Gasket	3	50 °C
with Flexible Graphite Flat Gasket	500	-600 °C
Vessel Details		
Mounting Style	Moveable	or Fixed Head
Closure	Split-Ring with 1	6 Compression Bolts
Head Opening Options	One head opening (rupture disc assembly port) Ports for valve, gage, rupture disc, and thermowell with thermocouple (VGR) Custom port option (Contact our technical sales department to discuss options)	
Maximum Head Openings 12 (Dependent on opening size, and required fittings)		ig size, and required fittings)
Valve Connections	1/4" N	PT Female
Pressure Gage, Size	4.5 inches	
Range	0-3000 psi (207 bar)	
Temperature Measurement	Thermowell	
Cooling Coil (optional) Serpentine		pentine
Bottom Drain Valve (BDV) (optional)	1" NPS (Floor Stand or Moveable Cart Only)	
Stand and Heater (optional)		
Stand Type*	Stand Type* Floor Cart with heater (moveable head only) Floor Stand with Pneumatic Lift and heater (fixed head)	
Heater Style	Cerai	nic Fiber
Heater Power, Watts		
Maximum Load, amps, 230V		12
*Heater watts and amps vary between differer	nt heater styles, contact Parr Technica	al Support for more information.
Weights & Dimensions		
Cylinder I.D. x Depth, inches	5.5 x 9.75	5.5 x 15.0
Vessel Assembly w/ Moveable Head Weight, pounds*	113	131
Vessel Assembly w/ Fixed Head Weight, pounds*	124	143
Cylinder Weight, pounds	3653	72
*Vessel weight is Vessel weight will vary base	based on a vessel with VGR head fittir ad on head configuration choice and fi	igs. tting options.

Series 4670 Ordering Guide

8

An example order number for a vessel in this series is: 4670-1G-FG-SS-VGR-5000

A composite identification number to be used when ordering a 4670 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

-P

A Vessel Volu	ume			
-1G	1 Gallon, 5.5" ID			
-1.5G	1.5 Gallon, 5.5" ID			
\sim				
B Head Moun				
-No Symbol	Moveable Head (Standard configuration)			
-FH	Fixed Head			
C Gasket / Maximum Temperature				
-T	PTFE Flat Gasket, 350 °C			
-FG	Flexible Graphite Flat Gasket, 500 or 600 °C			
D Material of Construction				
-SS	T316 Stainless Steel			
-C20	Alloy 20*			
-M0	Alloy 400*			
-IN600	Alloy 600			
-IN625	Alloy 625			
-HB	Alloy B-2 / B-3*			
-HC	Alloy C-276			
-NI	Nickel 200			
-TI2	Titanium Grade 2*			
-TI4	Titanium Grade 4*			
-ZR702	Zirconium Grade 702*			
-ZR705	Zirconium Grade 705*			
	option not available. for complete list of available alloys.			
E Head Confi				
-RD	Rupture disc port with rupture disc assembly & disc			
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple and thermowell			
-NS Non-standard, custom head opening				
F Optional Fittings & Custom Options (List all needed)				
F Optional Fi				
F Optional Fi				
	ttings & Custom Options (List all needed)			
-VD	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required)			
-VD -DVD	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required)			
-VD -DVD -SC	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil			
-VD -DVD -SC -BDV	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS)			
-VD -DVD -SC -BDV -CAD	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device			
-VD -DVD -SC -BDV -CAD -XCAD	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device			
-VD -DVD -SC -BDV -CAD -XCAD -SCP	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve)			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RC -RTC	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser			
-VD -DVD -SC -BDV -CAD -CAD -XCAD -SCP -RC -RTC G Pressure G	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser age			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RTC G Pressure G -No Symbol	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux/Take-Off Condenser age None			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RTC G Pressure G -No Symbol -5000	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser age None 5000 psi / 345 bar			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RC -RTC G Pressure G -No Symbol -5000 -3000 -2000	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser age None 5000 psi / 345 bar 3000 psi / 207 bar 2000 psi / 138 bar			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RC -RTC G Pressure G -No Symbol -5000 -3000 -2000 H Stand and	the second seco			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RC -RTC G Pressure G -No Symbol -5000 -3000 -2000	ttings & Custom Options (List all needed) Single Valve on head w/ dip tube (VGR Required) Double Valve Assembly on head w/ dip tube (VGR Required) Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser age None 5000 psi / 345 bar 3000 psi / 207 bar 2000 psi / 138 bar Heater Options None			
-VD -DVD -SC -BDV -CAD -XCAD -SCP -RC -RC -RTC G Pressure G -No Symbol -5000 -3000 -2000 H Stand and -No Symbol	the second seco			

Controller			
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with one additional display module.		
-4848	PID Control, Ramp & Soak Programming, and Data logging with Software. (RS-485 to USB cable not included). For use with up to three additional display modules.		
See <u>Chapter 6</u> for complete list of Controller options.			
J Controller Options (Check all needed)			
-PDM	Pressure Display Module		
-HTM	High Temperature Cut Off Module		
-ETLM	External Temperature Limit Module		
-SVM	Solenoid Valve Module (for cooling control)		
A 100FF4	DC 405 to UCD Conventor (no pulse of few data la police)		

-SVM	Solenoid Valve Module (for cooling control)		
-A1925E4	RS-485 to USB Converter (required for data logging)		
-A1925E6	RS-485 to USB Converter, isolated, 30-ft		
-A2208E	RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)		
-A3504HC	SpecView Software Package for 4838/4848 Controller		
K Electrical Supply for Hester and/or Controller			
K Electric	al Supply for Heater and/or Controller		
	al Supply for Heater and/or Controller		
K Electric -No Symbol	al Supply for Heater and/or Controller None		
\sim			
-No Symbol -230	None 230V		
-No Symbol	None 230V		
-No Symbol -230	None 230V		
-No Symbol -230	None 230V ations		

Parr Certification



Model 4670-1.5G Moveable Head Vessel on a Floor Cart with Heater and 4848 Controller



Type: General Purpose

Stand: Floor Stand

Mounting Style: Moveable or Fixed Head

Vessel Sizes, Gallons (Liters): **2.6 (10) and 5 (18.75)**

Standard Pressure MAWP Rating, psi (bar): **1900 (131)**

Standard Maximum Operating Temperature, °C: 225 w/ FKM 0-ring 300 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

Series 4675 General Purpose Vessel Systems, 2.6 & 5 Gallon


The Series 4675 General Purpose Vessels are the largest vessels

offered by Parr Instrument Company. Vessel styles are offered in a moveable head or fixed head design. These vessels are available with an FKM O-ring seal for operating temperatures to 225 °C, and FFKM O-ring for temperatures to 300 °C, or with a flat, PTFE gasket for operating temperatures up to 350 °C maximum.

In the moveable head design, the vessel is held in a support system which minimizes the physical effort required to handle these heavy components. The hoist is attached to a support column which provides a convenient means for lifting the head and cylinder out of the stand. Vessels equipped with a bottom drain valve will probably remain in the heater most of the time, but can be lifted out when necessary.

The fixed head support stand features hinged split-rings that swing to either side allowing the head to remain fixed to the stand while a pneumatic lift allows the cylinder to be raised and lowered. When lowered, the cylinder can be slid forward for cleaning and servicing.

These vessels are generally used for pilot plant or for custom chemical production purposes, usually with a variety of attachments added to the basic units.



4675 Cross Section



Model 4675-10L General Purpose Vessel

Model 4675-18.5L General Purpose Vessel with heater, and 4848M Temperature Controller



4675

8

	Shaded bar indicates specifications that change within	series.		
Model Number	46	75		
Approximate Volume, Gallon (Liter)	5 (18.75) 2.6 (10)			
Maximum Pressure (MAWP) 1900 psi (131 bar)		(131 bar)		
Maximum Temperature				
with FKM 0-Ring	225 °C			
with FFKM 0-Ring	300 °C			
with PTFE Flat Gasket	350	О° (
with Flexible Graphite Flat Gasket	350	3° (
Vessel Details				
Mounting Style	Moveable o	r Fixed Head		
Closure	Split-Ring with 12 Compression Bolts	Split-Ring with 12 Compression Bolts		
Head Opening Options	Ports for valve, gage, rupture disc, and	IPT, or rupture disc assembly port) d thermowell with thermocouple (VGR) cal sales department to discuss options)		
Maximum Head Openings	12 (Dependent on opening	size, and required fittings)		
Valve Connections	1/4" NP	T Female		
Pressure Gage, Size	4.5 ir	iches		
Range	0-2000 psi (138 bar)			
Temperature Measurement	Thern	Thermowell		
Cooling Coil (optional)	Serpe	entine		
Bottom Drain Valve (BDV) (optional)	1" N	NPT		
Stand and Heater (optional)				
Stand Type		nual hoist (moveable head) matic Lift (fixed head)		
Heater Options	Band Heater, 3-zone Ceramic, 3-zone (moveable head) Flexible Mantle Heater (low temperatures <200 °C)	Band Heater, 3-zone Flexible Mantle Heater (low temperatures <200 °C)		
Heater Power Requirements	Typical power requirements for Parr's large, electrically heated pressure vessels are 40 Amp single phase or 3-Phase power sources. Users are advised to have a qualified electrician determine and install an appropriate mains power supply for the large pressure vessel system. Large pressure vessel systems with lower electrical power requirements, such as low temperature applications are available for use with typical 20 amp 230 volt sources. Contact Parr Technical Sales staff for assistance with electrical specifications.			
Weights & Dimensions				
Cylinder I.D. x Depth, inches	9.5 x 16.25 7.75 x 12.2			
Moveable Vessel Assembly Weight, pounds*	335	194		
Fixed Vessel Assembly Weight, pounds*	336	195		
Cylinder Weight, pounds	157	97		
Vessel	*Vessel weight is based on a vessel with VGR head fir weight will vary based on head configuration choice and	ttings. I fitting options.		
Other options available. See Ordering Guide, visit <u>www.parrinst.com</u> , or call for more information.				

Series 4675 Ordering Guide

An example order number for a vessel in this series is: 4675-5G-FG-SS-VGR-5000

A composite identification number to be used when ordering a 4675 Series Pressure Vessel System can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Vessel Volume				
-2.6G 2.6 Gallon, 7.75" ID				
-5G	5 Gallon, 9.5" ID			
B Head Mounting Style				
-No Symbol	Moveable Head (Standard configuration)			
-FH	Fixed Head			
C Gasket / Maximum Temperature				
-T PTFE Flat Gasket, 350 °C				
-0V	FKM 0-ring, 225 °C FFKM 0-ring, 300 °C			
-0K				
-FG Flexible Graphite Flat Gasket, 350 °C				
D Material of Construction				
-SS T316 Stainless Steel				
-C20	Alloy 20			
-M0	Alloy 400			
-IN600	Alloy 600			
-IN625	Alloy 625			
-HB	Alloy B-2 / B-3			
-HC	Alloy C-276			
-NI	Nickel 200			
-TI2	Titanium Grade 2			
-TI4	Titanium Grade 4			
-ZR702	Zirconium Grade 702			
-ZR705	Zirconium Grade 705			
See <u>page 10</u> or <u>24</u> for complete list of available alloys.				
E Head Cont				
-1/4	One head opening, 1/4" NPT port, plugged			
-3/8	One head opening, 3/8" NPT port, plugged			
-RD Rupture disc port with rupture disc assembly & disc				
-VGR	Valve on gage adapter, pressure gage, rupture disc assembly, & thermocouple and thermowell			
-NS	Non-standard, custom head opening			
F Optional F	ittings & Custom Options (List all needed)			
-VD	Single Valve on head w/ dip tube (VGR Required)			
-DVD	Double Valve Assembly on head w/ dip tube			
	(VGR Required)			
-SC	(VGR Required) Serpentine Cooling Coil			
-SC	Serpentine Cooling Coil			
-SC -BDV	Serpentine Cooling Coil Bottom Drain Valve (1" NPS)			
-SC -BDV -CAD	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device			
-SC -BDV -CAD -XCAD	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device			
-SC -BDV -CAD -XCAD -SCP	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve)			
-SC -BDV -CAD -XCAD -SCP -RC -RTC	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser			
-SC -BDV -CAD -XCAD -SCP -RC -RTC	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser			
-SC -BDV -CAD -XCAD -SCP -RC -RTC G Pressure (Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser Gage (VGR Required)			
-SC -BDV -CAD -XCAD -SCP -RC -RTC G Pressure C -No Symbol	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser Gage (VGR Required) None			
-SC -BDV -CAD -XCAD -SCP -RC -RTC G Pressure -No Symbol -2000 -1000 -600	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser Gage (VGR Required) None 2000 psi / 138 bar 1000 psi / 69 bar 600 psi / 40 bar			
-SC -BDV -CAD -XCAD -SCP -RC -RTC G Pressure (-No Symbol -2000 -1000 -600 -200	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser Sage (VGR Required) None 2000 psi / 138 bar 1000 psi / 69 bar 600 psi / 40 bar 200 psi / 14 bar			
-SC -BDV -CAD -XCAD -SCP -RC -RTC C Pressure C -No Symbol -2000 -1000 -600	Serpentine Cooling Coil Bottom Drain Valve (1" NPS) Internal Catalyst Addition Device External Catalyst Addition Device Solids Charging Port (Ball Valve) Reflux Condenser Reflux/Take-Off Condenser Gage (VGR Required) None 2000 psi / 138 bar 1000 psi / 69 bar 600 psi / 40 bar			

H Stand - -No Symbol	None		
-NU Symbol			
-FS	Floor Stand, Table Support with Manual Lift (Moveable Head)		
-FSEH	Floor Stand, Table Support with Electric Hoist (115 VAC) (Moveable Head)		
-FSP	SP Floor Stand with Pneumatic Lift (Fixed Head)		
Heater Options (require stand)			
-No Symbol	None		
-C3	Ceramic, 3-zone (moveable head, 5 Gallon only)		
-B3	Band Heaters, 3-zone		
-WJ	Welded Jacket		
-FM	Flexible Mantle Heater, 1-zone (200 °C Max. Temp.)		
J Contro	ller		
-No Symbol	None		
-4838	PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included) For use with one additional display module. (Available only for 10 L/18.75 L with Flexible Mantle Heate or 10L with Band Heaters		
-4848M	With Slave Box, PID Control, Ramp & Soak Programming, and Data logging with Software (RS-485 to USB cable not included). For use with up to three additional display mod- ules. (Recommended for 10 L with Band Heaters. Required for 18.75 L with Ceramic or Band Heaters		
See <u>Chapter 6</u> for complete list of Controller options.			
K Contro	ller Options (Check all needed)		
-PDM Pressure Display Module			
-HTM	High Temperature Cut Off Module		
-ETLM	External Temperature Limit Module		
-SVM	Solenoid Valve Module (for cooling control)		
-A1925E4	RS-485 to USB Converter (required for data logging)		
-A1925E6 RS-485 to USB Converter, isolated, 30-ft			
-A2208E RS-485 Daisy Chain for Multiple Controllers (Must be used with A1925E6)			
-A3504HC	SpecView Software Package for 4838/4848 Controller		
L Electrical Supply for Heater and/or Controller			
	None		
-No Symbol	None		
\sim	230 VAC, Single Phase (Typically 40A supply required)		
-No Symbol			
-No Symbol -230	230 VAC, Single Phase (Typically 40A supply required) 400-415V 3-Phase "Y"		
-No Symbol -230 -400	230 VAC, Single Phase (Typically 40A supply required) 400-415V 3-Phase "Y"		
-No Symbol -230 -400 M Certific	230 VAC, Single Phase (Typically 40A supply required) 400-415V 3-Phase "Y" cations		
-No Symbol -230 -400 M Certific -No Symbol	230 VAC, Single Phase (Typically 40A supply required) 400-415V 3-Phase "Y" cations No Certification		

4310A Gage

Block Assembly

4310B Gage

Block Assembly

Gage Block Assemblies

Parr gage block assemblies combine the

functions of an Inlet / Outlet Valve, a Pressure Gage, and a Safety Rupture Disc into a compact assembly which can be attached to the head of any small vessel with limited head space using a connecting tube. There is a threaded socket in the block for a gas connection with a pressure hose or tubing using a type "A" coned pressure fitting. Coned fittings are also used on the thick-walled tube which connects the block to the pressure vessel. The valve in this assembly controls the gas flow into the vessel and the gage shows the pressure in the vessel when the valve is closed. A rupture disc with a burst rating matched to the gage is installed in the block.


For small vessels with a 1-inch I.D. that can only accommodate a singe head port, we offer a version of the 4310A with a 1/16" O.D. thermocouple.

Two Sizes

These gage block assemblies are made in two styles which differ only in the size of the pressure gage and the size of the tube connection to the pressure vessel. The smaller 4310A Gage Block Assembly is usually furnished with a 3-1/2 inch diameter pressure gage and includes a 3-inch long connecting tube with type "A" coned fitting. This assembly is normally installed on pressure vessels up to 600 mL.

The 4310B Gage Block Assembly has a 4-1/2 inch diameter gage and type "B" coned fitting on a 3-inch long connecting tube. This assembly is normally installed on vessels of 1000 mL and larger.

The gage block, connecting tube and pressure gage on all of these assemblies are normally made of type 316 Stainless Steel, but they can also be made of Alloy 400 on special order. Longer tubes can be furnished on special order.



Available Gages


Any of the gages shown in this table can be installed on these gage block assemblies.

Gages for Parr Pressure Vessels			
Pressure Range		4-1/2" Diameter	3-1/2" Diameter
psi	bar	Gage Number	Gage Number
0-100	0-7	56HCPA	593HCP1AD
0-200	0-14	56HCPB	593HCP2AD
0-600	0-40	56HCPC	593HCP6AD
0-1000	0-69	56HCPD	593HCPD
0-2000	0-138	56HCPF	593HCPF
0-3000	0-207	56HCPG	593HCPG
0-5000	0-345	56HCPH	593HCP50AD
0-7500	0-517	56HCP75A	N/A
0-10000	0-690	56HCPK	N/A

When ordering, specify the catalog number, pressure gage diameter, range, and if a thermocouple is required.



Gage Blocks, pressure hoses and other parts that are frequently removed from a pressure vessel are attached with a coned, socket type connector. The male segment of a coned fitting consists of a sleeve with a left-hand thread screwed onto a thick-walled tube, plus a compression nut to complete the assembly. When screwed into a mating socket these parts produce a rigid joint which will remain tight over a wide temperature and pressure range, yet the joint can be made and broken repeatedly without destroying the sealing faces. No gasket or sealing tape is required. These fittings have the added advantage of behaving like union fitting,




allowing the gage block, connecting tube or hose to remain stationary while the nut is tightened.


Parr coned pressure fittings are made in two sizes, identified as types "A" and "B", for 1/4-inch and 3/8-inch tubing, respectively.

Coned Pressure Fittings			
Туре	Α	В	
Tube Size	1/4"	3/8"	
Thread (D)	9/16 -18	3/4 -16	
Bore (E)	3/32"	1/8"	
Nut No.	35HC	326HC	
Sleeve No.	40HC	366HC	

Notes



Notes



8

Notes

